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- 1.Burtis CA, Ashwood ER. Textbook of clinical chemistry. Philadelphia: W.B. Saunders Co., 1994:58-69 (techniques and procedures to minimize laboratory infections), 75-76 (specimen collection and

- storage recommendations).
2.Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
3.Clinical and Laboratory Standards Institute /NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4.Wu A, et al. National Academy of Clinical Biochemistry Standards of Laboratory Practice: Recommendations for the Use of Cardiac Markers in Coronary Artery Diseases. Clin Chem 1999; 45(7): 1104-1121.
5.Larue C, et al. Cardiac-specific immunoenzymometric assay of troponin I in the early phase of acute myocardial infarction. Clinical Chemistry 1993;39(6):972-979.
6.Bodor GS. Cardiac troponin-I: A highly specific marker for myocardial infarction. Journal of Clinical Immunoassay 1994;17(1):40-44.
7.Adams JE, et al. Cardiac troponin-I: A marker with high specificity for cardiac injury. 1993;Circulation 88(1):101-106.
8.Adams JE, et al. Diagnosis of perioperative myocardial infarction with measurement of cardiac troponin-I. New England Journal of Medicine 1994;330(10):67-74.
9.Antman EM, et al. Cardiac specific troponin-I levels to predict the risk of mortality in patients with acute coronary syndromes. New England Journal of Medicine 1996;335(18):1342-1349.
10.Braunwald E, et al. ACC/AHA Practice Guidelines: ACC/AHA Guideline Update for the Management of Patients With Unstable Angina and Non–ST-Segment Elevation Myocardial Infarction. American College of Cardiology and American Heart Association 2002.
11.Heeschen C, et al. Evaluation of a Rapid Whole Blood ELISA for Quantification of Troponin I in Patients with Acute Chest Pain. Clin Chem 1999; 45(10):1789-1796.
12.Newby K, et al. Bedside Multimarker Testing for Risk Stratification in Chest Pain Units Circulation 2001; 103(14):1832-1837.
13.Caragher T, et al. Evaluation of Qualitative Cardiac Biomarker Point-of-Care Testing in the Emergency Department. Journal of Emergency Medicine 2002; 22(1).
14.Chapelle JP, et al. Comparison of Cardiac Troponin I Measurements on Whole Blood and Plasma on the Stratus® CS Analyzer and Comparison with AxSYM. Clin Chem 2000; 46(11):1864-1866.
15.Cummins B, Auckland MS, Cummins P. Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction. American Heart Journal 1987;113(6):1333-1344.
16.Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in suspected cases of myocardial infarction. Clinical Chemistry 1992;38(11):2203-2214.
17.World Health Organization. Report on the Joint International Society and Federation of Cardiology-World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria for diagnosis of ischemic heart disease. Circulation 1979;59:607-609.
18.Dade Stratus® Cardiac Troponin-I Fluorometric Enzyme Immunoassay product insert, rev E. 5/97.
19.Myocardial Infarction Redefined – A Consensus Document of The Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction.
20.Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. Clin Chem 1999; 45:7:942-956.
21.Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal IgG. Clin Chem 1992; 38:1737-1742.
22.1. Zot AS, Potter JD, Structural aspects of troponin-tropomyosin regulation of skeletal muscle contraction. Ann. Rev. Biophys.1987;Chem. 16:535-559.
23.Protection of laboratory workers from instrument biohazards and infectious disease transmitted by blood. Approved Guideline; NCCLS Document M29-A, 1997.

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