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- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition.* CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Tietz, NW, Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1986, pp. 52-53 (techniques and procedures to minimize laboratory infections) and pp 478-497 (specimen collection and storage recommendations).
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition.* CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Tietz NM. Clinical Guide to Laboratory Tests, 2nd Edition, W.B. Saunders Co., Philadelphia, PA 1990; pp 554.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition.* CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Kaplan AK, Amadeo JP, Clinical Chemistry: Theory, Analysis, and Correlation, 3rd Edition, Mosby, St. Louis Missouri, 1996, pp. 680.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- National Cholesterol Education Program: Third Report of the Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). Journal of American Medical Association, May 16, 2001; p 18.
- Rifai, N., Russell Warnick, G., Dominiczak, M.H., Handbook of Lipoprotein Testing, AACC Press, Washington, DC 1997; p. 115.
- Burtis, C.A., Ashwood, E. R., Tietz Textbook of Clinical Chemistry, 2nd Edition, W.B. Saunders Co., Philadelphia, PA 1994; p 1017.
- Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline- Second Edition.* CLSI/NCCLS Document EP5-A2 [ISBN 1- 56238-542-9] CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline- Second Edition.* CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

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