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- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
- Clinical and Laboratory Standards Institute. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2000. CLSI Document C28-A2.
- Inloes R, Clark D, Drobenies A. Interference of fluorescein, used in retinal angiography, with certain clinical laboratory tests. *Clin Chem*. 1987;33 (11) :2126-2127.
- Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45 (7) :942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F (ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38 (9) :1737-1742.
- Mackey SE, Creasman WT. Ovarian cancer screening. *J Clin Oncol*. 1995;13 (3) :783-793.
- Schwartz MK. Cancer markers. In: DeVita VT, Hellman S, Rosenberg SA, eds. *Cancer: Principles and Practice of Oncology*. Philadelphia: JB Lippincott Company, 1993;534.
- Lavin PT, Knapp RC, Malkasian G, et al. CA 125 for the monitoring of ovarian carcinoma during primary therapy. *Obstet Gynecol*. 1987;69 (2) :223-227.
- O'Brien TJ, Hardin JW, Bannon GA, et al. CA 125 antigen in human amniotic fluid and fetal membranes. *Am J Obstet Gynecol*. 1986;155 (1) :50-55.
- Kenemans P, Verstraeten AA, van Kamp GJ, von Mensdorff-Pouilly S. The second generation CA 125 assays. *Ann Med*. 1995;27 (1) :107-113.

- Bonfrer JM, Baan AW, Jansen E, Lentfer D, Kenemans P. Technical evaluation of three second generation CA 125 assays. *Eur J Clin Chem Clin Biochem*. 1994;32 (3) :201-207.
- Kenemans P, van Kamp GJ, Oehr P, Verstraeten RA. Heterologous double-determinant immunoradiometric assay CA 125 II: reliable second-generation immunoassay for determining CA 125 in serum. *Clin Chem*. 1993;39 (12) :2509-2513.
- Parker D, Bradley C, Bogle SM, et al. Serum albumin and CA 125 are powerful predictors of survival in epithelial cancer. *Br J Obstet Gynaecol*. 1994;101 (10) :888-893.
- Niloff JM, Knapp RC, Lavin PT, et al. The CA 125 assay as a predictor of clinical recurrence in epithelial ovarian cancer. *Am J Obstet Gynecol*. 1986;155 (1) :56-60.
- Seveldal P, Schemper M, Spona J. CA 125 as an independent prognostic factor for survival in patients with epithelial ovarian cancer. *Am J Obstet Gynecol*. 1989;161 (5) :1213-1216.
- Berek JS, Knapp RC, Malkasian GD, et al. CA 125 serum levels correlated with second-look operations among ovarian cancer patients. *Obstet Gynecol*. 1986;67 (5) :685-689.
- Niloff JM, Bast RC, Schaetzl EM, et al. Predictive value of CA 125 antigen levels in secondlook procedures for ovarian cancer. *Am J Obstet Gynecol*. 1985;151 (7) :981-986.
- Helzlsouer KJ, Bush TL, Alberg AJ, et al. Prospective study of serum CA-125 levels as markers of ovarian cancer. *JAMA*. 1993;269 (9) :1123-1126.
- Mogensen O. Prognostic value of CA 125 in advanced ovarian cancer. *Gynecol Oncol*. 1992;44 (3) :207-212.
- Malkasian GD, Knapp RC, Lavin PT, et al. Preoperative evaluation of serum CA 125 levels in premenopausal and postmenopausal patients with pelvic masses: discrimination of benign from malignant disease. *Am J Obstet Gynecol*. 1988;159 (2) :341-346.
- Chou CY, Chang CH, Yao BL, et al. Color Doppler ultrasonography and serum CA 125 in the differentiation of benign and malignant ovarian tumors. *J Clin Ultrasound*. 1994;22 (8) :491-496.
- Cannistra SA. Cancer of the ovary. *N Engl J Med*. 1993;329 (21) :1550-1559.
- Bast RC, Klug TL, St. John E, et al. A radioimmunoassay using a monoclonal antibody to monitor the course of epithelial ovarian cancer. *N Engl J Med*. 1983;309 (15) :883-887.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
- Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP07-ed3.
- Wild D, ed. *The Immunoassay Handbook*. 4th ed. Oxford, UK: Elsevier Science; 2013:842.

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