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- Tietz NW, Textbook of Clinical Chemistry, W.B. Saunders Co. Philadelphia, PA 1999: pp. 42-72.
- Woodrum D, French C, Shamel LB. Stability of free Prostate-Specific Antigen in serum samples under a variety of sample collection and sample storage conditions. Urology 1996; 48(6A): 33-39.
- Jung K, et al. Preanalytical Determinants of Total and Free Prostate-Specific Antigen and Their Ratio: Blood Collection and Storage Conditions. Clin Chem 1998; 44:685-688.
- Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition. CLSI/NCCLS document H1-A5 (ISBN 1-56238-519-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 (ISBN 1-56238-555-0). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')2 conjugate and polyclonal mouse IgG. Clin Chem 1992; 38:1737-42.
- Fricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. Clin Chem 1999; 45:7,942-956.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the "hook effect in two-site immunometric assays". Anal Biochem 1982; 127:308-15.
- Carroll P, et al: Prostate-Specific Antigen best practice policy- Part I: Early detection and diagnosis of prostate cancer. Urology 2001; 57: 217-224.
- Coley CM, Barry MJ, Fleming C, et al: Early detection of prostate cancer: I. Prior probability and effectiveness of tests. ANN Intern Med 1997; 126:394-406.
- Lilja H. A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. J. Clin Invest 1985;76:1899-903.
- Watt KWK, Lee P-J, Timkulu T, Chan W-P, Loor R. Human prostate-specific antigen: structural and functional similarities with serine protease. Proc Natl Acad Sci. 1986; 83:3166-3170.
- Stenman V-H, Leinonen J, Alfthan H, et al. A complex between prostate-specific antigen and α 1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer: Assay of the complex improves clinical sensitivity for cancer. Cancer Res 1991;51:222-6.
- Lilja H, Christensson A, Dahlen V, et al. Prostate-specific antigen in human serum occurs predominantly in complex with α 1-antichymotrypsin. Clin Chem 1991;37:1618-25.
- Catalona WJ, et al. Use of the Percentage of Free Prostate-Specific Antigen to Enhance Differentiation of Prostate Cancer from benign Prostatic Disease. JAMA 1998; 279: 1542-1547.
- Potter SR, Partin AW. The clinical utility, benefits, and limitations of percent free PSA testing. Contemporary Urology 1999;15-27.
- Clinical and Laboratory Standards Institute /NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

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