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SIEMENS

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- 1.Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline* – Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- 2.Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard* – Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- 3.Clinical and Laboratory Standards Institute/NCCLS. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline* – Second Edition. CLSI document C28-A2 [ISBN 1-56238-406-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA 2000.
- 4.van Beek EJR, et al. A comparative analysis of D-dimer assays in patients with clinically suspected pulmonary embolism. *Thromb. Haemost.* 1993; 70:408.
- 5.Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999; 45:7;942-956.

- 6.Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and Polyclonal mouse IgG. *Clin Chem* 1992; 38:1737-1742.
- 7.Budzynski AZ, et al. Antigenic markers on fragment DD, a unique plasmonic derivative of human crosslinked fibrin. *Blood* 1979, 54: 794.
- 8.Bongard O, et al. D-dimer plasma measurement in patients undergoing major hip surgery: use in the prediction and diagnosis of postoperative proximal vein thrombosis. *Thromb. Res.* 1994; 74:487.
- 9.Bounameaux H, et al. D-dimer testing in suspected venous thromboembolism: an update. *Q J Med* 1997; 90: 437.
- 10.Bovill EG. Disseminated intravascular coagulation: pathophysiology and laboratory diagnosis. *Fibrinolysis* 1993; 7:17.
- 11.Houdijk WP. Utility of D-Dimer Assays in Strategies to Diagnose Venous Thromboembolism (VTE): Selection Criteria for Optimal Performance. *Clinical Hemostasis Review.* 2000: 14 (6):6.
- 12.Jensen R. The Clinical Utility of the D-Dimer Assay. *Clinical Hemostasis Review.* 2000: 14 (6):1.
- 13.Brotman DJ, et al. Limitations of D-dimer Testing in Unselected Inpatients with Suspected Venous Thromboembolism. *The American Journal of Medicine.* 2003; 276-282
- 14.Ver Elst K, et al. Plasma D-Dimer Concentrations in Different Clinical Conditions. *Acta Clinica Belgica.* 2002 ;57 : 325-330.
- 15.Schreengost E, et al. Comparison of Diagnostic Accuracies in Outpatients and Hospitalized Patients of D-Dimer Testing for the Evaluation of Suspected Pulmonary Embolism. *Clinical Chemistry.* 2003 49:9 1483-1490.
- 16.P.S. Wells, et al. Value of Assessment of Pretest probability of Deep Vein Thrombosis in Clinical Management. *Lancet.* 1997; 350:1795 – 1798.
- 17.P.S. Wells, et al. Derivation of a Simple Clinical Model to Categorize Patients Probability of Pulmonary Embolism: Increasing the Models Utility with the SimpliRED D-Dimer. *Thromb Haemost.* 2000; 83: 416-420.
- 18.Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline. CLSI document M29-A [ISBN 1-56238-567-4]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 1997.
- 19.Burtis CA, Ashwood ER. Textbook of clinical chemistry. Philadelphia: W.B. Saunders Co., 1994:58–69 (techniques and procedures to minimize laboratory infections), 75-76 (specimen collection and storage recommendations).
- 20.Wells PS. The role of qualitative D-Dimer assays, clinical probability, and non invasive imaging tests for the diagnosis of deep vein thrombosis and pulmonary embolism. *Semin Vasc Med* 2005; 5: 340-50.
- 21.Kelly J, Hunt BJ. The utility of pretest probability assessment in patients with clinically suspected venous thrombo-embolism. *J Thromb Haemost* 2003; 1: 1888-96.
- 22.Taylor FB, Toh CH, Hoots WK, et al. Towards definition, clinical and laboratory criteria, and a scoring system for disseminated intravascular coagulation. *Thromb Haemost* 2001;86:1327-30.
- 23.Palareti G, Cosmi B. Predicting the risk of recurrence of venous thromboembolism. *Curr Opin Hematol* 2004;11:192-7.
- 24.Lowe GDO. Fibrin D-Dimer and cardiovascular risk. *Semin Vasc Med* 2005; 5: 387-98.
- 25.Eichinger S. D-Dimer in pregnancy. *Semin Vasc Med* 2005; 5: 375-8.

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