L2KPP2

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- 1) Logue FC, Fraser WD, O'Reilly DSJ, Beastall GH. The circadian rhythm of intact parathyroid hormone (1-84) and nephrogenous cyclic adenosine monophosphate in normal men. J Endocrinol 1989;121:R1-R3.
- 2) UK National External Quality Assurance Scheme (NEQAS) for PTH, ACTH, and Calcitonin. Distribution 65, 2001.
- 3) Bowers GN, Brassard C, Sena SF. Measurement of ionized calcium in serum with ion-selective electrodes. Clin Chem 1986;32:1437-47.
- 4) Halvorsen JF, Gautvik VT, Gautvik KM. Improved diagnosis of primary hyperparathyroidism by defining the inverse relationship between serum immunoreactive parathyroid hormone and calcium. Scand J Clin Lab Invest 1986;46:435-42.
- 5) Hawker CD, Clark SW, Martin KJ, Slatopolsky E, Di Bella FP. Radioimmunoassay of parathyroid hormone: clinical utility and interpretation. In: Thompson NW, Vinik AL, editors. Endocrine surgery update. New York: Grune & Statton, 1983:321-40.
- 6) Boyd JC, Ladenson JH. Value of laboratory tests in the differential diagnosis of hypercalcemia. Am J Med 1984; 77:863-72.
- 7) Groth TL, Ljunghall S, De Verdier C-H. Optimal screening for patients with hyperparathyroidism with use of serum calcium observation; a decision-theoretical analysis. Scand J Clin Lab Invest 1983;43:699-707.
- 8) Kleeman K, Kleeman CR. Parathyroid hormone and calcitonin. Contemp Endocrinol 1985;2:247-344.
- 9) Orwoll ES, Meier D. Alterations in calcium, vitamin D, and parathyroid hormone physiology in normal men with aging: relationship to the development of senile osteopenia. J Clin Endocrinol Metab 1986;63:1262-9.
- 10) Kao PC. Parathyroid hormone assay. Mayo Clin Proc 1982;57:596-7.
- 11) Heath H. Tests of parathyroid function: utility and limitations. American Association for Clinical Chemistry: Continuing Education in Endocrinology and Metabolism. Endo (Feb) 1984;2 (8):10 pages.
- 12) Johnson DR, Howarth AT. Differential diagnosis of hypercalcemia. CRC Crit Rev Clin Lab Sci 1984;21:51-97.
- 13) Audran M, Kumar R. The physiology and pathophysiology of vitamin D. Mayo Clin Proc 1985;60:851-66.
- 14) Boscatto LM, Stuart MC. Heterophilic antibodies : a problem for all immunoassays. Clin Chem 1998;34:27-33.
- 15) 社内資料
- * 16) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8.

- * 17) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
- * 18) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

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