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1. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et.al., Luminescent oxygen channeling assay (LOCI): sensitive,broadly applicable homogeneous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
2. Ullman EF, Kirakossian H, Sharat S, Ping Wu,Z, Irvin BR, et.al Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl. Acad.Sci.USA, Vol 91, pp.5426-5430, June 1994 Biochemistry
3. Ullman EF. Homogeneous Immunoassays The immunoassay Handbook,2nd ed., David Wild editor,2001 192-194.
4. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition*. CLSI/NCCLS document H3-A5[ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition*. CLSI /NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400 Wayne, PA 19087-1898 USA. 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898vUSA, 2004.
7. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, Second Edition, AACC Press, Washington, DC 1997; PP.3-386.
8. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
9. Brogan GX, et al. Evaluation of a new rapid quantitative immunoassay for serum myoglobin versus ck-mb for ruling out acute myocardial infarction in the emergency department. Annals of Emerg. Med.1994;4:665-671.
10. Carraro P, et al. A new immunoassay for the measurement of myoglobin in serum. Journal of Clinical Laboratory Analysis 1994;8:70-75.
11. Braunwald E, Antman EM, etal. ACC/AHA guidelines for the management of patients with unstable angina and non-ST-segment elevation myocardial infarction: executive summary and recommendations: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on Management of Patients with Unstable Angina). Circulation 2000;102: 1193-1209.
12. de Winter RJ, et al. Value of Myoglobin, Tropt, CK-MB mass in ruling out an acute myocardial infarction in the emergency room. Circulation 1995;3401-3407.
13. Newby LK, Storrow AB, Gibler B, et al. Bedside multimarker testing for risk stratification in chest pain units: The chest pain evaluation by Creatine Kinase-MB, myoglobin, and Troponin I (CHECKMATE) Study. Circulation 2001;103:1832-1837.
14. Montague C, Kircher T. Myoglobin in the early evaluation of acute chest pain. A.J.C.P. 1995;4:472-476.
15. Wu HB, Apple FS, Warshaw MM. Recommendations for the use of cardiac markers in coronary artery diseases. National Academy of Clinical Biochemistry Standard of Laboratory Practice1999; (5):17-18.
16. Thygesen K, Alpert JS, et al. Myocardial Infarction Redefined—A Consensus Document of The Joint European Society of Cardiology/ American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Journal of the American College of Cardiology 2000;36 (3):959-969.
17. Kricka LJ, Human anti-animal antibody interference in immunological assays. ClinChem 1999; 45 (7): 942-956.
18. Ryall RG, Story CJ, Turner DR. Reappraisal of the causes of the "hook effect" in two-site immunoradiometric assays. Anal Biochem 1982:127:308.
19. Tucker JF, et al. Value of serial myoglobin levels in the early diagnosis of patients admitted for acute myocardial infarction. Annals of Emerg Med 1994;4: 704-708.
20. Woo J, et al. Is myoglobin useful in the diagnosis of acute myocardial infarction in the emergency department setting? Am J Clin.Pathol. 1995;103: No.6;725-729.
21. Clinical and Laboratory Standards Institute /NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2[ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
22. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

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