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1. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et.al., Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogenous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
2. Ullman EF, Kirakossian H, Sharat S, Ping Wu,Z, Irvin BR, et.al Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl. Acad.Sci.USA, Vol 91, pp.5426-5430, June 1994 Biochemistry
3. Ullman,EF Homogenous Immunoassays The immunoassay Handbook,2nd ed., David Wild editor,2001 192-194.
4. Tietz, NM. Textbook of Clinical Chemistry, 3rd ed.(edited by).Carl A. Burtis, Edward R. Ashwood, W.B. Saunders Co., Philadelphia, PA 1999:662.
5. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition. NCCLS document H3-A5 [ISBN 1-56238-515-1]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2003.
6. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2004.
7. Clinical and Laboratory Standards Institute /NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition; NCCLS document H1-A5 [ISBN 1-56238-519-4]. NCCLS, 940 West Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898 USA. 2003.
8. Clinical and Laboratory Standards Institute /NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI/ NCCLS document EP7-A2 [ISBN 1-56238-584-4]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898, 2005.
9. Dimension® clinical chemistry system, MMB Instructions for Use, 2002-10 E.
10. Skogen W, Osborne D, Price G. Guidelines for Interpreting the CK-MB Relative Index Value. Clin Chem 1991; 37:916.
11. Wu HB, Apple FS, Warshaw MM. Recommendations for the use of cardiac markers in coronary artery diseases. National Academy of Clinical Biochemistry Standard of Laboratory Practice1999; (5):17-18.
12. Thygesen K, Alpert JS, et al. Myocardial Infarction Redefined—A Consensus Document of The Joint European Society of Cardiology/ American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Journal of the American College of Cardiology 2000; 36 (3):959-969.
13. Marin MM, Teichman S L. Use of rapid serial sampling of creatine kinase MB for very early detection of myocardial infarction in patients with acute chest pain. Am Heart 1992; 123:354-361.
14. Sharkey SW, Brunette DD, Ruiz E, Hession WT et al. An Analysis of Time Delays Preceding Thrombolysis for Acute Myocardial Infarction. J Am Med Assoc 1989;262:3171-3174
15. Apple FS. Creatine Kinase-MB. Lab Med 1992; 23:298-302.
16. Apple FS, Sharkey SW, Werdick M, Elsperger KJ, Tilbury RT. Analysis of Creatine Kinase Isoenzymes and Isoforms in Serum to Detect Reperfusion after Acute Myocardial Infarction. Clin Chem 1987; 33:507-511.
17. Galen RS, Gambino SR. Isoenzymes of CPK and LDH in Myocardial Infarction and Certain Other Diseases. Pathobiol Ann 1975; 5:283-315.
18. Kricka LJ, Human anti-animal antibody interference in immunological assays. Clin Chem 1999; 45 (7): 942-956.
19. Lott JA. Serum Enzyme Determinations in the Diagnosis of Acute Myocardial Infarction: An Update. Hum Pathol 1984; 15:706-716.
20. Gulbis B. Mass Concentration of Creatine Kinase MB Isoenzyme and Lactate Dehydrogenase Isoenzyme 1 in Diagnosis of Perioperative Myocardial Infarction after Coronary Bypass Surgery. Clin Chem 1990; 36:1784-1788.
21. Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline -- Second Edition. CLSI/NCCLS Document EP5-A2 [ISBN 1-56238-542-9]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
22. Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI/NCCLS document EP9-A2 (ISBN 1-56238-472-4). NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
23. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline*. CLSI Document EP17-A2 [ISBN 1-56238-796-0].CLSI,940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2012.
- **24. Grimsey P, Frey N, Bendig G, *et al*. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet*. 2017;2(4):247-256.
- **25. Piketty ML, Prie D, Sedel F, *et al*. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med*. 2017;55(6):817-825.

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