Dimension® EXL™ integrated chemistry system
LOCI® Module

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- Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et. al., Luminescent oxygen channeling assay (LOC®): sensitive, broadly applicable homogenous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
- Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc. Natl. Acad. Sci. USA, Vol 91 pp.5426-5430, June 1994 Biochemistry.
- Ullman, EF Homogenous Immunoassays The immunoassay Handbook, 2nd ed., David Wild editor, 2001 192-194.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1- 56238-555-0]. CLSI, 940 West Valley Road, Suite1400, Wayne, PA 19087-1898, 2004.
- Stability of Human Cardiac Troponin-I in Refrigerated Plasma Samples. Dade GLS '04 – 48.
- Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Cummins B, Auckland MS, Cummins P, Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction. American Heart Journal 1987; 113 (6): 1333-1344.
- Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in suspected cases of myocardial infarction. Clinical Chemistry 1992; 38(11):2203-2214.
- Larue C, et al., Cardiac-specific immunoenzymometric assay of troponin-I in the early phase of acute myocardial infarction. Clinical Chemistry 1993; 39(6): 972-979.
- National Academy of Clinical Biochemistry Laboratory Medicine Practice: clinical characteristics and utilization of biochemical markers in acute coronary syndromes. Clin Chem 2007; 53: 552-574.
- Bodor GS, Cardiac troponin-I: A highly specific marker for myocardial infarction. Journal of Clinical Immunoassay 1994; 17(1): 40-44.
- Adams JE, et al., Cardiac troponin-I: A marker with high specificity for cardiac injury. 1993; Circulation 88(1): 101-106.
- Adams JE, et al., Diagnosis of perioperative myocardial infarction with measurement of cardiac troponin-I. New England Journal of Medicine 1994; 330 (10): 67-674.
- Antman EM, et al., Cardiac specific troponin-I levels to predict the risk of mortality in patients with acute coronary syndromes. New England Journal of Medicine 1996; 335(18): 1342-1349.
- Braunwald, E. et al. ACC/AHA Practice Guidelines: ACC/AHA Guideline Update for the Management of Patients With Unstable Angina and Non-ST Segment Elevation Myocardial Infarction. American College of Cardiology and American Heart Association 2002.
- National Academy of Clinical Biochemistry and IFCC Committee for standardization of markers of cardiac damage laboratory medicine practice guidelines: analytical issues for biochemical markers of acute coronary syndromes. Clin Chem 2007; 53: 547-551.
- Kricka LJ, Human anti-animal antibody interference in immunological assays. Clin Chem 1999; 45 (7): 942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')2 conjugate and polyclonal IgG. Clin Chem 1992; 38:1737-1742.
- Ryall RG, Story CJ, Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays. Anal Biochem 1982:127:308.
- Zot AS, Potter JD, Structural aspects of troponin-tropomyosin regulation of skeletal muscle contraction. Ann Rev Biophys 1987; Chem. 16: 535-559.
- Myocardial Infarction redefined – A consensus Document of the Joint European Society of Cardiology/American College of Cardiology committee for the Redefinition of Myocardial Infarction. J Am Coll Cardiol 2000, 36 (3): 959 – 69.
- Braunwald E et al. ACC/AHA 2002 guideline update for the management of patients with unstable angina and non-ST-segment elevation myocardial infarction – summary article: a report of the American College of Cardiology/American Heart Association task force on practice guidelines (Committee on the Management of Patients with Unstable Angina). J Am Coll Cardiol 2002; 40: 1366-1374.
- Jaffe AS et al. Biomarkers in acute cardiac disease – the present and the future. J Am Coll Cardiol 2006; 48: 1-11.
- Waxman DA et al. A model for troponin-I as a quantitative predictor of in-hospital mortality. J Am Coll Cardiol 2006; 48: 1755-1762.
- Kontos MC et al. Implication of different cardiac troponin-I levels for clinical outcomes and prognosis of acute chest pain patients. J Am Coll Cardiol 2004; 43: 958-965.
- Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition. CLSI/NCCLS Document EP5-A2 [ISBN 1- 56238-542-9] CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. Int. J. Pharmacokinet. 2017;2(4):247–256.
- Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles:validation of a simple method to overcome biotin interference. Clin Chem Lab Med. 2017;55(6):817-825.

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