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- Clinical and Laboratory Standards Institute (CLSI). Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard-Fifth Edition. CLSI document H3-A5 (ISBN 1-56238-515-1). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute (CLSI). Tubes and additives for venous blood specimen collection; Approved Standard-Fifth Edition. CLSI document H1-A5 (ISBN 1-56238-519-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute (CLSI). Procedures for the handling and processing of blood specimens; Approved Guideline-Third Edition. CLSI document H18-A3 (ISBN 1-56238-555-0). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Schreiber W, Rausch S, Lammers M. Choice of specimen for immunonephelometric protein assays. Clin Chim Acta 2005;355:S407.
- Clinical and Laboratory Standards Institute (CLSI). Interference testing in clinical chemistry; Approved Guideline-Second Edition. CLSI document EP7-A2 (ISBN 1-56238-584-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- National Kidney Foundation. K/DOQI clinical practice guidelines for chronic kidney disease: evaluation, classification and stratification. Am J Kidney Dis 2002;39 (2 Supp 1):S1-266.
- Mathew TH; Australasian Creatinine Consensus Working Group. Chronic kidney disease and automatic reporting of estimated glomerularfiltration rate: a position statement. Med J Aust 2005;183:138-41.
- Levey AS, Eckardt KU, Tsukamoto Y, et al. Definition and classification of chronic kidney disease: a position statement from Kidney Disease: Improving Global Outcomes (KDIGO). Kidney Int 2005;67:2089-100.
- Hoek FJ, Kemperman FAW, Krediet RT. A comparison between cystatin C, plasma creatinine and the Cockcroft and Gault formula for the estimation of glomerular filtration rate. Nephrol Dial Transplant 2003;18:2024-31.
- Zegers I, Bilrup-Jensen S, Althaus H, et al. Certification of cystatin C in the human serum reference material ERM®-DA471/IFCC. 2010 report EUR 24408 EN (ISSN 1018-5593). 1-34.
- Simonsen O, Grubb A, Thysell H. The blood serum concentration of cystatin C (γ -trace) as a measure of the glomerular filtration rate. Scand J Clin Lab Invest 1985;45:97-101.
- Randers E, Erlandsen EJ. Serum cystatin C as an endogenous marker of the renal function -a review. Clin Chem Lab Med 1999;37:389-95.
- Grubb A. Diagnostic value of analysis of cystatin C and protein HC in biological fluids. ClinNephrol 1992;38 Suppl 1:S20-7.
- Coll E, Botey A, Alvarez L, et al. Serum cystatin C as a new marker for noninvasiveestimation of glomerular filtration rate and as a marker for early renal impairment. Am J Kidney Dis 2000;36:29-34.
- Clinical and Laboratory Standards Institute (CLSI). Evaluation of precision performance of quantitative measurement methods; Approved Guideline-Second Edition. CLSI document EP5-A2 (ISBN 1-56238-542-9). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute (CLSI). Method comparison and bias estimation using patient samples; Approved Guideline-Second Edition. CLSI document EP9-A2 (ISBN 1-56238-472-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

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