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1. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et. al., Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogeneous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
2. Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al. Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl. Acad. Sci. USA, Vol 91, pp 5426 – 5430, June 1994 Biochemistry.
3. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition. CLSINCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
4. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, Second Edition. AACC Press Washington DC 1997.
5. DiMagno EP, et.al. Effect of Long-term Freezer Storage, Thawing, and Refreezing on Selected Constituents of Serum. Mayo Clin Proc, October 1989, Vol 64.
6. Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
7. Fraser, Callum G. PhD, Biological Variation: From Principles to Practice, Washington, DC: AACC Press, 2001 [pp. 71 – 78].
8. Trapé J, Botargés JM, Porta F, Ricós C, Badal JM, Salinas R, Sala M, Roca A. Reference change value for alpha-fetoprotein and its application in early detection of hepatocellular carcinoma in patients with hepatic disease. Clin Chem, 2003; 49:1209-1211.
9. Davison AC, Hinkley DV. Bootstrap methods and their application. Cambridge, NY : Cambridge University Press, 1997.
10. Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays, Anal Biochem 1982;127:308-315.
11. Kricka LJ Human Anti-Animal Antibody Interferences in Immunological Assays. Clin Chem 1999; 45:7; 942-956.
12. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')2 conjugate and polymouse IgG. Clin Chem 1992; 38:1737-1742.
13. Abelev GI, Perova SD, Kharmkova NI et al. Production of Embryonal α-Globulin by Transplantable Mouse Hepatomas. Transplantation 1963; 1: 174-180.
14. Abelev GI. Production of Embryonal Serum α-Globulin by Hepatomas: Review of Experimental and Clinical Data. Cancer Research 1968; 28: 1344-1350.
15. Bergstrand CG, Czar B. Demonstration of a New Protein Fraction in Serum from the Human Fetus. Scand. J. Clin Lab Invest. 1956; 8: 74.
16. Chan DW, Booth RA, Diamandis EP. Tumor Markers in Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, Fourth Edition 2006; 767-768.
17. Bosl GJ, Bajorin DF, Shienfeld J et al. Cancer of the Testes in DeVita VT Jr, Hellman S, Rosenberg SA eds Cancer: Principles and Practices of Oncology, Seventh Edition, Philadelphia PA, Lippincott Williams and Wilkins 2005; 1269-90.
18. Bartlett NL, Freiha, FF, Torti FM. Serum Markers in Germ Cell Neoplasms. Hem./Onc. Clinics of N.A. 1991; 5: 1245-1260.
19. Clinical and Laboratory Standards Institute /NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI/NCCLS document EP5-A2 [ISBN1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
20. Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

21. Clinical and Laboratory Standards Institute/NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- *22. Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet.* 2017;2(4):247–256.
- *23. Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med.* 2017;55(6):817-825.

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