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- Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al. Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogenous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
- Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et al. Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl. Acad. Sci. USA, Vol 91, pp 5426 – 5430, June 1994 Biochemistry.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition. CLSI/NCCLS document H3-A5 [ISBN1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, Second Edition. AACC Press Washington DC 1997.
- Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Fraser Callum G. PhD, Biological Variation: From Principles to Practice, Washington, DC: AACC Press, 2001 [pp. 71 – 78].
- Ricos C, Alvarez V, Cava F, Garcia-Lario JV, Hernandez A, Jimenez CV, Minchinela J, Perich C, Simon M. “Current databases on biologic variation: pros, cons and progress.” Scand J Clin Lab Invest 1999; 59: 491-500.<http://www.westgard.com/biodatabase1.html>
- B. Emir, S. Wieand, John Q.S., and S. Cha. “Analysis of repeated markers used to predict progression of cancer”. Statistics in Medicine. 1998 [2563-2578].
- Davison AC, Hinkley DV. Bootstrap methods and their application. Cambridge, NY : Cambridge University Press, 1997.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays, Anal Biochem 1982;127:308-315.
- Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays Clin Chem 1999; 45:7; 942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')2 conjugate and polyclonal IgG. Clin Chem 1992; 38:1737-1742.
- Gold P, Freedman SO. Demonstration of Tumor Specific Antigens in Human Colonic Carcinomata by Immunological Tolerance and Absorption Techniques. J Exp Med.1965; 121: 439-462.
- Chan DW, Booth RA, Diamandis EP. Tumor Markers in Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, Fourth Edition. 2006; 768-769.
- Paus E, Almasbak H et al. A Single-chain Fv-based Immunofluorometric Assay Specific for the CEA Variant NCA-2. J Immunol Methods. 2003; 286: 125-139.
- Tietz Clinical Guide to Laboratory Tests, Fourth Edition. 2006; 220-222.
- NACB: Practice Guidelines And Recommendations For Use Of Tumor Markers In The Clinic Quality Requirements [Section 2]. National Academy of Clinical Biochemistry Guidelines on Quality Requirements for the Use of Tumor Markers. Catharine Sturgeon, Elizabeth Hammond, Soo-Ling Ch'ng, György Sölétormos, Daniel F Hayes.
- Perkins GL, Slater ED et al. Serum Tumor Markers. Am Family Physician.2003; 68: 1075-1082.
- Goldstein MJ, Mitchell EP. Carcinoembryonic Antigen in the Staging and Follow-up of Patients with Colorectal Cancer. Cancer Invest.2005; 23: 338-351.
- Fakih MG. CEA Monitoring in Colorectal Cancer. Oncology. 2006; 20.
- Clinical and Laboratory Standards Institute /NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI/NCCLS document EP5-A2 [ISBN1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Clinical and Laboratory Standards Institute/NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8].CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- **24. Grimsey P, Frey N, Bendig G, *et al.* Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet.* 2017;2(4):247–256.
- **25. Piketty ML, Prie D, Sedel F, *et al.* High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. *Clin Chem Lab Med.* 2017;55(6):817-825.

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