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- Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et. al., Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogenous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
- Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl.Acad.Sci. USA, Vol 91, pp.5426-5430, June 1994 Biochemistry.
- Tietz NW, Textbook of Clinical Chemistry, W.B. Saunders Co. Philadelphia, PA 1999: pp. 42-72.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 (ISBN 1-56238-519-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Woodrum D, French C, Shamel LB. Stability of free Prostate-Specific Antigen in serum samples under a variety of sample collection and sample storage conditions. Urology 1996; 48(6A): 33-39.
- Jung K, et al. Preanalytical Determinants of Total and Free Prostate-Specific Antigen and Their Ratio: Blood Collection and Storage Conditions. Clin Chem 1998; 44: 685-688.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline-Third Edition*. CLSI/NCCLS document H18-A3 (ISBN 1-56238-555-0). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook”effect in two-site immunometric assays. Anal Biochem 1982; 127: 308-15.
- Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. Clin Chem 1999; 45:7; 942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab)2 conjugate and polyclonal mouse IgG. Clin Chem 1992; 38: 1737-42.
- www.cancer.org
- Coley CM, Barry MJ, Fleming C, et al: Early detection of prostate cancer: I. Prior probability and effectiveness of tests. ANN Intern Med 1997; 126: 394-406.
- Lilja H. A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. J Clin Invest 1985;76: 1899-903.
- Watt KWK, Lee P-J, Timkulu T, Chan W-P, Loor R. Human prostate-specific antigen: structural and functional similarities with serine protease. Proc Natl Acad Sci. 1986; 83: 3166-3170.
- Stenman V-H, Leinonen J, Alifthan H, et al. A complex between prostate-specific antigen and α1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer: Assay of the complex improves clinical sensitivity for cancer. Cancer Res 1991;51: 222-6.
- Lilja H, Christensson A, Dahlen V, et al. Prostate-specific antigen in human serum occurs predominantly in complex with α1-antichymotrypsin. Clin Chem 1991; 37: 1618-25.
- Catalona WJ, et al. Use of the Percentage of Free Prostate-Specific Antigen to Enhance Differentiation of Prostate Cancer from benign Prostatic Disease. JAMA 1998; 279: 1542-1547.
- Potter SR, Partin AW. The clinical utility, benefits, and limitations of percent free PSA testing. Contemporary Urology 1999; 15-27.
- Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

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