SIEMENS

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- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimen; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite1400, Wayne, PA 19087-1898 USA, 2004.
- Tietz NW. *Textbook of Clinical Chemistry*. 4th ed., W.B. Saunders Co. Philadelphia, PA 2006: pp 602.
- Wu A, ed. *Tietz Clinical Guide to Laboratory Tests*, 4th ed. Philadelphia; W.B. Saunders Company, 2006: p. 648.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Klein G, Berger A, Bertholf R et al. Abstract: Multicenter evaluation of liquid reagents for CK, CK-MB and LDH with determination of reference intervals on Hitachi Systems. Clin Chem 2001; 47:Suppl. A30.
- Siekmann L, Ceriotti F, Ferard G, Kanno T, Schumann G, et. al. IFCC Primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37° C. Part 3. Reference procedure for the Measurement of Catalytic Concentration of Creatine Kinase. Clin Chem Lab Med. 2002; 40(6):635-642.
- McPherson RA, Pincus MR, eds. *Henry's clinical diagnosis and management by laboratory methods* 21st ed. St. Louis: Saunders Elsevier, 2007: 221, 260.
- Anderson JL, Adams CD, Antman EM, Bridges CR, Califf RM, Casey DE, et. al. ACC/AHA 2007 Guidelines for the management of patients with unstable angina/non-ST-elevation myocardial infarction-executive summary. JACC. 50(7), 2007:652-726.
- Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline- Second Edition*. CLSI/NCCLS Document EP5-A2 [ISBN 1- 56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline- Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

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