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10995644

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1. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
2. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
3. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
4. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
5. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
6. Clinical and Laboratory Standards Institute. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2000. CLSI Document C28-A2.
7. Catalona WJ, Smith DS, Ratliff TL, et al. Measurement of prostate-specific antigen in serum as a screening test for prostate cancer. *N Engl J Med*. 1991;324 (17) :1156-1161.
8. Gormley GJ, Stoner E, Bruskewitz RC, et al. The effect of finasteride in men with benign prostatic hyperplasia. The Finasteride Study Group. *N Engl J Med*. 1992;327 (17) :1185-1191.
9. Graves HC, Wehner N, Stamey TA. Comparison of a polyclonal and monoclonal immunoassay for PSA: need for an international antigen standard. *J Urol*. 1990;144 (6) : 1516-1522.
10. Lilja H, Christensson A, Dahlen U, et al. Prostate-specific antigen in human serum occurs predominantly in complex with alpha 1-antichymotrypsin. *Clin Chem*. 1991;37 (9) :1618-1625.
11. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45 (7) :942-956.
12. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F (ab') 2 conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38 (9) :1737-1742.
13. Wang MC, Valenzuela LA, Murphy GP, Chu TM. Purification of a human prostate specific antigen. *Invest Urol*. 1979;17 (2) :159-163.
14. Watt KW, Lee PJ, M'Timkulu T, et al. Human prostate-specific antigen: structural and functional similarity with serine proteases. *Proc Natl Acad Sci USA*. 1986;83 (10) :3166-3170.
15. Lilja H. A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. *J Clin Invest*. 1985;76 (5) :1899-1903.
16. Christensson A, Laurell CB, Lilja H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine proteinase inhibitors. *Eur J Biochem*. 1990;194 (3) :755-763.
17. Lilja H. Significance of different molecular forms of serum PSA. *Urol Clin North Am*. 1993;20 (4) :681-686.
18. McCormack RT, Rittenhouse HG, Finlay JA, et al. Molecular forms of prostate-specific antigen and the human kallikrein gene family: a new era. *Urology*. 1995;45 (5) :729-744.
19. Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and characterization of different molecular forms of prostate-specific antigen in human seminal fluid. *Clin Chem*. 1995;41 (11) :1567-1573.
20. Chen Z, Chen H, Stamey TA. Prostate specific antigen in benign prostatic hyperplasia: purification and characterization. *J Urol*. 1997;157 (6) :2166-2170.
21. Mikolajczyk SD, Grauer LS, Millar LS, et al. A precursor form of PSA (pPSA) is a component of the free PSA in prostate cancer serum. *Urology*. 1997;50 (5) :710-714.
22. Mikolajczyk SD, Marks, LS, Partin AW, Rittenhouse HG. Free prostate-specific antigen in serum is becoming more complex. *Urology*. 2002;59 (6) :797-802.
23. Stamey TA, Yang N, Hay AR, et al. Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. *N Engl J Med*. 1987;317 (15) :909-916.

24. Schellhammer PF, Schlossberg SM, el-Mahdi AM, et al. Prostate specific antigen levels after definitive irradiation for carcinoma of the prostate. *J Urol.* 1991;145 (5) :1008-1010.
25. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem.* 1988;34 (1) :27-33.
26. Stenman UH, Leinonen J, Alfthan H, et al. A complex between prostate-specific antigen and alpha 1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostatic cancer: assay of the complex improves clinical sensitivity for cancer. *Cancer Res.* 1991;51 (1) :222-226.
27. Zhou Z, Armstrong EG, Belenky A, et al. Equivalent recognition of free and ACT-complexed PSA in a monoclonal-polyclonal sandwich assay is conferred by binding specificity of the monoclonal antibody. *J Clin Lab Anal.* 1998;12 (4) :242-249.
28. Christensson A, Bjork T, Nilsson O, et al. Serum prostate-specific antigen complexed to α -1-antichymotrypsin as an indicator of prostate cancer. *J Urol.* 1993;150 (1) :100-105.
29. Catalona WJ, Smith DS, Wolfert RL, et al. Evaluation of percentage of free serum prostatespecific antigen to improve specificity of prostate cancer screening. *JAMA.* 1995; 274 (5) :1214-1220.
30. Luderer AA, Chen YT, Soriano TF, et al. Measurement of the proportion of free to total prostate-specific antigen improves diagnostic performance of prostate-specific antigen in the diagnostic gray zone of the total prostate-specific antigen. *Urology.* 1995; 46 (2) :187-194.
31. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
32. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
33. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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