Dimension® EXL™ integrated chemistry system LOCI® Module

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1.Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al. *Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogeneous immunoassay method.* Clin Chem 42:9 1996, 1518-1526.

2.Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al. *Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence.* Proc. Natl. Acad. Sci. USA, Vol 91, pp.5426-5430, June 1994 Biochemistry.

3.Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture: Approved Standard—Fifth Edition.* CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.

4.Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition.* CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

5.Shimizu H., Aono K., Masuta K., et al. *Stability of brain natriuretic peptide (BNP) in human blood samples.* Clinica Chim Acta 1999;285:169-72.

6.Shimizu H, Aono K, Masuta K, Asada H, Misaki A, Teraoka H. *Degradation of human brain natriuretic peptide (BNP) by contact activation of blood coagulation system.* Clin Chim Acta 2001;305:181-6.

7.Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition.* CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.

8.Pereira M., Azevedo A., Severo M., Barros H. *Long-term stability of endogenous B-type natriuretic peptide after storage at -20°C or -80°C.* Clin Chem Lab Med 2008;46(8):1171-1174.

9.Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.

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11.Kricka LJ. *Human Anti-Animal Antibody Interferences in Immunological Assays.* Clin Chem 1999; 45:7:942-956.

12.Vaidya HC, Beatty BG. *Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using Fabγ: conjugate and polyclonal mouse IgG.* Clin Chem 1992; 38:1737-1742.

13.Ryall RG, Story CJ, and Turner DR. *Reappraisal of the causes of the “hook effect” in two-site immunoradiometric assays,* Anal Biochem 1982;127:308-315.

14.Remme WJ, Swedberg K. *Task Force Report: guidelines for the diagnosis and treatment of chronic heart failure.* Eur Heart J 2001;22:1527-60.

15.Hunt SA, Baker DW, Chin MH, Cinquegrani MP, Feldman AM, Francis GS, Ganiats TG, Goldstein S, Gregoratos G, Jessup ML, Noble RJ, Packer M, Silver MA, Stevenson LW, Gibbons RJ, Antman EM, Alpert JS, Faxon DP, Fuster V, Jacobs AK, Hiratzka LF, Russell RO, Smith SC Jr. *ACC/AHA guidelines for the evaluation and management of chronic heart failure in the adult: executive summary. A report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines.* J Am Coll Cardiol. 2001;38:2101-13.

16.Stein BC, Levin RJ. *Natriuretic peptides: physiology, therapeutic potential, and risk stratification in ischemic heart disease.* Am Heart J 1998;135:914-23.

17.Sagnella GA. *Measurement and significance of circulating natriuretic peptides in cardiovascular disease.* Clin Science 1998;95:519-29.

18.Suzuki T, Yamazaki T, Yazaki Y. *The role of natriuretic peptides in the cardiovascular system.* Cardiovasc Res 2001;51:489-94.

19.Levin ER, Gardner DG, Samson WK. Natriuretic peptides. N Engl J Med 1998;339:321-8.

20.Wilkins M, Redondo J, Brown L. *The natriuretic peptide family.* The Lancet 1997;349:1307-10.

21.Dickstein K. Natriuretic peptides in detection of heart failure. The Lancet 1998;35:3-4.

22.Chen HH, Burnett JC Jr. *The natriuretic peptides in heart failure: diagnostic and therapeutic potentials.* Proc Assoc Am Physicians 1999;111:406-16.

23.Hirata Y, Matsumoto A, Aoyagi T, Ymaoki K, Kumoro I, Suzuki T, Ashida T, Sugiyama T, Hada Y, Kuwajima I, Nishinaga M, Akioka H, Nakajima O, Nagai R, Yazaki Y. *Measurement of plasma brain natriuretic peptide level as a guide for cardiac overload.* Cardiovascular Res 2001;51:585-91.

24.Mair J, Hammerer-Lercher, Puschendorf B. *The impact of cardiac natriuretic peptide determination on the diagnosis and management of heart failure.* Clin Chem Lab Med 2001;39:571-88.

25.Clerico A, Del Ry S, Giannessi D. *Measurement of cardiac natriuretic hormones (atrial natriuretic peptide, brain natriuretic peptide, and related peptides) in clinical practice: the need for a new generation of immunoassay methods.* Clin Chem 2000;46:1529-34.

26.Yan W, Wu F, Morser J, Wu Q, Corin, A. *transmembrane cardiac serine protease, acts as a pro-atrial natriuretic peptide-converting enzyme.* Proc Natl Acad Sci USA 2000;97:8525-9.

27.Kelly R, Struthers A. Are natriuretic peptides clinically useful as markers of heart failure? Ann Clin Biochem 2001;38:94-102.

28.Daly C, Fox K, Henein M. *Natriuretic peptides in the diagnosis of heart disease – first amongst equals?* Intl J Cardiol 2002;84:107-113.

29.Tabbibizar R, Maisel A. The impact of B-type natriuretic peptide levels on the diagnoses and management of congestive heart failure. Curr Op Cardiol 2002;17:340-45.

30.Lee SC, Stevens TL, Sandberg SM, Heublein BM, Nelson SM, Jougasaki M, Redfield MM, Burnett JC Jr. *The potential of brain natriuretic peptide as a biomarker for New York Heart Association class during the outpatient treatment of heart failure.* J Cardiac Failure 2002;8:149-54.

31.Lubien E, DeMaria A, Krishnaswamy P, Clopton P, Koon J, Kazanegra R, Gardetto N, Wanner E, Maisel AS. *Utility of B-natriuretic peptide in detecting diastolic dysfunction: comparison with doppler velocity recordings.* Circulation 2002;105:595-601.

32.Fisher Y, Filzmaier K, Stiegler H, Graf J, Fuhs S, Franke A, Janssens U, Gressner AM. *Evaluation of a new, rapid bedside test for quantitative determination of B-type natriuretic peptide.* Clin Chem 2001;47:591-4.

33.Maisel AS, Krishnaswamy P, Nowak RM, McCord J, Hollander JE, Duc P, Omland T, Storrow AB, Abraham WT, Wu AHB, Clopton P, Steg PG, Westheim A, Wold Knudsen C, Perez A, Kazanegra R, Herrmann H, McCullough PA. *Rapid measurement of B-type natriuretic peptide in the emergency diagnosis of heart failure.* New Engl J Med. 2002;347:161-7.

34.Dao Q, Krishnaswamy P, Kazanegra R, Harrison A, Amirnovin R, Lenert L, Clopton P, Alberto J, Hlavin P, Maisel AS. *Utility of B-type natriuretic peptide in the diagnosis of congestive heart failure in an urgent-care setting.* J Am Col Cardiol 2001;37:379-85.

35.Morrison LK, Harrison A, Krishnaswamy P, Kazanegra R, Clopton P, Maisel AS. *Utility of a rapid B-type natriuretic peptide assay in differentiating congestive heart failure from lung disease in patients presenting with dyspnea.* J Am Col Cardiol 2002;39:202-9.

36.McCullough PA, Nowak RM, McCord J, Hollander JE, Herrmann H, Steg PG, Duc P, Westheim A, Omland T, Wold Knudsen C, Storrow AB, Abraham WT, Lamba S, Wu AHB, Perez A, Clopton P, Krishnaswamy P, Kazanegra R, Maisel AS. *B-type natriuretic peptide and clinical judgment in emergency diagnosis of heart failure: analysis from the breathing not properly (BNP) multinational study.* Circulation 2002;106:416-22.

37.Wieczorek SJ, Wu AHB, Christenson R, Krishnaswamy P, Gottlieb S, Rosano T, Hager D, Gradetto N, Chiu A, Bailly KR, Maisel AS. *A rapid B-type natriuretic peptide assay accurately diagnoses left ventricular dysfunction and heart failure: a multicenter evaluation.* Am Heart J 2002;144:834-9.

38.Maisel AS, Koon J, Krishnaswamy P, Kazanegra R, Clopton P, Gradetto N, Morrisey NC, Chiu A, DeMaria A. *Utility of B-natriuretic peptide as a rapid, point-of-care test for screening patients undergoing echocardiography to determine left ventricular dysfunction.* Am Heart J 2001;141:367-74.

39.Valli N, Georges A, Corcuff JB, Barat JL, Bordenave L. *Assessment of brain natriuretic peptide in patients with suspected heart failure: comparison with radionuclide ventriculography data.* Clin Chim Acta 2001; 306:19-26.

40.McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, Tunstall-Pedoe H, McMurray JJV, Dargie HJ. *Biochemical detection of left-ventricular systolic dysfunction.* The Lancet 1998;351:9-13.

41.Berger R, Huelsman M, Strecker K, Bojic A, Moser P, Stanek B, Pacher R. *B-type natriuretic peptide predicts sudden death in patients with chronic heart failure.* Circulation 2002;105:2392-97.

42.Koglin J, Pehlivanli S, Schwaiblmair M, Vogeser M, Cremer P, von Scheidt W. *Role of brain natriuretic peptide in risk stratification of patients with congestive heart failure.* J Am Coll Cardiol 2001;38:1934-41.

43.Cheng V, Kazanegra R, Garcia A, Lenert L, Krishnaswamy P, Gardetto N, Clopton P, Maisel A. *A rapid bedside test for B-type natriuretic peptide predicts treatment outcomes in patients admitted for decompensated heart failure: a pilot study.* J Am Coll Cardiol 2001;37:386-91.

44.Stanek B, Frey B, Hulsmann M, Berger R, Sturm B, Strametz-Juranek J, Bergler-Klein J, Moser P, Bjic A, Hartter E, Pacher R. *Prognostic evaluation of neurohumoral plasma levels before and during beta-blockertherapy in advanced left ventricular dysfunction.* J Am Col Cardiol 2001;38:436-42.

45.DeLemos JA, Morrow DA, Bentley JH, Omland T, Sabatine MS, McCabe CH, Hall C, Cannon CP, Braunwald E. *The prognostic value of B-type natriuretic peptide in patients with acute coronary syndromes.* N Engl J Med 2001;345:1014-21.

46.Omland T, Aakvaag A, Bonarjee VVS, Caidahl K, Lie RT, Nilsen DWT, Sundsfjord JA, Dickstein K. *Plasma brain natriuretic peptide as an indicator of left ventricular systolic function and long-term survival after acute myocardial infarction.* Circulation 1996;93:1963-69.

47.Mega JL, Morrow DA, deLemos JA, Sabatine MS, Murphy SA, Rifai N, Gibson CM, Antman EM, Braunwald E. *BNP at Presentation and Prognosis in Patients with ST-Elevation Myocardial Infarction: An ENTIRE-TIMI 23 Substudy.* J Am Col Cardiol. 2004;44:335-9.

48.Morrow DA, de Lemos JA, Blazing MA, Sabatine MS, Murphy SA, Jarolim P, White HD, Fox KA, Califf RM, Braunwald E. *Prognostic value of serial B-type natriuretic peptide testing during follow-up of patients with unstable coronary artery disease.* JAMA. 2005;294(22):2866-71.

49.Richards AM, Lainchbury JG, Nicholls MG, Troughton RW, Yandel TG. *BNP in hormone-guided treatment of heart failure.* Trends Endocrinol Metab 2002;13:151-5.

50.Johnson W, Omland T, Hall C, Lucas C, Myking OL, Collins C, Pfeffer M, Rouleau JL, Stevenson LW. *Neurohormonal activation rapidly decreases after intravenous therapy with diuretics and vasodilators for class IV heart failure.* J Am Coll Cardiol 2002;39:1623-9.

51.Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline—Second Edition.* CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

52.Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples;Approved Guideline—Second Edition.* CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

53.Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.* CLSI/NCCLS document EP6-A [ISBN 1-56238-498-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.

54.Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline.* CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

55.Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline- Second Edition.* CLSI document EP17-A2 [ISBN 1-56238-795-2]. CLSI, 950 West Valley Road, Suite 2500, Wayne, PA 19087 USA, 2012.

※※56.Grimsey P, Frey N, Bendig G, et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. Int. J. Pharmacokinet. 2017;2(4):247-256.

※※57.Piketty ML, Prie D, Sedel F, et al. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. Clin Chem Lab Med. 2017;55(6):817-825

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