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1. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Sixth Edition*. CLSI Document GP41-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
2. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard-Sixth Edition*. CLSI Document GP39-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
3. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline-Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
4. Savukoski T, Engström E, Engblom J, Ristiniemi N, et al. Troponin-specific autoantibody interference in different cardiac troponin I assay configurations. *Clin Chem*. 2012;58(6): 1040-1048.
5. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition*. CLSI Document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
6. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline-Third Edition*. CLSI Document EP28-A3c. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
7. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.
8. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Clin Chem* 1988;34:261-264.

9. Mair J, Wagner I, Puschendorf B, et al. Cardiac troponin I to diagnose myocardial injury. *Lancet* 1993;341:838–839.
10. Adams JE, 3rd, Bodor GS, Davila-Roman VG, et al. Cardiac troponin I. A marker with high specificity for cardiac injury. *Circulation* 1993;88:101–106.
11. Corin SJ, Juhasz O, Zhu L, Conley P, Kedes L, Wade R. Structure and expression of the human slow twitch skeletal muscle troponin I gene. *J Biol Chem* 1994;269:10651–10659.
12. Mair J, Larue C, Mair P, Balogh D, Calzolari C, Puschendorf B. Use of cardiac troponin I to diagnose perioperative myocardial infarction in coronary artery bypass grafting. *Clin Chem* 1994;40:2066–2070.
13. Adams JE, 3rd, Sicard GA, Allen BT, et al. Diagnosis of perioperative myocardial infarction with measurement of cardiac troponin I. *N Engl J Med* 1994;330:670–674.
14. Adams JE, 3rd, Schechtman KB, Landt Y, Ladenson JH, Jaffe AS. Comparable detection of acute myocardial infarction by creatine kinase MB isoenzyme and cardiac troponin I. *Clin Chem* 1994;40:1291–1295.
15. Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in suspected cases of myocardial infarction. *Clin Chem* 1992;38:2203–2214.
16. Apple FS. Acute myocardial infarction and coronary reperfusion. Serum cardiac markers for the 1990s. *Am J Clin Pathol* 1992;97:217–226.
17. Thygesen K, Alpert JS, Jaffe AS. Joint ESC/ACCF/AHA/WHF Task Force for the Universal Definition of Myocardial Infarction. *J Am Coll Cardiol* 2012;60:1581–1598.
18. Wildi K, Twerenbold R, Mueller C. How acute changes in cardiac troponin concentrations help to handle the challenges posed by troponin elevations in non-ACS-patients. *Clin Biochem.* 48:2015;218–222.
19. IFCC Task Force on Clinical Applications of Cardiac Biomarkers (June 2014). Analytical characteristics of high-sensitivity cardiac troponin assays. *Clin Chem.* 2012;58:1. 54–61.
20. Amsterdam EA, Wenger NK, Brindis RG, Casey DE, et al. 2014 AHA/ACC Guideline for the Management of Patients With Non-ST-Elevation Acute Coronary Syndromes. *Circulation.* 2014;130:e344–e426. <https://doi.org/10.1161/CIR.0000000000000134>. Accessed 12/15/2016.
21. Roffi M, Patrono C, Collet JP, Mueller C, et al. 2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation. *European Heart Journal.* 2016;37: 267–315.
22. Antman E, Tanasijevic M, Thompson B, Schactman M, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *N Engl J Med.* 1996;Oct 31;335(18):1342–1349.
23. Chin C, Shah A, McAllister D, Cowell S, et al. High-sensitivity troponin I concentrations are a marker of an advanced hypertrophic response and adverse outcomes in patients with aortic stenosis. *Eur Heart J.* 2014 Sep 7;35:2312–2321.
24. Youssef A, Hassan A, El-Ghamry R, Ahmed A. Serum troponin-I as a prognostic marker in acute exacerbated chronic obstructive pulmonary disease patients. *Egypt J Chest Dis Tuberc.* 2013;62:549–555.
25. Hijazi Z, Siegbahn A, Andersson U, Granger C, et al. High-sensitivity troponin I for risk assessment in patients with atrial fibrillation: insights from the Apixaban for Reduction in Stroke and other Thromboembolic Events in Atrial Fibrillation (ARISTOTLE) trial. *Circulation.* 2014;Feb 11;129(6):625–34.
26. Omland T. Cardiac Troponins: a tool for a personalized medicine strategy in stable coronary artery disease? *J Am Coll Cardiol.* 2014;63 (4):355–357.
27. Avouac J, Meune C, Chenevier-Gobeaux C, Borderie D, et al. Cardiac biomarkers in systemic sclerosis: contribution of high-sensitivity cardiac troponin in addition to N-terminal pro-brain natriuretic peptide. *Arthritis Care & Res.* 2015;67(7):1022–1030.
28. White HD, Tonkin A, Simes J, Stewart R, et al. Association of Contemporary Sensitive Troponin I Levels at Baseline and Change at 1 Year With Long-Term Coronary Events Following Myocardial Infarction or Unstable Angina: Results from the LIPID Study (Long-Term Intervention With Pravastatin in Ischaemic Disease). *J Am Coll Cardiol.* 2014;63(4):345–354.
29. Yeon JL, Lee H, Park JS, Kim SJ, et al. Cardiac troponin I as a prognostic factor in critically ill pneumonia patients in the absence of acute coronary syndrome. *J Crit Care.* 2015;30(2):390–394.
30. Clinical and Laboratory Standards Institute. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.* CLSI Document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
31. Clinical and Laboratory Standards Institute. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Third Edition.* CLSI Document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
32. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.* CLSI Document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
33. Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR.* 1988. 37:377–382, 387–388.
34. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline-Fourth Edition.* CLSI Document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
35. Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.
36. Apple FS, Sandoval Y, Jaffe AS, Ordonez-Llanos J. IFCC Task Force on Clinical Applications of Cardiac Bio-Markers. Cardiac troponin assays: Guide to understanding analytical characteristics and their impact on clinical care. *Clin Chem.* 2017;63(1):73–81.
- **37. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

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