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- Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al., Luminescent oxygen channeling assay (LOCI[®]): sensitive, broadly applicable homogenous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
- Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc. Natl. Acad. Sci. USA, Vol 91, pp.5426-5430, June 1994 Biochemistry.
- Ullman, EF, Homogenous Immunoassays. The immunoassay Handbook, 2nd ed., David Wild editor, 2001;192-194.
- Clinical and Laboratory Standards Institute (CLSI). Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. CLSI document GP41-A6 (ISBN 1-56238-650-6). Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2007.
- Clinical and Laboratory Standards Institute (CLSI). Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. CLSI document GP39-A6 (ISBN 1-56238-740-5). Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2010.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 (ISBN 1- 56238-555-0). CLSI,940 West Valley Road, Suite1400, Wayne, PA. 19087-1898, 2004.
- Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. CLSI Document EP07-A2. Wayne, PA. Clinical and Laboratory Standards Institute; 2005.
- Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. CLSI Document EP28-A3c. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- Fred S. Apple, Yader Sandoval, Allan S. Jaffe, and Jordi Ordoñez-Llanos, for the IFCC Task Force on Clinical Applications of Cardiac Bio-MarkersCardiac Troponin Assays: Guide to Understanding Analytical Characteristics And Their Impact On Critical Care. Clinical Chemistry 63:1. 73-81;2017.
- Thygesen K, Alpert JS, Jaffe AS. Joint ESC/ACCF/AHA/WHF Task force for the Universal Definition of Myocardial Infarction. J AM Coll Cardiol 2012;60:1581-1598.
- Wildt K, Twerenbold R, Mueller C. How accurate are cardiac troponin concentrations help to handle the challenges posed by troponin elevations in non-ACS-patients. Clin Biochem. 48:2015:218-222.
- IFCC Task Force on Clinical Applications of Cardiac Biomarkers (June 2014) Analytical Characteristics of high-sensitivity cardiac troponin assays. Clin Chem 58:1:54-61;2012
- Amsterdam EA, Wenger NK, Brindis RG, Casey DE, et al. 2014 AHA/ACC Guideline for the Management of Patients with Non-ST_Elevation Acute Coronary Syndromes.Circulation. 2014;130:e344-e426.
- Roffi M, Patrono C, Collet JP, Mueller C, et al. 2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation. European Heart Journal. 2016;37:267-315.
- Antman E, Tanasejivic M, Thompson B, Schactman M, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med. 1996;Oct 31;335(18):1342–1349.
- Chin C, Shah A, McAllister D, Cowell S, et al. High-sensitivity troponin I concentrations are a marker of an advanced hypertrophic response and adverse outcomes in patients with aortic stenosis. Eur Heart J. 2014 Sep 7;35:2312–2321.
- Youssef A, Hassan A, El-Ghannry R, Ahmed A. Serum troponin-I as a prognostic marker in acute exacerbated chronic obstructive pulmonary disease patients. Egypt J Chest Dis Tuberc. 2013;62:549–555.
- Hijazi Z, Sieghain A, Andersson U, Granger C, et al. High-sensitivity troponin I for risk assessment in patients with atrial fibrillation: insights from the Apixan for Reduction in Stroke and other Thromboembolic Events in Atrial Fibrillation (ARISTOTLE) trial. Circulation 2014;Feb 11;129(6):625-34.
- Omland T. Cardiac Troponins: a tool for a personalized medicine strategy in stable coronary artery disease. J Am Coll Cardiol 2014;63(4):355–357.
- Avouac J, Meune C, Chenavier-Gobeaux C, Borderie D, et al. Cardiac biomarkers in systemic sclerosis: contribution of high-sensitivity cardiac troponin in addition to N-terminal pro-brain natriuretic peptide. Arthritis Care & Res. 2015;67(7):1022–1030.
- White HD, Tonkin A, Simes J, Stewart R, et al. Association of Contemporary Sensitive Troponin I Levels at Baseline and Change at 1 Year With Long-Term Coronary Events Following Myocardial Infarction or Unstable Angina: Results from the LIPID Study (Long-Term Intervention With Pravastatin in Ischaemic Disease). J Am Coll Cardiol 2014;63(4):345–354.
- Yeon JL, Lee H, Park JS, Kim SJ, et al. Cardiac troponin I as a prognostic factor in critically ill pneumonia patients in the absence of acute coronary syndrome. J Crit Care 2015;30(2):390–394.
- Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999; 45 (7): 942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creative Kinane MB by using F(ab')2 conjugate and polymouse IgG. Clin Chem 1992; 38:1737-1742.
- Savukoski T, Engström E, Engblom J, Ristiemi N, et al. Troponin-specific autoantibody interference in different cardiac troponin I assay configurations. Clin Chem 2015;61(6):1040-1048.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. Clin Chem 1988;34:261–264.
- Mair J, Wagner J, Puschendorf B, et al. Cardiac troponin I to diagnose myocardial injury. Lancet 1993;341:838-839.
- Adams JE, 3rd, Bodor GS, Davila-Roman VG, et al. Cardiac Troponin I. A marker with high specificity for cardiac injury. Circulation 1993;88:101-106.
- Corin SJ, Jubelt Q, Zhu L, Conley P, Kedes L, Wade R. Structure and expression of the human slow twitch skeletal muscle troponin I gene. J Biol Chem 1994;269:10651-10659.
- Mair J, Larue C, Mair P, Balogh D, Catzolari C, Puschendorf B. Use of Cardiac troponin I to diagnose perioperative myocardial infarction in coronary artery bypass grafting. Clin Chem 1994;40:2066-2270.
- Adams JE, 3rd, Sicard GA, Allen BT, et al. Diagnosis of perioperative myocardial infarction with measurement of cardiac troponin I. N Engl J Med 1994;330:670-674.
- Adams JE, 3rd,Schechtman KB, Landt Y, Ladenson JH, Jaffe AS. Comparable detection of acute myocardial infarction by creatine kinase MB Isoenzyme and cardiac troponin I. Clin Chem 1994;40:1291-1295.
- Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibodies for an assay of cardiac troponin I and preliminary results in suspected cases of myocardial infarction. Clin Chem 1992;38:2203-2214.
- Apple FS. Acute myocardial infarction and coronary reperfusion. Serum cardiac markers for the 1990s. AM J Clin Pathol 1992;97:217-226.
- Clinical and Laboratory Standards Institute. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Third Edition. CLSI Document EP05-A3.Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. CLSI Document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- C