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- CLSI. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. CLSI document GP-41-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
- CLSI. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. CLSI document GP39-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI document GP44-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
- Piketty ML, Prie D, Sedel F, *et al*. High-dose biotin therapy leading to false biochemical endocrine profiles:validation of a simple method to overcome biotin interference. *Clin Chem Lab Med*. 2017;55(6):817-825
- Bischoff-Ferrari HA, Giovannucci E, Willett WC, *et al*. Estimation of optimal serum concentrations of 25-hydroxyvitamin D for multiple health outcomes. *Am J Clin Nutr* 2006; 84:18-28.
- Smith DC, Johnson CS, *et.al*. A Phase I Trial of Calcitriol (1,25 dihydroxycholecalciferol) in patients

with Advanced Malignancy. *Clin Cancer Res* June 1999; 5; 1339.

8. Holick MF et al. Evaluation, Treatment, and Prevention of Vitamin D Deficiency: An Endocrine Society Clinical Practice Guideline, *J Clin Endocrinol Metab* 2011; 96 (7): 1911-1930.

9. Grimsey P, Frey N, Bendig G, *et al.* Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and *in vitro* immunoassay interference. *Int. J. Pharmacokinet.* 2017;2(4):247–256.

10. Holick MF. Vitamin D status: measurement, interpretation, and clinical application. *Ann Epidemiol.* 2009; 19(2): 73–78.

11. KDOQI Clinical Practice Guidelines for Bone Metabolism and Disease in Children With Chronic Kidney Disease. http://www.kidney.org/PROFESSIONALS/kdoqi/guidelines_pedbone/guide8.htm

12. Souberbielle JC, Body JJ, Lappe JM, et al. Vitamin D and musculoskeletal health, cardiovascular disease, autoimmunity and cancer: Recommendations for clinical practice. *Autoimmun Rev* 2010; 9: 709-715.

13. Dawson-Hughes B, Heaney RP, Holick MF. Estimates of optimal vitamin D status. *Osteoporos Int.* 2005; 16(7): 713-716.

14. Vieth R. Why the minimum desirable serum 25-hydroxyvitamin D level should be 75 nmol/L (30 ng/mL). *Best Pract Res Clin Endocrinol Metab.* 2011; 25(4): 681-91.

15. Holick MF. The role of vitamin D for bone health and fracture prevention. *Curr Osteoporosis Rep.* 2006; 4(3): 96–102.

16. Ozkan B, Hatun, Ş, Bereket A. Vitamin D intoxication. *Turk J Pediatr.* 2012; 54: 93–98.

17. CLSI. *Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline— Third Edition.* CLSI Document EP28-A3c Wayne, PA: Clinical and Laboratory Standards Institute; 2008.

18. Holick MF. Vitamin D Deficiency. *N Engl J Med.* 2007; 357: 266–81.

19. Holick MF. Vitamin D deficiency. *N Engl J Med* 2007 July 19; 357(3): 266-81.

20. Lai JK, Lucas RM, Clements MS, Harrison SL, Banks E. Assessing vitamin D status: pitfalls for the unwary. *Mol Nutr Food Res* 2010 August; 54(8): 1062-71.

21. Rollins G. Vitamin D testing—what's the right answer? Labs grapple with confusing analytics, evidence. *Clin Lab News* 2009 July; 35(7): 1-9.

22. Hart GR. Medical Conditions Associated with Vitamin D Deficiency and the Clinical Consequences. Review Series Volume 1, Sept 2004; 1–8.

23. CLSI. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline— Second Edition.* CLSI document EP05-A2. CLSI, Wayne, PA: Clinical and Laboratory Standards Institute; 2004.

24. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline— Third Edition.* CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

25. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition.* CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.

26. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.* CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.

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