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- 1) Sack, R.B., et al.: Enterotoxigenic *Escherichia coli* Isolated from Patients with Severe Cholera-like Disease, J. Infect. Dis., 123, 378(1971).
- 2) Smith, H.W. & Gyles, C.L.: The Relationship Between Two Apparently Different Enterotoxins Produced by Enteropathogenic Strains of *Escherichia coli* of Porcine Origin, J. Med. Microbiol., 3, 387(1970).
- 3) Gyles, C.L. & Barnum, D.A.: A Heat-Labile Enterotoxin from Strains of *Escherichia coli* Enteropathogenic for Pigs, J. Infect. Dis., 120, 419(1969).
- 4) Gyles, C.L.: Immunological Study of the Heat-Labile Enterotoxins of *Escherichia coli* and *Vibrio cholerae*, Infect. Immun., 9, 564(1974).
- 5) Clements, J.D. & Finkelstein, R.A.: Immunological Cross-Reactivity Between a Heat-Labile Enterotoxin(s) of *Escherichia coli* and Subunits of *Vibrio cholerae* Enterotoxin, Infect. Immun., 21, 1036(1978).
- 6) Clements, J.D. & Finkelstein, R.A.: Demonstration of Shared and Unique Immunological Determinants in Enterotoxins from *Vibrio cholerae* and *Escherichia coli*, Infect. Immun., 22, 709(1978).
- 7) Kudoh, Y., et al.: Detection of Heat-Labile Enterotoxin of *Escherichia coli* by Reversed Passive Hemagglutination Test with Specific Immunoglobulin Against Cholera Toxin, Proceedings of the 14th Joint Conf., U. S.-Japan Coop. Med. Sci. Program, Cholera Panel, Toho Univ., Tokyo, 266(1979).
- 8) Finkelstein, R.A., et al.: Pathogenesis of Experimental Cholera: Biologic Activities of Purified Procholeragen A, J. Immunol., 96, 440(1966).
- 9) Evans, D.J., et al.: Production of Vascular Permeability Factor by Enterotoxigenic *Escherichia coli* Isolated from Man, Infect. Immun., 8, 725(1973).
- 10) 日本細菌学会教育委員会編: 逆受身赤血球凝集反応(Reversed Passive Hemagglutination, RPHA)コレラ菌と毒素原性大腸菌の検査法, 菜根出版, 75(1981).

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