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1. Fitzgerald-Bocarsly P, Herberman R, Hercend T, et al. A definition of natural killer cells, In: Ades E, Lopez C, eds. Natural Killer Cells and Host Defense. Basel: Karger; 1989:1.
2. Jackson A. Basic phenotyping of lymphocytes: Selection and testing of reagents and interpretation of data. Clin Immunol Newslett. 1990;10:43-55.
3. Kidd P, Vogt R. Report of the Workshop on the evaluation of T-cell subsets during HIV infection and AIDS. Clin Immunol Immunopathol. 1989;52:3-9.
4. Landay A, Muirhead K. Procedural guidelines for performing immunophenotyping by flow cytometry. Clin Immunol Immunopathol. 1989;52:48-60.
5. Schmidt R. Monoclonal antibodies for diagnosis of immunodeficiencies. Blut. 1989;59:200-206.
6. Nicholson J. Use of flow cytometry in the evaluation and diagnosis of primary and secondary immunodeficiency diseases. Arch Pathol Lab Med. 1989;113:598-605

7. Foucar K, Gueken Jr. Clinical application of immunologic techniques to the diagnosis of lymphoproliferative and immunodeficiency disorders. *Lab Medicine*. 1982;13:403.
8. Cohen S, Weetman A. Activated interstitial and intraepithelial thyroid lymphocytes in autoimmune thyroid disease. *Acta Endocr*. 1988;119:161-166.
9. Smolen J, Chused T, Leiserson W, Reeves J, Alling D, Steinberg A. Heterogeneity of immunoregulatory T-cell subsets in systemic lupus erythematosus. *Am J Med*. 1982;72:783-790.
10. Giorgi J, Hultin L. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. *Clin Immunol Newslett*. 1990;10:55-62.
11. Antel J, Bania M, Noronha A, Neely S. Defective suppressor cell function mediated by T8+ cell lines from patients with progressive multiple sclerosis. *J Immunol*. 1986;137:3436-3439.
12. Gratama J, Naipal A, Oljans P, et al. T lymphocyte repopulation and differentiation after bone marrow transplantation: Early shifts in the ratio between T4+ and T8+ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. *Blood*. 1984;63(6):1416-1423.
13. Bishop G, Hall B, Duggin G, Horvath J, Sheil A, Tiller D. Immunopathology of renal allograft rejection analyzed with monoclonal antibodies to mononuclear cell markers. *Kidney Internat*. 1986;29:708-717.
14. Wolde-Mariam W, Peter J. Recent diagnostic advances in cellular immunology. *Diagnost Med*. 1984;7:25-32.
15. Landay A, Ohlsson-Wilhelm B, Giorgi J. Application of flow cytometry to the study of HIV infection. *AIDS*. 1990;4:479-497.
16. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: New and previously defined clusters. In: McMichael A, ed. *Leucocyte Typing III: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1986:788-803.
17. Beverley P. Production and use of monoclonal antibodies in transplantation immunology. In: Touraine J, Trager J, Betuel H, eds. *Transplantation and Clinical Immunology XI*. Amsterdam: Excerpta Medica; 1980:87-94.
18. Rowe D, Beverley P. Characterization of breast cancer infiltrates using monoclonal antibodies to human leucocyte antigens. *Br J Cancer*. 1984;49:149-159.
19. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman S, eds. *Leucocyte Typing*. Berlin: Springer-Verlag; 1984:25-60.
20. Dimitriu-Bona A, Burmester G, Waters S, Winchester R. Human mononuclear phagocyte differentiation antigens. I. Patterns of antigenic expression on the surface of human monocytes and macrophages defined by monoclonal antibodies. *J Immunol*. 1983;130(1):145-152.
21. Herrmann F, Komischke B, Odenwald E, Ludwig W. Use of monoclonal antibodies as a diagnostic tool in human leukemia. I. Acute myeloid leukemia and acute phase of chronic myeloid leukemia. *Blut*. 1983;47:157-163.
22. Schwinzer R. Cluster report: CD45/CD45R. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford University Press; 1989:628-634.
23. Goyert S, Tesio L, Ashman L, et al. Report on the CD14 Cluster Workshop. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:789-794.
24. Bernstein I, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leukocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. *Leukocyte Typing II: Human Myeloid and Hematopoietic Cells*. New York: Springer-Verlag; 1984;3:1-25.
25. Jayaram Y, Hogg N. Surface expression of CD14 molecules on human neutrophils. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford University Press; 1989:796-797.
26. Haynes B. Summary of T cell studies performed during the Second International Workshop and Conference on Human Leukocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York: Springer-Verlag; 1986:1:3-30.
27. Ledbetter J, Evans R, Lipinski M, Cunningham-Rundles C, Good R, Herzenberg L. Evolutionary conservation of surface molecules that distinguish T lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. *J Exp med*. 1981;153 (February):310-323.
28. Kan E, Wang C, Wang L, Evans R. Non-covalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with anti-Leu-4 antibody. *J Immunol*. 1983;131(2):536-539.
29. Knowles R. Immunochemical analysis of the T cell-specific antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York: Springer-Verlag; 1986:259-288.
30. Nadler L, B cell/leukemia panel workshop: Summary and comments. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. *Leukocyte Typing II: Human B Lymphocytes*. New York: Springer-Verlag; 1986:2:3-43.
31. van Dongen J, Krissansen G, Wolvers-Tettero I, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. *Blood* 1988;71(3):603-612.
32. Brenner M, Groh V, Porcelli F, Hochstenbach H, et al. Structure and distribution of the human gd T-cell receptor. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:1049-1053.
33. Duken B, Moller P, Pezzutto A, Schwartz-Albiez R, Molderhauer G. B cell antigens: CD19. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1986:34-36.
34. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. *Annual Rev Immunol*. 1988;(6):629.

35. Evans R, Wall D, Platsoucas C, et al. Thymus-dependent membrane antigens in man: Inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. *Proc Natl Acad Sci USA*. 1981;78(1):544-548.
36. Maddon P, Dalgleish A, McDougal J, Clapham P, Weiss R, Axel R. The T4 gene encodes the AIDS virus receptor and is expressed in the immune system and the brain. *Cell*. 1986;47:333-348.
37. Dalgleish A, Beverly P, Clapham P, Crawford D, Greaves M, Weiss R. The CD4 (T4) antigen is an essential component of the receptor for the AIDS virus. *Nature*; 1984;312(December):763- 767.
38. Kurrle R. Cluster report: CD3, In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:290-293.
39. Wood F, Warner N, Warnke R. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. *J Immunol*. 1983;131(1):212-216.
40. Lanier L, Le A, Phillips J, Warner N, Babcock G. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. *J Immunol*. 1983; 131(4):1789-1796.
41. Anderson P, Blue M-L, Morimoto C, Schlossman S. Cross-linking of T3 (CD3) with T4 (CD4) enhances the proliferation of resting T lymphocytes. *J Immunol*. 1987;139:679-682.
42. Eichmann K, Johnson J-I, Falk I, Emmrich F. Effective activation of resting mouse T lymphocytes by cross-linking submitogenic concentration of the T cell antigen receptor with either Lyt-2 or L3T4. *Eur J Immunol*. 1987;17:643-650.
43. Gallagher P, Fazekas de St. Groth B, Miller J. CD4 and CD8 molecules can physically associate with the same T-cell receptor. *Proc Natl Acad Sci USA*. 1989;86:10044-10048.
44. Rudd C, Burgess K, Barber E, Schlossman S. Monoclonal antibodies to the CD4 and CD8 antigens precipitate variable amounts of CD4/CD8-associated p56-lck activity. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:326-327.
45. Terry L, Disanto J, Small T, Flomenberg N. Differential expression of the CD8 and Lyt-3 antigens on a subset of human T-cell receptor g/d-bearing lymphocytes. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989; 345-346.
46. Moebius U. Cluster report: CD8. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:342-343.
47. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1 a monoclonal antibody blocking Fc receptor functions. I. Characterization of the lymphocyte subset reactive with B73.1. *J Immunol*. 1983;130(5):2133-2141.
48. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions, II. Studies of B73.1 antibody-antigen interaction on the lymphocyte membrane. *J Immunol*. 1983;130(5):2142-2148.
49. Perussia B, Trinchieri G, Jackson A, et al. The Fc receptor for IgG on human natural killer cells: Phenotypic, functional, and comparative studies with monoclonal antibodies, *J Immunol*. 1984;133(1):180-189.
50. Schmidt R. Non-lineage/natural killer section report: New and previously defined clusters. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:517-542.
51. Schubert J, Lanier L, Schmidt R. Cluster report: CD56. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:699-702.
52. Lanier L, Testi R, Bindl J, Phillips J. Identity of Leu-19 (CD56) leukocyte differentiation antigen and neural cell adhesion molecule. *J Exp Med*. 1989;169(June):2233-2238.
53. Lanier L, Chang C, Azuma M, Ruitenberg J, Hemperly J, Phillips J. Molecular and functional analysis of human natural killer cell-associated neural cell adhesion molecule (N-CAM/CD56). *J Immunol*. 1991;146:4421-4426.
54. Lanier L, Le A, Civin C, Loken M, Phillips J. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. *J Immunol*. 1986;136(12):4480-4486.
55. Phillips J, Lanier L. Dissection of the lymphokine-activated killer phenomenon: Relative contribution of peripheral blood natural killer cells and T lymphocytes to cytotoxicity. *J Exp Med*. 1986;164:814-825.
56. Giorgi J. Lymphocyte subset measurements: Significance in clinical medicine. In: Rose N, Friedman H, Fahey J, eds. *Manual of Clinical laboratory Immunology*. 3rd ed. Washington DC: American Society of Microbiology. 1986:236-246.
57. Angadi C. Lack of Leu-3a epitope on T-helper (CD4) lymphocytes. *J Clin Lab Anal*. 1990;4:193-195.
58. Prince H, Hirji K, Waldbeser L, Plaeger-Marshall S, Kleinman S, Lanier L. Influence of racial background on the distribution of T cell subsets and Leu-11-positive lymphocytes in healthy blood donors. *Diag Immunol*. 1985;3:33-37.

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