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Bp	設定標準化団体	VME	Very Major Error
EA	Essential Agreement	ME	Major Error
CA	Category Agreement	mE	minor Error

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= US Food and Drug Administration 510(k) cleared
CLSI = Clinical and Laboratory Standards Institute
CASFM = Comite de l'Antibiogramme de la Societe Francaise de Microbiologie
E = External performance data
I = Internal performance data
- = Not available
N/A = Not applicable（非該当）
Ref. = Reference method for clinical performance study.

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1. Barry, AL The Antimicrobial Susceptibility Test, Principles and Practices, Lea and Febiger, Philadelphia, PA. 1976.
2. Clinical Laboratory Standards Institute (CLSI), Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically, M7- A7, Wayne, Pennsylvania, January 2006.
3. Clinical and Laboratory Standards Institute (CLSI), Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement, M100-S18, Vol. 27, No. 1, January 2008.
4. Clinical and Laboratory Standards Institute (CLSI), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-third Informational Supplement, M100-S22, January 2012.
5. Clinical and Laboratory Standards Institute (CLSI), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fourth Informational Supplement; M100-S24, January 2014.
6. Clinical and Laboratory Standards Institute (CLSI), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Informational Supplement, M100-S25, January 2015.
7. Comité de l'Antibiogramme de la Société Française de Microbiologie. Communiqué 1996. Path Biol, 1996, 44, n° 8, I-VIII.
8. Comité de l'Antibiogramme de la Société Française de Microbiologie, Communiqué 2007.
9. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM), Recommendations 2012.
10. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2014.
11. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2015.
12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.
13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.
14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.
15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.
16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.
17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Tenover RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).
19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

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