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1. Balows A, Hausler WJ, Herrmann KL, Isenberg HD, Shadomy HJ, editors. Manual of Clinical Microbiology, 5th ed. American Society for Microbiology, Washington, D.C. 1991.
2. Balows A, Truper HG, Dworkin M, Harder W, and Schleifer K-H, editors. The Prokaryotes - a Handbook on the Biology of Bacteria: Exophysiology, Isolation, Identification, Applications, 2nd ed., Volume II. Springer-Verlag, New York 1992.
3. Brenner DJ, Krieg NR, Staley JT, Garrity GM, editors. Bergey’s Manual of Systematic Bacteriology, 2nd Edition. Springer, New York, NY. 2005.
4. Clinical and Laboratory Standards Institute, M50-A, Quality Control for Commercial Microbial Identification Systems; Approved Guideline, Vol. 28 No. 23.
5. Clinical Laboratory Improvement Amendments of 1988. 42 U.S.C 263a. PL 100-578. 1988.
6. Difco Manual Dehydrated Culture Media and Reagents for Microbiology. 10th ed. 1984.
7. Holmes B, Costas M, On SL, Vandamme P, Falsen E, Kersters K. *Neisseria weaveri* sp. nov. (formerly CDC group M-5), from dog bite wounds of humans. Int J Syst Bacteriol. 1993 Oct; 43(4):687-93.
8. Holt JG, Krieg NR, Sneath PHA, Staley JT, Williams ST, editors. Bergey’s Manual of Determinative Bacteriology, 9th ed. Williams and Wilkins, Baltimore 1994.
9. Kilian M, Nicolet J, Biberstein EL, Biochemical and Serological Characterization of *Haemophilus pleuropneumoniae* (Matthews and Pattison 1961) Shope 1964 and Proposal of a Neotype Strain. Int J Syst Bacteriol. 1978. 28:20-26.
10. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC, editors. Color Atlas and Textbook of Diagnostic Microbiology 4th ed. Lippincott, Philadelphia, PA 1992.
11. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC, editors. Color Atlas and Textbook of Diagnostic Microbiology 5th ed. Lippincott, Philadelphia, PA. 1997.
12. Koneman EW, Allen SD, Janda WM, Procop GW, Schreckenberger PC, Winn WC, Woods GL, editors. Color Atlas and Textbook of Diagnostic Microbiology 6th ed. Lippincott, Philadelphia, PA. 2006.

13. Krieg NR, Holt JG, editors. *Bergey's Manual of Systematic Bacteriology*, 9th ed. Williams and Wilkins, Baltimore 1984.
14. Lindqvist, K., A *Neisseria* Species Associated with Infectious Keratoconjunctivitis of Sheep: *Neisseria ovis* nov. spec. *The Journal of Infectious Diseases*. 1960. 106:162-165.
15. Manafi M, Kneifel W, Bascomb S. Fluorogenic and chromogenic substrates used in bacterial diagnostics. *Microbiol.* 1991; Rev. 55:335- 348.
16. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover RH, editors. *Manual of Clinical Microbiology*, 7th ed. American Society for Microbiology, Washington, D.C. 1999.
17. Murray PR, Baron EJ, Jorgensen JH, Pfaller MA and Tenover RH, editors. *Manual of Clinical Microbiology*, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue — Approved Guideline, 1997.
19. Nørskov-Lauritsen N, Kilian M. 2006. Reclassification of *Actinobacillus actinomycetemcomitans*, *Haemophilus aphrophilus*, *Haemophilus paraphrophilus*, and *Haemophilus segnis* as *Aggregatibacter actinomycetemcomitans* gen. nov., comb. nov., *Aggregatibacter aphrophilus* comb. nov., and *Aggregatibacter segnis* comb. nov., and emended description of *Aggregatibacter aphrophilus* to include V factor dependent and V factor-independent isolates. *IJSEM*. 2006. 56:2135- 2146.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.
21. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry and D.W. Warnock. 2011. *Manual of Clinical Microbiology*, 10th edition. American Society for Microbiology, Washington, D.C.
22. Weyant RS, Moss CW, Weaver RE, Hollis DG, Jordon JG, Cook EC, and Daneshvar MI. *Identification of Unusual Pathogenic and Gram-Negative Aerobic and Facultatively Anaerobic Bacteria* 2nd ed. Williams & Wilkins, Philadelphia, PA 1996.
- * 23. Tonjum, T. 2005. *Neisseria*. In *Bergey's Manual of Systematic Bacteriology*, 2nd ed., Volume 2 Part C. Garrity, G.M., Brenner, D.J., Krieg, N.R., Staley, J.T. (eds), Springer, New York, p. 778
- * 24. Vandamme, P, F.E. Dewhirst, B.J. Pasterand S.L.W. On. 2005. *Campylobacter*. In *Bergey's Manual of Systematic Bacteriology*, 2nd ed., Volume 2 Part C. Garrity, G.M., Brenner, D.J., Krieg, N.R., Staley, J.T. (eds), Springer, New York, p. 1151.
- * 25. Vandamme, P., P. Segers, M. Vancanneyt, K. Van Hove, R. Mouton, J. Hommez, F. Dewhirst, B. Paster, K. Kersters, E. Falsen, L.A. Devriese, M. Bisgaard, K.-H. Hinze and W. Mannheim. *Ornithobacterium rhinotracheale* gen. nov., sp. nov., Isolated from the Avian Respiratory Tract. *Int. J. Syst. Bacteriol.* 1994. 44:24-37.
- * 26. Dickerman, A., A.B. Bandara and T.J. Inzana. Phylogenomic analysis of *Haemophilus parasuis* and proposed reclassification to *Glaesserella parasuis*, gen. nov., comb. nov. *Int. J. Syst. Evol. Microbiol.* 2020. 70:180-186.

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