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1. Approved Guideline – Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests, GP44-A4. 2010. Clinical and Laboratory Standards Institute.
2. Cembrowski GS, Carey RN. Laboratory quality management: QC ↔ QA. ASCP Press, Chicago, IL, 1989.
3. Demers, LM, Spencer, CA. Laboratory medicine practice guidelines: laboratory support for the diagnosis and monitoring of thyroid disease. Clin Endocrinol (Oxf) 2003; 58(2): 138-140.
4. Lee, SL. When is the TSH normal? New criteria for diagnosis and management. 12th Annual Meeting of the American Association of Clinical Endocrinologists (AAEC), May, 2003.
5. Haugen, BR. When isn't TSH normal and why? Clinical implications and causes. 12th Annual Meeting of the American Association of Clinical Endocrinologists (AAEC), May, 2003.
6. Singer, PA. Now it's normal – now it's not: food for thought from a complex case. 12th Annual Meeting of the American Association of Clinical Endocrinologists (AAEC), May, 2003.
7. ベックマン・コールター株式会社 社内データ 386902B
8. ベックマン・コールター株式会社 社内データ B01902A
9. 市原清志・標準化 up-to-date 共有基準範囲設定国際プロジェクト・臨床化学 38: 416-423. 2009.
10. 市原清志・河口勝憲編: エビデンスに基づく検査診断実践マニュアル. 日本教育研究センター, 2011.
11. Ichihara K, Boyd J. An appraisal of statistical procedures used in derivation of reference intervals. Clin Chem Lab Med 2010; 48(11): 1537-1551.
12. Kricka, L. Interferences in immunoassays – still a threat. Clin Chem 2000; 46: 1037-1038
13. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
14. White, GH. Recent advances in routine thyroid function testing. CRC - Critical Reviews in Clinical Laboratory Sciences. 1987; 24: 315-362.
15. Liewendahl, K, Majuri, H, Helenius, T. Thyroid function tests in patients on long-term treatment with various anticonvulsant drugs. Clinical Endocrinology, 1978; 8: 187-191.
16. Wenzel, KW. Pharmacological interference with in vitro tests of thyroid function. Metabolism. 1981; 30(7): 717-732.
17. Wilke, TJ. Estimation of free thyroid hormone concentrations in the clinical laboratory. Clinical Chemistry. 1986; 32(4): 585-592.
18. Spencer, CA. Thyroid status: trends in testing – selective test use cuts cost. Clinical Chemistry News. November, 1989; 15(11): 9-14.

** 19. Piketty ML, et al. 2017. High-dose biotin therapy leading to false biochemical endocrine profiles: validation of a simple method to overcome biotin interference. Clin Chem Lab Med. 55(6): 817-825.

** 20. Grimsey P, et al. 2017. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. 2(4): 247-256.

21. Gornall, AG, Luxton, AW, Bhavnani, BR. Endocrine Disorders. In Applied biochemistry of clinical disorders. 1986; 305-318. Philadelphia, PA: J. B. Lippincott Co.
22. Watts, NB, Keffer, JH. The thyroid gland. In Practical Endocrine Diagnosis. 1982; 77-96. Philadelphia, PA: Lea & Febiger.
23. Tentative Guideline - Internal quality control: Principles and definitions, C24-T. 1987. National Committee for Clinical Laboratory Standards.
24. Krouwer, JS, Rabinowitz, R. How to improve estimates of imprecision. Clinical Chemistry. 1984; 30: 290-292.
25. DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.

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