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- 1. C A and Ashwood, Tietz Textbook of Clinical Chemistry, 2nd Edition. 1994. Edited by Burtis, E R Philadelphia, PA: W B Saunders Co., 1857-1863.
- 2. Yen, S S C , The Human Menstrual Cycle- Neuroendocrine regulation. In Reproductive Endocrinology. Edited by Yen, S S C and Jaffe, R B Philadelphia, PA: W B Saunders Co., 273-308, 1991.
- 3. Carr, B R , Disorders of the ovary and female reproductive tract. In Williams Textbook of Endocrinology, 8th edition. Edited by Wilson, J D and Foster, D W Philadelphia, PA: W B Saunders Co., 733-798, 1992.
- 4. Hall, J E Polycystic ovarian disease as a neuroendocrine disorder of the female reproductive axis. In Endocrinology and Metabolism Clinics of North America, Neuroendocrinology II. Edited by Veldhuis, J D Philadelphia, PA: W B Saunders Co., 75-92.
- 5. Reyes-Fuentes, A. and Veldhuis, J D 1993. Neuroendocrine physiology of the normal male gonadal axis. In Endocrinology and Metabolism Clinics of North America. Edited by Velduis, J D Philadelphia, PA: W.B. Saunders Co., 22:93-124.
- 6. Manual Guide – Safety Management, No. CDC-22, Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta GA: Centers for Disease Control.
- 7. Approved Standard – Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture – H3-A4. 1998. National Committee for Clinical Laboratory Standards, 4th edition.
- 8. Approved Guideline – Procedures for the Handling and Processing of Blood Specimens, H18-A2. 1999. National Committee for Clinical Standards.
- 9. Cembrowski GS, Carey RN. Laboratory Quality Management: QC & QA. ASCP Press, Chicago, IL, 1989.
- 10. Kricka, L. Interferences in Immunoassays – Still a Threat. Clin Chem 2000; 46: 1037.
- 11. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613-621.
- 12. Williams, Textbook of Endocrinology, 8th edition. Edited by Wilson, J D and Foster, D W Philadelphia, PA: W B Saunders Co., 978 - 995, 1992.
- 13. Tentative Guideline – User evaluation of precision performance of clinical chemistry devices, EP5-T. 1984. National Committee for Clinical Laboratory Standards, 4, N8.
- 14. Krouer, J S, Rabinowitz, R. 1984. How to improve estimates of imprecision. Clinical Chemistry, 30: 290-292.
- 15. HHS Publication No 93-8395, 3rd ed., May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: U.S. Government Printing Office.
- 16. ベックマン・コールター株式会社 社内データ

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