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1. Barry, A. L. 1970. The Antimicrobial Susceptibility Test, Principles and Practices. p. 26. Lea and Febiger, Philadelphia, PA.
2. Chitwood, L. A. 1969. Tube dilution anti-microbial susceptibility testing: efficacy of a microtechnique applicable to diagnostic laboratories. Appl. Microbiol. 17:707.
3. Jorgensen, J.H., M.A. Pfaller, K.C. Carroll, G. Funke, M.L. Landry, S.S. Richter, D.W. Warnock (eds), 2015. Manual of Clinical Microbiology, 11th ed. American Society for Microbiology, Washington D. C.
4. Ericsson, H. M., and J. C. Sherris. 1971. Antibiotic sensitivity testing. Report of an international collaborative study. Acta Pathol. Microbiol. Scand. Sect. B Suppl. 217: 1-90.
5. Gerlach, E. H. 1974. Microdilution I: A comparative study, p. 63. In A. Balows (ed.), Current Techniques for Antibiotic Susceptibility Testing. C. C. Thomas, Springfield, IL.
6. Harwick, J. H., P. Weiss, and F. Fekety, Jr. 1968. Application of microtitration techniques to bacteriostatic and bacteriocidal antibiotic susceptibility testing. J. Lab. Clin. Med. 72:511.
7. Jorgensen, J.H. 1994. Detection of antimicrobial resistance in *Streptococcus pneumoniae* by use of standardized susceptibility testing methods and recently developed interpretive criteria. Clin. Microbiol. Newsl. 16: 97-101.
8. Jorgensen, J.H. and G.V. Doern. 1995. Practical guidelines for in vitro susceptibility testing of selected gram-positive bacteria. Clin. Microbiol. Newsl. 17: 81-85.
9. Jorgensen, J.H., J.M. Swenson, F.C. Tenover, M.J. Ferraro, J.A. Hindler and P.R. Murray. 1994. Development of interpretive criteria and quality control limits for broth microdilution and disk diffusion antimicrobial susceptibility testing of *Streptococcus pneumoniae*. J. Clin. Microbiol. 32: 2448-2459.
10. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically. 2015. Approved Standard M07-A10. Clinical and Laboratory Standards Institute, Wayne, PA.
11. Performance Standards for Antimicrobial Susceptibility Testing. 2016. CLSI Document M100-S26. Clinical and Laboratory Standards Institute, Wayne, PA.
12. Tilton, R. C., and L. Newberg. 1974. Standardization of the microdilution susceptibility test, p. 77. In A. Balows (ed.), Current techniques for antimicrobial susceptibility testing. C. C. Thomas, Springfield, IL.
13. Washington, J. A., E. Warren, and A. G. Karlson. 1973. Stability of barium sulfate turbidity standards. Appl. Microbiol. 24:1013.
14. World Health Organization 1961. Standardization of methods for conducting microbic sensitivity tests. Second Report of the Expert Committee on Antibiotics. Technical Report Series #210, p. 1. Geneva.

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Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

© 2021 Beckman Coulter, Inc.
All rights reserved.

Eng. Ref.: C56630