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- Bernard A, Boumsell L Dausset J, Milstein C and Schlossman SF, eds:1984.Leukocyte Typing. New York:Springer-Verlag. p. 26, 28, 41-42, 196, 201, 215-216.
- McMichael AJ. ed:1987.Leukocyte Typing III .Oxford: Oxford University Press,p.38-40, 42, 106-107, 116, 199, 302-306, 308, 315, 475.
- Reinherz EL, Haynes BF, Nadler LM and Bernstein ID:1986. Leukocyte Typing II New York:Springer-Verlag. Vol.1.p.6, 16, 364-365.
- Reinherz EL,Haynes BF Nadler LM and Bernstein ID : 1986.Leukocyte Typing II New York:Springer-Verlag.Vol.2.p. 8, 15-25, 37.
- Reinherz EL and Schlossman SF:1980. The differentiation and function of human T lymphocytes.Cell: 821-827,
- Aiuta F, Cerottini J-C Coombs RRA, Cooper M, Dickler HB Froland S, Fudenberg HH, Greaves MF, Grey HM, Kunkel HG, Natvig J, Preud'homme J-L, Rabellino E, Ritts RE, Rowe DS, Seligmann M, Siegal FP, Stjernsward J, Terry WD and Wybran J:1975. Identification, enumeration and isolation of B and T lymphocytes from human peripheral blood, international Union of Immunoblgical Societies (IUIS) Report-July 1974.Clin Immunol and Immunopathol 3:584-597.
- Nadler LM, Anderson KC Marti. G, Bates M, Park E, Daley JF and Schlossman SF : 1983 . B4 , a human B lymphocyte-associated antigen expressed on normal mitogen-activated, and malignant B lymphocytes. J Immunol 131 : 244-250.
- Reinherz EL, Haynes BF, Nadler LM and Bernstein ID. 1986. Leukocyte Typing II . New York:Springer-Verlag.Vol.2, p.8. 20-21, 189
- Stashenko P. Nadler LM, Hardy R and Schlossman SF. 1980. Characterization of a human B lymphocyte-specific antigen. J Immunol 125 : 1678-1680,
- Nadler LM, Stashenko P, Hardy R van Agthoven A, Terhorst C and Schlossman SF : 1981. Characterization of a human B cell-specific antigen (B2) distinct from B1. J Immunol 126 : 1941.
- Benjamin, E and Leskowitz S:1991. Disorders of the Immune Response. In Immunology:A Short Course. Second Edition. New York:Wiley-Liss.p211-244.
- Reinherz EL, O'Brien C, Rosenthal P and Schlossman SF: 1980. The cellular basis for viral-induced immunodeficiency: Analysis by monoclonal antibodies. J Immunol 125: 1269-1274.
- Felsenstein D, Carney WP, Iacoviello VR and Hirsch MS:1985. Phenotypic properties of atypical lymphocytes in cytomegalovirus-induced mononucleosis. J Infect Dis 152 : 198-203.
- Rinaldo CR, Ho M, Hamoudi WH, Gui X and DeBiasio RL: 1983. Lymphocyte subsets and natural killer cell responses during cytomegalovirus mononucleosis. Infect Immun 40 : 472-477.
- Goldstein G, Litter J and Mittler R:1982. Immunoregulatory changes in human disease detected by monoclonal antibodies to T lymphocytes, In : Monoclonal Antibodies in Clinical Medicine. McMichael AJ and Fabre JW,eds. New York, NY: Academic Press, p. 39-70.
- Schmidt Re: 1989. Monoclonal antibodies for diagnosis of immunodeficiencies. Blut 59: 200-206.
- de Martini RM and Parker JW:1989 Immunoblgic alterations in human immunodeficiency virus infection:A review. J Clin Lab

Anal 3:56-70.

18. Fauci AS: 1988. The human deficiency virus: Infectivity and mechanism of pathogenesis. Science 239:617-622.
19. Taylor MGJ, Fahey JL, Detets R and Giorgi JV: 1989. CD4 Percentage, CD4 number and CD4:CD8ratio in HIV infection. How to choose and how to use. JAIDS 2:114-124.
20. Pedrazzini A, Freedman AS, Andersen J, Heflin L, Anderson K, Takvorian T, Canellos GP, Whitman J, Coral F, Ritz J and Nadler LM: 1989. Anti-B cell monoclonal antibody purged autologous bone marrow transplantation for B-cell non-Hodgkin's lymphoma. Phenotypic reconstitution and B-cell function: Blood 74: 2203-2211.
21. Ramos EL, Turka LA, Leggat JE, Wood IG, Milford EL and Carpenter CB: 1989. Decrease in phenotypically defined T helper inducer cells (T4+4B4+) and increase in T suppressor effector cells (T8+2H4+) in stable renal allograft recipients. Transplantation 47: 465-471.
22. National Comminee for Clinical Laboratory Standards. 1991. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture. 3rd Edition (H3-A3). National Committee for Clinical Laboratory Standards, Villanova, PA.
23. Loken MR, Brosnan JM, Bach BA and Ault KA: 1990. Establishing optimal lymphocyte gates for immunophenotyping by flow cytometry. Cytometry 11:453-459.
24. Koepke JA and Landay AL: 1989. Precision and accuracy of absolute lymphocyte counts. Clin Immunol Immunopathol 52: 19-27.

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