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2/3

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1. Wu AH, Apple FS, Gibler WB, Jesse RL, Warshaw MM, Valdes R Jr. National Academy of Clinical Biochemistry Standards of Laboratory Practice: recommendations for the use of cardiac markers in coronary artery disease. Clin Chem 1999; 45: 1104-1121.
2. Approved Guideline-Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests, H18-A4. 2010. Clinical and Laboratory Standards Institute.
3. Troponin: what laboratorians should know to manage elevated results. <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/TipsandArticle/sonDeviceSafety/ucm109362.htm>
4. Approved Standard-Sixth Edition, Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture, H3-A6. 2007. Clinical and Laboratory Standards Institute.
5. Cembrowski GS, Carey RN. Laboratory quality management: QC $\approx$ QA. ASCP Press, Chicago, IL, 1989.
6. ベックマン・コールター株式会社 社内データ
7. Kricka L. Interferences in immunoassays-still a threat. Clin Chem 2000; 46: 1037-1038.
8. Bjerner J, Nustad K, Norum LF, Olsen KH, Børmer OP. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
9. Lum G, Solarz D, Farney L. False positive cardiac troponin results in patients without acute myocardial infarction. Labmedicine 2006; 37(9): 546-550.
10. Lingwood D, Ballantyne JS. Alkaline phosphatase-immunoglobulin conjugate binds to lipids in vitro, independent of antibody selectivity. J Immunol Methods 2006; 311: 174-177.
11. Mitchell JRA. Back to the future: so what will fibrinolytic therapy offer your patients with myocardial infarction? Br Med J 1986; 292: 973-978.
12. Thygesen K, Alpert JS, Jaffe AS. Joint ESC/ACCF/AHA/WHF Task Force for the Universal Definition of Myocardial Infarction. Third universal definition of myocardial infarction. Eur Heart J 2012; 33: 2551-67. J Am Coll Cardiol 2012; 60: 1581-98. Circulation 2012; 126: 2020-35.
13. Perry SV. The regulation of contractile activity in muscle. Biochem Soc Trans 1979; 7: 593-617.
14. Vallins WJ, Brand NJ, Dabhade N, Butler-Browne G, Yacoub MH, Barton PJ. Molecular cloning of human cardiac troponin I using polymerase chain reaction. FEBS Lett 1990; 270: 57-61.
15. Perry SV. Troponin I: inhibitor or facilitator. Mol Cell Biochem 1999; 190: 9-32.
16. Katrukha A. Antibody selection strategies in cardiac troponin assay. In: Wu AHB, editor. Cardiac markers, 2nd ed. Totowa (NJ): Humana Press Inc., 2003, 173-185.
17. Larue C, Defacque-Lacquement H, Calzolari C, Le Nguyen D, Pau B. New monoclonal antibodies as probes for human cardiac troponin I: epitopic analysis with synthetic peptides. Mol Immunol 1992; 29: 271-278.
18. Mair J, Morandell D, Genser N, Lechleitner P, Dienstl F, Puschendorf B. Equivalent early sensitivities of myoglobin, creatine kinase MB mass, creatine kinase isoform ratios, and cardiac troponins I and T for acute myocardial infarction. Clin Chem 1995; 41: 1266-1272.
19. Mair J, Genser N, Morandell D, Maier J, Mair P, Lechleitner P, Calzolari C, Larue C, Ambach E, Dienstl F, Pau B, Puschendorf B. Cardiac troponin I in the diagnosis of myocardial injury and infarction. Clin Chim Acta. 1996; 245: 19-38.
20. Jaffe AS, Babuin L, Apple FS. Biomarkers in acute cardiac disease: the present and the future. J Am Coll Cardiol 2006; 48: 1-11.
21. Babuin L, Jaffe AS. Troponin: the biomarker of choice for the detection of cardiac injury. CMAJ 2005; 173: 1191-1202.
22. Jeremias A, Gibson CM. Narrative Review: alternative causes for elevated cardiac troponin levels when acute coronary syndromes are excluded. Ann Intern Med 2005; 142: 786-91.
23. Thygesen K, Alpert JS, White HD; Joint ESC/ACC/AHA/WHF Task Force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. Eur Heart J 2007; 28: 2525-38.
24. Myocardial infarction redefined-a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the redefinition of myocardial infarction. Eur. Heart J 2000; 21: 1502-1513. J Am Coll Cardiol 2000; 36: 959-969.
25. Hamm CW, Giannitsis E, Katus HA. Cardiac troponin elevations in patients without acute coronary syndrome. Circulation 2002; 106: 2871-2872.
26. FM, Decker WW, Diercks DB, Ghaemmaghami CA, Nazarian D, Brady WJ, Hahn S, Jagoda AS; American College of Emergency Physicians Clinical Policies Subcommittee (Writing Committee) on Non-ST-Segment Elevation Acute Coronary Syndromes. Clinical policy: critical issues in the evaluation and management of adult patients with non-ST-segment elevation acute coronary syndromes. Ann Emerg Med. 2006; 48: 270-301.
27. Jaffe AS, Apple FS; High Sensitivity Cardiac Troponin: Hype, Help, and Reality. Clin Chem 2010; 56:342-344.
28. Agewall S, Giannitsis E, Jernberg T, Katus H. Troponin elevation in coronary vs non-coronary disease. Eur Heart J. 2011; 32: 404-411.
29. Kelley WE, Januzzi JL, Christenson RH. Increases of cardiac troponin in conditions other than acute coronary syndrome and heart failure. Clin Chem 2009; 55: 2098-2122.DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13-Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.
30. DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13-Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.

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