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1. Cancer Facts and Figures - 1997. American Cancer Society Inc., 1997; 1-32.
2. Wang MC, Valenzuela LA, Murphy GP, Chu TM. Purification of human prostate specific antigen. Invest Urol 1979;17: 159-163.
3. McCormack RT, Rittenhouse HG, Finlay JA, Sokoloff RL, Wang TJ, Wolfert RL, Lilja H, Oesterling JE. Molecular forms of prostate-specific antigen and the human kallikrein gene family: A new era. Urology 1995; 45: 729-744.
4. Christensson A, Laurell CB, Lilja H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine proteinase inhibitors. Eur J Biochem 1990; 194: 755-763.
5. Lilja H, Christensson A, Dahlen V, Matikainen MT, Nilsson O, Pettersson K, Lovgren T. Prostate-specific antigen in human serum occurs predominately in complex with alpha 1 antichymotrypsin. Clin Chem 1991; 37: 1618-1625.
6. Stenman UH, Leinonen J, Alfthan H, Rannikko S, Tuhkanen K, Alfthan O. A complex between prostate-specific antigen and alpha 1 antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostatic cancer: Assay of the complex improves clinical sensitivity for cancer. Cancer Res 1991; 51: 222-226.
7. Christensson A, Bjork T, Nilsson O, Dahlen U, Matikainen MT, Cockett ATK, Abrahamsson PA, Lilja H. Serum prostate antigen complexed to alpha 1 antichymotrypsin as an indicator of prostate cancer. J Urol 1993; 150: 100-105.
8. Catalona WJ, Smith DS, Wolfert RL, Wang TJ, Rittenhouse HG, Ratliff TL, Nadler RB. Evaluation of percentage of free serum prostate-specific antigen to improve specificity of prostate cancer screening. JAMA 1995; 274: 1214-1220.
9. Partin AW, Catalona WJ, Southwick PC, Subong ENP, Gasior GH, Chan DW. Analysis of percent free prostate-specific antigen (PSA) for prostate cancer detection: Influence of total PSA, prostate volume, and age. Urology 1996; 48 (Suppl): 55-61.
10. Van Cangh PJ, De Nayer P, Sauvage P, Tombal B, Elsen M, Lorge F, Opsomer R, Wese F. Free to total prostate-specific antigen (PSA) ratio is superior to total PSA in differentiating benign prostate hypertrophy from prostate cancer. Prostate 1996; 7(Suppl): 30-34.
11. Woodrum DL, Brawer MK, Partin AW, Catalona WJ, Southwick PC. Interpretation of free prostate specific antigen clinical research studies for the detection of prostate cancer. J Urol 1998; 159: 5-12.
12. Catalona WJ, Partin AW, Slawin KM, Brawer MK, Flanigan RC, Patel A, Richie JP, deKernion JB, Walsh PC, Scardino PT, Lange PH, Subong ENP, Parson RE, Gasior GH, Loveland KG, Southwick PC. Use of the percentage of free prostate specific antigen to enhance differentiation of prostate cancer from benign prostatic disease: A prospective multicenter clinical trial. JAMA 1998; 279: 1542-1547.
13. Manual Guide - Safety Management, No. CDC-22, Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta GA: Centers for Disease Control.
14. Ornstein DK, Rao GS, Smith DS, Ratliff TL, Basler JW, Catalona WJ. Effect of digital rectal examination and needle biopsy on serum total and percentage of free prostate specific antigen levels. J Urol 1997; 157:195-198.
15. Yuan JJ, Coplen DE, Petros JA, Figenshau RS, Ratliff TL, Smith DS, Catalona WJ. Effects of rectal examination, prostatic massage, ultrasonography and needle biopsy on serum prostate specific antigen levels. J Urol 1992; 147: 810-814.

16. Woodrum DL, French C, Shamel LB. Stability of free PSA in serum samples under a variety of sample collection and sample storage conditions. Urology 1996; 48 (Suppl): 33–39.
17. Woodrum DL, French CF, Hill TM, Roman SJ, Slatore HL, Shaffer JL, York LG, Eure KL, Loveland KG, Gasior GH, Southwick PC, Shamel LB. Analytical performance of the Tandem®-R free PSA immunoassay measuring free prostate specific antigen. Clin Chem 1997; 43: 1203–1208.
18. Woodrum DL, York L. Two year stability of free and total PSA in frozen serum samples. Urology 1998; 52: 247–251.
19. Cembrowski GS, Carey RN. Laboratory Quality Management: QC & QA. ASCP Press, Chicago, IL, 1989.
20. Kricka, L. Interferences in Immunoassays – Still a Threat. Clin Chem 2000; 46: 1037.
21. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613–621.
22. Catalona WJ, Richie JP, Ahmann FR, Hudson MA, Scardino PT, Flanigan RC, deKernion JB, Ratliff TL, Kavoussi LR, Dalkin BL, Waters WB, MacFarland MT, Southwick PC. Comparison of digital rectal examination and serum prostate specific antigen in the early detection of prostate cancer. Results of a multicenter clinical trial of 6,630 men. J Urol 1994; 151: 1283–1290.
23. Keetch DW, Catalona WJ, Smith DS. Serial prostatic biopsies in men with persistently elevated serum prostate specific antigen values. J Urol 1994; 151: 1571–1574.
24. HHS Publication No 93–8395, 3rd ed., May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: U.S. Government Printing Office.
25. ベックマン・コールター株式会社 社内データ

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