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1. Whitley, RJ. Endocrinology. In Tietz textbook of clinical chemistry, 1660-1694. Edited by Burtis, CA, Ashwood, ER. WB Saunders, Philadelphia, PA,1994.
2. Ashby, CD. Prolactin. In Methods in clinical chemistry, 258-265. Edited by Pesce, AJ, Kaplan, LA. CV Mosby, St. Louis, MO, 1987.
3. Jeffcoate, SL. Assays for prolactin: guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem, 1986. 23:638-651.
4. Franks, S. Prolactin. In Hormones in blood, 280-331. Edited by Gray, CH, James, VHT. Academic Press, London, UK, 1979.
5. Frohman, LA. Diseases of the anterior pituitary. In Endocrinology and metabolism, 247-337. Edited by Felig, P. McGraw Hill, New York, NY, 1987.
6. Reichlin, S. Neuroendocrinology. In Textbook of endocrinology, 589-645. Edited by Williams RH. WB Saunders, Philadelphia, PA, 1981.
7. Kohler, PO. Diseases of the hypothalamus and anterior pituitary. In Harrison's Principles of internal medicine, 587-604. Edited by Pertersdorf, RG. McGraw Hill, New York, NY, 1983.
8. Howanitz, JH. Evaluation of endocrine function. In Clinical diagnosis and management by laboratory methods, 308-348. Edited by Henry, JB. WB Saunders, Philadelphia, PA, 1991.
9. Frantz, AG. The breasts. In Textbook of endocrinology, 400-411. Edited by Williams RH. WB Saunders, Philadelphia, PA, 1981.
10. Bardin, CW. The testes. In Textbook of endocrinology, 293-354. Edited by Williams RH. WB Saunders, Philadelphia, PA, 1981.
11. Liwicz BH. The hypothalamopituitary system. In Clinical chemistry theory, analysis, and correlation, 739-747. Edited by Kaplan, LA , Pesce, A J. CV Mosby, St. Louis, MO, 1984.
12. Kubasik, NP. Prolactin: its functional effects and measurement. Laboratory management 1987. 25(8):49-51.
13. Jacobs, LS. Prolactin. In Methods of hormone radioimmunoassay, 199-222. Edited by Jaffe, BM, Behrman HR : Academic Press, New York, NY, 1979.
14. Approved Guideline – Procedures for the handling and processing of blood specimens, H18-A2. 1999. National Committee for Clinical Laboratory Standards.
15. Cembrowski GS, Carey RN. Laboratory quality management: QC & QA. ASCP Press, Chicago, IL, 1989.
16. Kricka, L. Interferences in immunoassays – still a threat. Clin Chem 2000; 46: 1037.
17. Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613–621.
18. Tentative Guideline - User evaluation of precision performance of clinical chemistry devices, EP5-T. National Committee for Clinical Laboratory Standards, 4(8) 1984.
19. Krouwer, J.S., Rabinowitz, R. 1984. How to improve estimates of imprecision. Clinical Chemistry 30:290-292.T
20. Manual guide – Safety management, No. CDC-22, Decontamination of laboratory sink drains to remove azide salts. April 30, 1976. Atlanta GA: Centers for Disease Control.
21. ベックマン・コールター株式会社 社内データ

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