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- Approved Guideline – Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests, GP44-A4.2010. Clinical and Laboratory Standards Institute.
- Kricka L. Interferences in immunoassays – still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- Tsung SH. Creatine kinase isoenzyme patterns in human tissue obtained at surgery. Clin Chem 1976; 22: 173-174.
- Urdal P, Urdal K, Stromme JH. Cytoplasmic creatine kinase isoenzymes quantitated in tissue specimens obtained at surgery. Clin Chem 1983; 29: 310-313.
- Gouille JP, Mechard D, Laine G, Jeanmet A, Cramer J, Maitrot B, Fondimare A, Gruchy D, Letac B. Répartition isozymique de la créatine kinase dans différents organes humains. Intérêt en pathologie humaine. Ann Biol Clin 1979; 37: 303-307.
- El Allaf M, Chapelle JP, El Allaf D, Adam A, Faymonville ME, Laurent P, Heugheum C. Differentiating muscle damage from myocardial injury by means of the serum creatine kinase (CK) isoenzyme MB mass measurement/total CK activity ratio. Clin Chem 1986; 32: 291-295.
- Apple F, Preese L. 1 Creatine Kinase-MB: Detection of myocardial infarction and monitoring reperfusion. Journal of Clinical Immunoassay 1994; 17(1): 24-29.
- Adam JE, Abendschein DR, Jaffe AS. Biochemical markers of myocardial injury. Circulation 1993; 88(2): 750-763.
- Galen RS, Reiffel JA, Gambino SR. Diagnosis of acute myocardial infarction: relative efficiency of serum enzyme and isoenzyme measurements. J Am Med Assoc 1975; 232: 145-147.
- Mercer DW, Varat MA. Detection of cardiac specific creatine kinase isoenzyme in sera with normal or slightly increased total creatine-kinase activity. Clin Chem 1975; 21: 1088-1092.
- Strom S, Bendz R. Serum creatine kinase (CK) activity after M-subunit inhibition in patients with atypical CK isoenzymes – a comparison to findings in acute myocardial infarction. Clin Chem Acta 1983; 132: 73-81.
- Apple F, Preese L, Bennet R, Fredrickson A. Clinical and analytical evaluation of two immunoassays for direct measurement of creatine kinase MB with monoclonal anti-CK-MB antibodies. Clin Chem 1988; 4: 2364-2367.
- Schwartz JG, Brown RW, McMahan CA, Gage CL, Herber SA. Clinical and analytical evaluation of different methods for measurement of creatine-kinase isoenzyme MB. Clin Chem 1989; 35: 130-134.

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