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1. Croce MV, et al. Detection of Circulating Mammary Mucin (MUC1) and MUC1 Immune Complexes (MUC1-CIC) in Healthy Women. The International Journal of Biological Markers, 2001, Vol. 16 no. 2, 112–120.
2. Duffy MJ. CA 15-3 and related Mucins as Circulating Markers in Breast Cancer, Ann. Clin. Biochem, 1999, 36, 579–586.
3. von Mensdorff-Pouilly S, et al. Human MUC1 Mucin: A Multifaceted Glycoprotein. The International Journal of Biological Markers. 2000, Vol. 15 No. 4, 343–356.
4. Price MR, Rye PD, Finn OJ, Hilgers J. Summary Report on the ISOBM TD-4 Workshop: Analysis of 56 monoclonal antibodies against the MUC1 Mucin. Tumor Biol. 1998, 19 (suppl 1):1–20.
5. Clinton et al. A comparative study of four serological tumor markers for the detection of breast cancer. ISA Paper #2003-070 0067-8856/2003.
6. Colomer et al. Circulating CA 15-3 levels in the postsurgical follow-up of breast cancer patients and in non-malignant diseases. Breast Cancer Research and Treatment, 1989, Vol 13, 123–133.
7. American Cancer Society, Facts and Figures 2003.
8. Safi et al. The value of the tumor marker CA 15-3 in diagnosing and monitoring breast cancer. Cancer, Aug 1991, 574–582.
9. HHS Publication, 4th ed., April 1999. Biosafety in Microbiological and Biomedical Laboratories.
Available <http://www.cdc.gov/od/ohs/biosfty/bmbl4/bmbl4toc.htm>
10. DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.
11. Approved Guideline – Procedures for the Handling and Processing of Blood Specimens, H18-A3. 2004. Clinical and Laboratory Standards Institute.
12. Cembrowski GS, Carey RN. Laboratory Quality Management: QC⇌QA. ASCP Press, Chicago, IL, 1989.
13. Kricka, L. Interferences in Immunoassays – Still a Threat. Clin Chem 2000; 46: 1037-1038.
14. Bjerner J, et al. Immunometric Assay Interference: Incidence and Prevention. Clin Chem 2002; 48: 613–621.
15. Bon GG, et al. Maternal serum CA125 and CA15-3 antigen levels in normal and pathological pregnancy. Fetal Diagn Ther 2001; 16: 166-172.
16. Bombardieri E, et al. CA 15.3 Determination in patients with breast cancer: clinical utility for the detection of distant metastases. Eur J Cancer 1993; Vol 29A, No. 1: 144-146.
17. Krouwer JS, Rabinowitz R. How to improve estimates of imprecision. Clin Chem 1984; 30: 290-292.
18. Tentative Guideline – User evaluation of precision performance of clinical chemistry devices, EP5-T. 1984. National Committee for Clinical Laboratory Standards, 4: No. 8.
19. ベックマン・コールター株式会社 社内データ

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