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- Approved Guideline - Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests, GP44-A4. 2010. Clinical and Laboratory Standards Institute.
- Approved Guideline - Seventh Edition, Collection of Diagnostic Venous Blood Specimens, GP41. 2017. Clinical and Laboratory Standards Institute.
- Cembrowski GS, Carey RN. Laboratory quality management: QC ≡ QA. ASCP Press, Chicago, IL, 1989.
- Approved Guideline - Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory, EP28-A3c. October 2010. Clinical and Laboratory Standards Institute.
- Meisner M, Procalcitonin — Biochemistry and Clinical Diagnosis, ISBN 978-3-8374-1241-3, UNI-MED, Bremen 2010.
- Meisner M, Procalcitonin — Biochemistry and Clinical Diagnosis, ISBN 978-3-8374-1241-3, UNI-MED, Bremen 2010.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- Lingwood D, Ballantyne JS. Alkaline phosphatase-immunoglobulin conjugate binds to lipids in vitro, independent of antibody selectivity. Journal of Immunological Methods 2006; 311: 174-177.
- Kishnani PS, Rush ET, et al. Monitoring guidance for patients with hypophosphatasia treated with asfotase alfa. Molecular Genetics and Metabolism 2017; 122: 4-17.
- Meisner M. Update on Procalcitonin Measurements. Ann Lab Med 2014; 34:263-273.
- Approved Guideline - Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures, EP17-A2. June 2012. Clinical and Laboratory Standards Institute.
- Wacker C, Prkno A, Brunkhorst FM, Schlattmann P (2013) Procalcitonin as a diagnostic marker for sepsis: a systematic review and meta-analysis. Lancet Infect Dis. 13(5):426-435.
- Rhodes A, et al. Surviving Sepsis Campaign: International Guidelines for Management of Sepsis and Septic Shock:2016. Society of Critical Care Medicine and the European Society of Intensive Care Medicine. March 2017; 45, Number 3.
- Singer M, Deutschman CS, Seymour CW, et al: The Third International Consensus Definitions for Sepsis and Septic Shock (Sepsis-3). JAMA 2016; 315:801-810.
- Fleischmann C, Scherag A, Adhikari NK, et al. Assessment of Global Incidence and Mortality of Hospital-treated Sepsis. Current Estimates and Limitations. Am J Respir Crit Care Med 2016; 193(3): 259-72.
- Maruna P, Nedelnikova K, Gurlich R: Physiology and Genetics of Procalcitonin. Physiol Res 2000; 49(Suppl 1):S57-S61.
- Müller B, White JC, Nylén ES, et al. Ubiquitous Expression of the Calcitonin-Receptor-Like Receptor 1 Gene in Multiple Tissues in Response to Sepsis. J Clin Endocrinol & Metab 2001; 86(1):396-404.
- Brunkhorst, et al. Kinetics of procalcitonin iatrogenic sepsis. Intensive Care Med (1998) 24:888-892.
- Schuetz P, Bretscher C, Bernasconi L, and Mueller, B. Overview of procalcitonin assays and procalcitonin-guided protocols for the management of patients with infections and sepsis. EXPERT REVIEW OF MOLECULAR DIAGNOSTICS, 2017. VOL. 17, NO. 6, 593-601.
- Clementine YF, Yap, Tar Choon Aw. The Use of Procalcitonin in Clinical Practice. Proceedings of Singapore Healthcare, Volume 23, Number 1, 2014.

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