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- 1) Boyanton BL, Blick KE. Stability studies of Twenty-Four Analytes in Human Plasma and Serum. Clin Chem 2002;48/12:2242-2247.
- 2) O'Keane MP, Cunningham SK. Evaluation of three different specimen types (serum, plasma lithium heparin and serum gel separator) for analysis of certain analytes: clinical significance of differences in results and efficiency in use. Clin Chem Lab Med 2006;44(5):662-668.
- 3) Guder WG, da Fonseca-Wollheim F, Heil W, et al. Quality of Diagnostic Samples, in brochure: Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine, 3rd Ed. 2010.
- 4) Elfath D, Cooney J, McDaniel R, et al. Effect of frozen storage of serum on the level of 22 chemistry analytes. Clin Chem 1991;37:931.
- 5) Gambino SR, Schreiber H. The measurement of CO2 content with the autoanalyzer. A comparison with 3 standard methods and a description of a new method (alkalinization) for preventing loss of CO2 from open cups. Am J Clin Path 1966;45:406.
- 6) Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 7) Breuer, J. Report on the Symposium, Drug Effects in Clinical Chemistry Methods. Eur J Clin Chem Clin Biochem. 1996, 34, p.385～386.
- 8) Sonntag, O. et al. Drug interferences in clinical chemistry, recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem. 2001, 38, p.376～385.
- 9) Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243
- 10) Scott, MG. et al. Electrolytes and blood gases, in Tietz NW. Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA. WB Saunders Co. 1999, p.1,065～1,066.
- 11) Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.

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COBAS is a trademark of Roche.