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1. Braunstein GD, Rasor J *et al.* Serum human chorionic gonadotropin levels throughout normal pregnancy. Am J Obstet Gynecol 1976;126: 678-81
2. Brody S, Carlström G. Immunoassay of human chorionic gonadotropin in normal and pathologic pregnancy. J Clin Endocrinol Metab 1962; 22: 564-74
3. Cole LA, Kardana A, *et al.* The deactivation of hCG by nicking and disassociation. J Clin Endocrinol Metab 1993; 76: 704-13
4. Pastorfi de GB, *et al.* Serum chorionic gonadotropin activity after molar pregnancy, therapeutic abortion, and term delivery. Am J Obstet Gynecol 1974; 118: 293-94
5. Reyes FI, Winter JSD, *et al.* (Postpartum disappearance of chorionic gonadotropin from the maternal and neonatal circulations. Am J Obstet Gynecol 1985; 153: 486-99
6. Cole LA. Immunoassay of human chorionic gonadotropin, its free subunits, and metabolites. Clin Chem 1997; 43: 12: 223-43
7. Lövgren T, Pettersson K. Time-resolved fl uoroimmunoassay, advantages and limitations. In: Van Dyke K, Van Dyke R, eds. Luminescence immunoassay and molecular applications. Boca Raton: CRC Press, 1990: 233-53.
8. Lövgren T *et al.* One-step all-in-one dry reagent immunoassays with fl uorescent europium chelate label and time-resolved fl uorometry. Clin Chem 1996; 42, 1196-1201.
9. Burtis CA *et al.* Clinical Chemistry of Pregnancy. In: Tietz textbook of clinical chemistry and molecular diagnostics, fourth edition. St. Louis: Elsevier Saunders, 2006
10. Cammona F, Balasch J *et al.* Early hormonal markers of pregnancy outcome after in vitro fertilization and embryo transfer. Journal of Assisted Reproduction and Genetics 2003; 20, 12: 521-26.
11. Lambers MJ, van Weering HGI *et al.* Optimizing hCG cut-off values: A single determination on day 14 or 15 is sufficient for a reliable prediction of pregnancy outcome. European Journal of Obstetrics Gynecology and Reproductive Biology 2006; 127, 1: 94-8
12. Sugantha SE, Webster S, *et al.* Predictive value of plasma human chorionic gonadotropin following assisted conception treatment. Human Reproduction 2000; 15, 2: 469-73
13. Lode PV, Rainaho J, Petterson K. Quantitative, widerange, 5-minute point-of-care immunoassay for total hCG in whole blood. Clinical Chemistry 2004; 50, 6: 1025-35.
14. Azzazy HME, Romero LF, Hall L *et al.* Two-center clinical evaluation of a new automated fl uorometric immunoassay for the quantitative analysis of total betahuman chorionic gonadotropin. Clinical Biochemistry 2003; 36: 523-28
15. Bauer R *et al.* Stratus CS Stat Fluorometric Analyzer. Point-of-Care 2002; 1, 1: 54-7.

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