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1. Filatov VL et al. Troponin: Structure, properties, and mechanism of functioning. Biochemistry (Moscow) 1999; 64(9):969-85.
2. Farah CS, Reinach FC. The troponin complex and regulation of muscle contraction. FASEB J 1995; 9:755-67.
3. David C Gaze et al. Multiple molecular forms of circulating cardiac troponin: analytical and clinical significance, Ann Clin Biochem 2008; 45: 349 -355.
4. Wu AHB. Cardiac markers. In: Wild D et al. The immunoassay handbook. London: Nature Publishing Group, 2001: 623-34.
5. Collinson PO, Boa FG, Gaze DC. Measurement of cardiac troponins. Ann Clin Biochem 2001; 38(5): 423-49.
6. Apple FS, et al. Improved detection of minor ischemic myocardial injury with measurement of serum cardiac troponin I. Clin Chem 1997; 43: 2047-51.
7. Manuel Del Rey JM, Madrid AH et al. Cardiac troponin I and minor cardiac damage: biochemical markers in a clinical model of myocardial lesions. Clin Chem 1998; 44: 2270-76.
8. Wu AHB, Apple FS, Gibler WB, Jesse RL et al. National Academy of Clinical Biochemistry Standards of Laboratory Practice:Recommendations for the use of cardiac markers in coronary artery diseases. Clin Chem 1999; 45:1104-21.
9. Thygesen K, Alpert JS et al. Myocardial infarction redefined – a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the redefinition of myocardial infarction. JACC 2000; 36(3):959-69.
10. Lövgren T, Pettersson K. Time-resolved fluoroimmunoassay, advantages and limitations. Van Dyke K, Van Dyke R, eds. In Luminescence Immunoassay and Molecular Applications. Boca Raton, USA: CRC Press, 1990; 233-53.
11. Lövgren T, Meriö L, et al. One-step all-in-one dry reagent immunoassays with fluorescent europium chelate label and time-resolved fluorometry. Clin Chem 1996; 42: 1196-201.
12. Katus HA, et al. Enzyme linked immunoassay of cardiac troponin T for the detection of acute myocardial infarction in patients. Mol Cell Cardiol 1989;21(7):1349-1353.
13. Katus HA, et al. Proteins of the troponin complex. Laboratory Medicine 1992;23(5):311-317.
14. Hamm CW, Ravkilde J, Gerhardt W et al. The prognostic value of serum troponin T in unstable angina. N Engl J Med 1992;327(3):146-150.
15. Ohmann EM, et al. Risk stratification with admission cardiac troponin T levels in acute myocardial ischemia. N Engl J Med 1996;335:1333-1334.
16. Christenson RH, Duh SH, Newby LK, Ohman EM et al. Cardiac troponin T and cardiac troponin I: relative values in short-term risk stratification of patients with acute coronary syndromes. Clin Chem 1998;44(3):494-501.
17. Lindahl B, et al. Mechanisms behind the prognostic value of Troponin T in unstable coronary artery disease: a FRISC II substudy. J Am Coll Card 2001;38:979-986.
18. Aviles RJ et al. Troponin T levels in patients with acute coronary syndromes, with or without renal dysfunction. N Engl J Med 2002; 346: 2047-2052

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