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- (1) Biosafety in Microbiological and Biomedical Laboratories, U.S. Department of Health and Human Services, Washington, 5th Edition, 2009.
- (2) CLSI Document H3-A6. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture. Vol. 27, No. 26, ISBN 1-56238-650-0, 2007.

(3) CLSI Document H21-A5. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays. Vol. 28, No. 5, ISBN 1-56238-657-3, 2008.

(4) Hanke AA, Roberg K, Monaca E, Sellmann T, Weber CF, Rahe-Meyer N, Görlinger K.

Impact of platelet count on results obtained from multiple electrode platelet aggregometry (Multiplate). Eur J Med Res. 2010 May 18;15(5):214-9.

(5) Stissing T, Dridi NP, Ostrowski SR, Bochsén L, Johansson PI. The influence of low platelet count on whole blood aggregometry assessed by Multiplate. Clin Appl Thromb Hemost. 2011 Nov-Dec;17(6):E211-7.

(6) Cardinal DC, Flower RJ. The electronic aggregometer: a novel device for assessing platelet behavior in blood. J Pharmacol Methods. 1980 Feb;3(2):135-58.

(7) Tóth O, Calatzis A, Penz S, Losonczy H, Siess W. Multiple electrode aggregometry: a new device to measure platelet aggregation in whole blood. Thromb Haemost. 2006 Dec;96(6):781-8.

(8) Görlinger K, Jambor C, Hanke AA, Dirkmann D, Adamzik M, Hartmann M, Rahe-Meyer N. Perioperative coagulation management and control of platelet transfusion by point-of-care platelet function analysis. Transfus Med Hemother. 2007; 34: 396-411.

(9) ROTEM platelet operating Manual Version 1.0.0.01.EN and higher.

(10) Uchiyama S. Clopidogrel resistance: identifying and overcoming a barrier to effective antiplatelet treatment. Cardiovasc Ther. Epub 2011 Jul 31.

(11) Woo KS, Kim BR, Kim JE, Goh RY, Yu LH, Kim MH, Han JY. Determination of the prevalence of aspirin and clopidogrel resistances in patients with coronary artery disease by using various platelet-function tests. Korean J Lab Med. 2010 Oct;30(5):460-8.

(12) Hummel T, Meves SH, Rüdiger K, Mügge A, Mumme A, Burkert B, Mühlberger D, Neubauer H. Prevalence of acetylsalicylic acid (ASA) - low response in vascular surgery. Chirurg. 2016 May;87(5):446-54.

(13) Chapman MP, Moore EE, Moore HB, Gonzalez E, Morton AP, Silliman CC, Sauaia A, Banerjee A.

Early TRAP pathway platelet inhibition predicts coagulopathic hemorrhage in trauma. Shock. 2015 Jun; 43 (6 Suppl 1):33

(14) Ortmann E, Klein AA, Sharples LD, Walsh R, Jenkins DP, Luddington RJ, Besser MW. Point-of-care assessment of hypothermia and protamine-induced platelet function with multiple electrode aggregometry (Multiplate) in patients undergoing cardiopulmonary bypass. Anesth Analg. 2013 Mar;116(3):533-40

(15) Scharf RE. Drugs that affect platelet function. Semin Thromb Hemost 2012; 38:865-883

(16) Flöck A, Zobel A, Bauriedel G, Tuleta I, Hammerstingl C, Höfels S, Schuhmacher A, Maier W, Nickenig G, Skowasch D. Antiplatelet effects of antidepressant treatment: a randomized comparison between escitalopram and nortriptyline. Thromb Res. 2010 Aug;126(2):E83-7.

(17) Velik-Salchner C, Maier S, Innerhofer P, Streif W, Klingler A, Kolbitsch C, Fries D. Point-of-care whole blood impedance aggregometry versus classical light transmission aggregometry for detecting acetylsalicylic acid and clopidogrel: the results of a pilot study. Anesth Analg 2008 Dec;107(6):1798-806.

(18) Solomon C, Trautinger S, Ziegler B, Hanke A, Rahe-Meyer N, Voelckel W, Schöchl H. Platelet function following trauma. A multiple electrode aggregometry study. Thromb Haemost. 2011 Aug;106(2):322-30.

(19) Rahe-Meyer N, Winterhalter M, Boden A, Froemke C, Piepenbrock S, Calatzis A, Solomon C. Platelet concentrates transfusion in cardiac surgery and platelet function assessment by multiple electrode aggregometry Acta Anaesthesiol Scand. 2009 Feb;53(2):168-75.

(20) Adamzik M, Görlinger K, Peters J, Hartmann M. Whole blood impedance aggregometry as a biomarker for the diagnosis and prognosis of severe sepsis. Crit Care. 2012 Oct 22;16(5):R204

(21) Petricevic M, Konosic S, Biocina B, Dirkmann D, White A, Mihaljevic MZ, Ivancan V, Konosic L, Svetina L, Görlinger K. Bleeding risk assessment in patients undergoing elective cardiac surgery using ROTEM and Multiplate impedance aggregometry. Anaesthesia. 2016 Jun;71(6):636-47

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