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Nose Spray

General Information

Refer to the label on the product or package before use.

Indications for Use

Nose Spay is an aerosol spray with a manual pump that enables users to administer a medicine by spraying into the ear, nose, and throat.

Contraindications

- Do not use if any damage or abnormality is observed before use or during use.
- Do not make any modification to this product.
- Do not reuse.

Instructions for Use

Preparation

1. Remove the cap.
2. Twist gently to detach the pump and the container.
Note: Pump and container are fixed with screw type.
3. Fill necessary quantity of medicine into the container under the instruction of pharmacists.
4. Twist the pump gently to attach to the container.
5. Confirm that the pump and the container are firmly fixed and then the cap is attached and firmly fixed.

Instructions for Use

1. Remove the cap.
2. Press the pump 2 or 3 times until medicine is sprayed in a mist.
3. Spray into the nasal cavity while inhaling slowly following prescribed usage and dosage.

Precautions

- Do not use if the package or the product is damaged.
- Do not make any modification to this product.
- Do not apply excessive force.
- Press the pump strongly. Spray may not be misty when pressing slowly.
- This instruction is only suggestion. The final decision on the technique and procedure used is the responsibility of the user.

Storage Instructions and Expiration Date

Before use: Store in a clean and dry condition away from sunlight, heat and humidity.

After use: Do not reuse and dispose the product after use following prescribed usage and dosage.

Symbols on the Label

	Do not reuse
	Batch code
	Re-order number