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HOYA INJECTOR SYSTEM (Handpiece & Cartridge)

Instructions for Use

PRESCRIPTION DEVICES

Caution: This product is for sale by or on the order of a physician, only.

DESCRIPTION

HOYA INJECTOR SYSTEM (hereinafter "system") is specially designed for reliably implanting an intraocular lens (IOL) in the capsular bag to correct aphakia following cataract surgery. The system consists of HOYA INJECTOR SYSTEM Cartridge (hereinafter "cartridge"; Fig.1) and HOYA INJECTOR SYSTEM Handpiece (hereinafter "handpiece"; Fig.2), both of which are used in combination. The polypropylene cartridge is supplied sterile (sterilized by ethylene oxide), and only for single use. The titanium alloy handpiece is supplied non-sterile and reusable, and must be sterilized prior to each use.

Fig.1 Cartridge (Type-E1)


Fig.2 Handpiece (ISH001)

INDICATIONS

The system is intended for use as an accessory for implanting a foldable soft acrylic IOL into the capsular bag. It is used only for qualified soft acrylic IOLs manufactured by HOYA CORPORATION, and must not be used for any other lens. Please refer to the "QUALIFIED MODEL LIST" affixed to the cartridge's outer box, for correct combination of IOL and handpiece models etc.

WARNINGS

1. The combination of the system is strictly limited. Never attempt to use any other combinations or systems from other manufactures.
2. Carefully check the "QUALIFIED MODEL LIST" affixed to the cartridge's outer box. Never attempt to use IOLs other than those applicable lenses.
3. Do not use the cartridge if its sterility seems to be compromised due to rupture in the sterilization pouch, etc.
4. Do not reuse or resterilize the cartridge. Reuse or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury or illness. Reuse or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to, the transmission of infectious diseases(s) from one patient to another. Contamination of the device may lead to injury or illness of the patient.
5. Do not use the cartridge after the indicated expiration date.

PRECAUTIONS

1. Before use, please consult this "INSTRUCTIONS FOR USE," and that attached to the handpiece. Use the devices only after training of surgical staff has been completed.
2. The handpiece must be cleaned and sterilized prior to each use.
3. Before use, confirm that there is no deformation, damage, foreign bodies or other abnormalities in the handpiece or cartridge. Never use the handpiece or cartridge if any such abnormalities are detected.
4. The cartridge is for single use only; after use, discard properly as medical waste. Handle the cartridge aseptically after opening the pouch. If implantation surgery is discontinued or suspended after opening the pouch, the IOL and cartridge must be discarded as medical waste.
5. Do not fold an IOL in methods other than using the provided lens case.
6. Do not use the IOL if the IOL is removed or dislodged from the lens case prior to folding.
7. Use ophthalmic viscosurgical device (OVD) (e.g. sodium hyaluronate for ophthalmic use) before loading the IOL into the cartridge. The OVD should be allowed to attain room temperature prior to use (for details, consult instructions for use attached to the OVD).
8. Do not attempt to turn the plunger to the reverse direction until the IOL is released from the cartridge.
9. Do not proceed if strong stiff resistance is felt, when implanting the IOL.
10. Do not attempt to disassemble or modify the handpiece or cartridge.
11. Store the cartridge and handpiece in a dry place. Do not store the cartridge in direct sunlight or at temperatures over 25°C.

CLEANING AND STERILIZATION OF HANDPIECE

Cleaning of handpiece

1. Remove and discard the used cartridge before cleaning the handpiece.
2. The handpiece must be cleaned by appropriate method just after use.
3. Flush the handpiece with sufficient water to rinse off blood, tissue fluid, cellular debris, intraocular irrigation solution, OVD, etc.
4. Rub and wash with a mild detergent and a soft brush.
5. Ultrasonic cleaning in water is also recommended.
6. After rinsing off the detergent under running water, rinse the handpiece thoroughly with distilled or deionized water.
7. After checking that the handpiece has been thoroughly cleaned and rinsed, dry it.

Notes

- Do not use hydrogen peroxide, bleaching agents, polishing agents, saline, etc. for cleaning of the handpiece.
- Use only a soft brush for rubbing and washing. Do not use hard brushes such as those made of metal, etc.

Sterilization of handpiece


1. The handpiece must be sterilized after clearing and drying.
2. Confirm that the rod does not stick out of the main body of the handpiece.
3. Put the handpiece in the sterilization pouch and seal it.
4. Sterilize the handpiece by steam (moist heat) according to one of the following general methods:
 - Gravity displacement cycle at 134 °C for 15 minutes
 - Prevacuum cycle at 134 °C for 4 minutes
 or
 - Sterilization at 121 °C for 30 minutes.

In all case, the above recommendations do not replace the need for each institution to identify the appropriate sterilization parameters and conduct validations to assure the sterility of this device.


5. Store the sterilized handpiece in clean and hygienic condition.

INSTRUCTIONS FOR USE


1. Clean and sterilize the handpiece in accordance with the "CLEANING AND STERILIZATION OF HANDPIECE."
2. Make sure that the handpiece works normally: push the plunger forward until the threads contact the main body, engage the threads and rotate the plunger clockwise, and then confirm that the plunger moves smoothly toward the end.
3. Confirm that there are no cracks, dirt, foreign bodies, etc., at the tip of the rod.
4. Pull back the plunger completely, so that the cartridge can be set.
5. Take the cartridge out of the sterilization pouch under aseptic condition, and infuse the OVD (e.g. sodium hyaluronate for ophthalmic use) into the cartridge (Fig.3). Use sufficient amounts of the OVD so that inside of cartridge from tapered part to the end can be filled with the OVD (Fig.4).


6. Fold the IOL by pushing the lever on the upper end of the provided lens case (Fig.5), and pick up with forceps (Fig.6).


7. Hold the cartridge with the lens diagram facing upward, and then load the folded IOL into the cartridge. Note: The direction of the intraocular lens of Fig.7 is inverted as compared to Fig.6.


8. Push the lens forward with leading haptic extending so that the lens may be placed into position as described in Fig.8.

9. Set the cartridge with the lens diagram facing upward on the handpiece (Fig.9).


10. Slowly push the plunger forward until the tip of rod touches the back edge of optic. At this time, the rod passes in front of the trailing haptic and pushes the edge of the optic (Fig.10). Continue pushing the plunger forward until the threads contact the main body of the handpiece. Rotate the plunger clockwise to engage the threads. Then rotate the plunger about one-quarter turn.

11. Carefully observe to make sure that the tip of the rod pushes the back edge of the optic and does not run over the optic. If the rod runs over the optic, discard the IOL and the cartridge as medical waste.
12. Insert the tip of the cartridge into the eye through incision, keeping the slit of the cartridge in a downward position.
13. Hold the main body of the handpiece in one hand, rotate the plunger clockwise in the other hand, and start releasing the lens into the eye. N.B.: The leading haptic should be extending to the left as seen from surgeon's side (Fig.11).

14. After the leading haptic goes beneath the capsulorhexis, rotate the main body of the handpiece clockwise to release the IOL (Fig.12) as monitoring the lens motion recovering from its folded state, while carefully controlling the handpiece to prevent capsular or corneal endothelium damage by the tip of the cartridge, rod or haptic. The rotation direction is also indicated by an arrow on the main body (Fig.13).

15. Continuously rotate the plunger clockwise to completely release the optical portion of the lens from the cartridge. Then, pull out the cartridge from the incision. N.B.: The trailing haptic should be remained outside of the eye through the incision.
16. As carefully observing that the lens has recovered from its folded state in the eye without turning over, adjust the position of the lens to place the trailing haptic in the capsular bag, using a hook or other devices (Fig.14).

SYMBOLS

Symbols	Description
STERILE EO	Sterilized using ethylene oxide
⌘	Do not reuse
⌘	Do not resterilize
📖	Consult instructions for use
🔥	Upper limit of temperature
🏭	Manufacturer
EC REP	Authorized representative in the European Community
🕒	Used by (YYYY-MM : year-month))
SN	Serial number
LOT	Lot number

Manufacturer

HOYA CORPORATION
6-10-1 Nishi-Shinjuku Shinjuku-ku Tokyo 160-0023 Japan

Made by

HOYA INJECTOR SYSTEM Handpiece Nichiyu Electronics Co., Ltd., Japan
HOYA INJECTOR SYSTEM Cartridge HOYA MEDICAL SINGAPORE PTE., LTD., Singapore

Authorized Representative in the European Community

HOYA SURGICAL OPTICS GMBH
De-Saint-Exupéry-Strasse 8, 60549 Frankfurt am Main, Germany
Tel: +49-(0)69-664-268-0