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1. Bore, P. et al. 1987. Pharmacokinetics of Methotrexate and 7-Hydroxy-Methotrexate After Methotrexate Infusions. Cancer Drug Delivery 4:177 – 183.
2. Collier, C. P. et al. 1982. Analysis of methotrexate and 7-hydroxymethotrexate by high-performance liquid chromatography and preliminary clinical studies. Ther Drug Monit 4:371 – 380.
3. Erttmann, R. et al. 1985. 7-Hydroxy-Methotrexate and Clinical Toxicity Following High-Dose Methotrexate Therapy. J Cancer Res Clin Oncol 109:86 – 88.
4. Jacobs, S. A. et al. 1976. 7-Hydroxymethotrexate as a urinary metabolite in human subjects and rhesus monkeys receiving high dose methotrexate. J Clin Invest 57:534 – 538.
5. Belz, S. et al. 1994. High-performance liquid chromatographic determination of methotrexate, 7-hydroxymethotrexate, 5-methyltetrahydrofolic acid and folinic acid in serum and cerebrospinal fluid. J Chromatogr B Biomed Appl 661:109 – 118.
6. Breithaupt, H. and Kuenzlen, E. 1982. Pharmacokinetics of methotrexate and 7-hydroxy-methotrexate following infusions of high dose methotrexate. Cancer Treat Rep 66:1733 – 1741.
7. Wolfrom, C. et al. 1990. Pharmacokinetic study of methotrexate, folinic acid and their serum metabolites in children treated with high-dose methotrexate and leucovorin rescue. Eur J Clin Pharmacol 39:377 – 383.
8. Al-Turkmani, M. R. et al., 2010. Difficulty Measuring Methotrexate in a Patient with High-Dose Methotrexate–Induced Nephrotoxicity. Clin Chem 56:1792 – 1796.
9. Schwarz, S. et al. 2007. Glucarpidase (Carboxypeptidase G2) intervention in adult and elderly cancer patients with renal dysfunction and delayed methotrexate elimination after high-dose methotrexate therapy. The Oncologist 12:1299 – 1308.
10. Chabner, B. A. et al. 1972. Enzymatic cleavage of methotrexate provides a method for prevention of drug toxicity. Nature 239:395 – 397.
11. Widemann, B. C. et al. 1995. Carboxypeptidase-G2 rescue in a patient with high dose methotrexate-induced nephrotoxicity. Cancer 76:521 – 526.
12. Buchen, S. et al. 2005. Carboxypeptidase G2 rescue in patients with methotrexate intoxication and renal failure. British Journal of Cancer 92:480 – 487.
13. Prescribing information. 2008. Methotrexate Injection, USP. Hospira, Inc. Lake Forest, IL.
14. Bleyer, W. A. 1978. The clinical pharmacology of methotrexate: New applications of an old drug. Cancer 41:36 – 51.
15. Jaffe, N. and Gorlick, R. 2008. High-Dose Methotrexate in Osteosarcoma: Let the Questions Surcease—Time for Final Acceptance. J Clin Oncol 26:4365 – 4366.
16. Colom, H. et al. 2009. Population Pharmacokinetics of High-Dose Methotrexate After Intravenous Administration in Pediatric Patients With Osteosarcoma. Ther Drug Monit 31:76 – 85.
17. Leucovorin (Fusilev) Prescribing Information. 2008. Spectrum Pharmaceuticals, Inc. Irvine, CA.
18. Dombrowsky, E. et al. 2011. Evaluating performance of a decision support system to improve methotrexate pharmacotherapy in children and young adults with cancer. Ther Drug Monit 33:99 – 107.
19. Widemann, B. C. and Adamson, P. C. 2006. Understanding and managing methotrexate nephrotoxicity. Oncologist 11:694 – 703.
20. Blum, R. et al. 2002. Significant impairment of high-dose methotrexate clearance following vancomycin administration in the absence of overt renal impairment. Annals of Oncology 13:327 – 330.
21. Martelli, N. et al. 2011. Methotrexate pharmacokinetics in childhood acute lymphoblastic leukaemia: a prognostic value? J Clin Pharm Ther 36:237 – 245.
22. Mazanec, D. J. and Grisanti, J. M. 1989. Drug-induced osteoporosis. Cleve Clin J of Med 56:297 – 303.
23. Chessells, J. M. et al. 1990. Neurotoxicity in lymphoblastic leukaemia: Comparison of oral and intramuscular methotrexate and two doses of radiation. Archives of Disease in Childhood 65:416 – 422.
24. Allen, J. C. et al. 1980. Leukoencephalopathy following high-dose IV methotrexate chemotherapy with leucovorin rescue. Cancer Treat Rep 64:1261 – 1273.
25. Jacobs, P. et al. 1991. Methotrexate encephalopathy. Eur J Cancer 27:1061 – 1062.
26. Flombaum, C. D. and Meyers, P. A. 1999. High-Dose Leucovorin as Sole Therapy for Methotrexate Toxicity. J Clin Oncol 17:1589 – 1594.
27. Widemann, B. C. et al. 2010. Glucarpidase, leucovorin, and thymidine for high-dose methotrexate-induced renal dysfunction: clinical and pharmacologic factors affecting outcome. J Clin Oncol 28:3979 – 3986.
28. Jonsson, O. G. and Kamen, B. A. 1991. Methotrexate and childhood leukemia. Cancer Investigation 9:53 – 60.
29. Saeter, G. et al. 1991. Treatment of osteosarcoma of the extremities with the T-10 protocol, with emphasis on the effects of preoperative chemotherapy with single-agent high-dose methotrexate: A scandinavian sarcoma group study. Journal of Clinical Oncology 9:1766 – 1775.
30. Abromowitch, M. et al. 1988. High-dose methotrexate improves clinical outcome in children with acute lymphoblastic leukemia: St. Jude total therapy study X. Medical and Pediatric Oncology 16:297 – 303.
31. Hann, I. M. et al. 1990. 'MACHO' chemotherapy for stage IV B cell lymphoma and B cell acute lymphoblastic leukaemia of childhood. British Journal of Haematology 76:359 – 364.
32. Wheeler, C. A. et al. 1991. Cisplatin, continuous infusion 5-fluorouracil, and intermediate dose methotrexate in the treatment of unresectable non-small cell carcinoma of the lung. Cancer 67:892 – 895.
33. Powles, T. J. et al. 1991. A randomized trial comparing combination chemotherapy using mitomycin C, mitozantrone and methotrexate (3M) with vincristine, anthracycline and cyclophosphamide (VAC) in advanced breast cancer. Br J Cancer 64:406 – 410.
34. Bablok, W. et al. 1988. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry. Part III. J. Clin Chem Clin Biochem 26:783 – 790.

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