Atellica™ IM
Analyzer

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1. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline−Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29‑A4.
2. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard−Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41‑A6.
3. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard −Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39‑A6.
4. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline−Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44‑A4.
5. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline−Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7‑A2.
6. 大野明美他：「ケンタウルスにおけるComplexed PSA 測定試薬「ケミルミACS-cPSA」の基礎的、臨床的有用性についての検討」、医学と薬学49（2）, 2003.
7. Clinical and Laboratory Standards Institute. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline−Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2000. CLSI Document C28‑A2.
8. Graves HC, Wehner N, Stamey TA. Comparison of a polyclonal and monoclonal immunoassay for PSA: need for an international antigen standard. *J Urol*. 1990;144（6）: 1516‑1522.
9. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45（7）:942‑956.

10. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F（ab'）2 conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38（9）:1737-1742.
11. Watt KW, Lee PJ, M'Timkulu T, et al. Human prostate-specific antigen: structural and functional similarity with serine proteases. *Proc Natl Acad Sci USA*. 1986;83（10）:3166-3170.
12. Christensson A, Laurell CB, Lilja H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine proteinase inhibitors. *Eur J Biochem*. 1990:194（3）:755-763.
13. Lilja H. Significance of different molecular forms of serum PSA. *Urol Clin North Am*. 1993;20（4）:681-686.
14. McCormack RT, Rittenhouse HG, Finlay JA, et al. Molecular forms of prostate-specific antigen and the human kallikrein gene family: a new era. *Urology*. 1995;45（5）:729-744.
15. Lilja H, Christensson A, Dahlen U, et al. Prostate-specific antigen in human serum occurs predominantly in complex with alpha 1-antichymotrypsin. *Clin Chem*. 1991;37（9）:1618-1625.
16. Stenman UH, Leinonen J, Alfthan H, et al. A complex between prostate-specific antigen and alpha 1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostatic cancer: assay of the complex improves clinical sensitivity for cancer. *Cancer Res*. 1991;51（1）:222-226.
17. Brawer MK, Meyer GE, Letran JL, et al. Measurement of complexed PSA improves specificity for early detection of prostate cancer. *Urology*. 1998;52（3）:372-378.
18. Brawer MK, Cheli CD, Neaman IE, et al. Complexed prostate specific antigen provides significant enhancement of specificity compared with total prostate specific antigen for detecting prostate cancer. *J Urol*. 2000;163（5）:1476-1480.
19. Okegawa T, Noda H, Nutahara K, Higashihara E. Comparison of two investigative assays for the complexed prostate-specific antigen in total prostate-specific antigen between 4.1 and 10.0 ng/mL. *Urology*. 2000;55（5）:700-704.
20. Mitchell ID, Croal BL, Dickie A, et al. A prospective study to evaluate the role of complexed prostate specific antigen and free/total prostate specific antigen ratio for the diagnosis of prostate cancer. *J Urol*. 2001;165（5）:1549-1553.
21. Miller MC, O'Dowd GJ, Partin AW, Veltri RW. Contemporary use of complexed PSA and calculated percent free PSA for early detection of prostate cancer: impact of changing disease demographics. *Urology*. 2001;57（6）:1105-1111.
22. Okihara K, Cheli CD, Partin AW. Comparative analysis of complexed prostate specific antigen, free prostate specific antigen and their ratio in detecting prostate cancer. *J Urol*. 2002;167（5）:2017-2023.
23. Stamey TA, Yemoto CE. Examination of the 3 molecular forms of serum prostate specific antigen for distinguishing negative from positive biopsy: relationship to transition zone volume. *J Urol*. 2000;163（1）:119-126.
24. Filella X, Alcover J, Molina R, et al. Measurement of complexed PSA in the differential diagnosis between prostate cancer and benign prostate hyperplasia. *Prostate*. 2000;42（3）: 181-185.
25. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline−Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
26. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline−Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
27. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline−Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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