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1. Clinical and Laboratory Standards Institute (CLSI). Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard-Sixth Edition. CLSI document H3-A6 (ISBN 1-56238-650-6). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2007.
2. Clinical and Laboratory Standards Institute (CLSI). Tubes and additives for venous and capillary blood specimen collection; Approved Standard-Sixth Edition. CLSI document H1-A6 (ISBN1-56238-740-5). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2010.
3. Clinical and Laboratory Standards Institute (CLSI). Procedures for the handling and processing of blood specimens for common laboratory tests; Approved Guideline-Fourth Edition. CLSI document H18-A4 (ISBN 1- 56238-724-3). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2010.
4. Schreiber W, Rausch S, Lammers M. Choice of specimen for immunonephelometric protein assays. Clin Chim Acta 2005;355:S407.
5. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP7-A2 (ISBN 1-56238-584-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
6. Steinmetz J, Tarallo P, Fournier B, et al. Reference limits of apolipoprotein A-I and apolipoprotein B using an IFCC standardized immunonephelometric method. Eur J Clin Chem Clin Biochem 1995;33:337-42.
7. Naito HK. The clinical significance of apolipoprotein measurements. J Clin Immunoassay 1986; 9:11-20.
8. Kottke BA, Zinsmeister AR, Holmes DR Jr, et al. Apolipoproteins and coronary artery disease. Mayo Clin Proc 1986;61:313-20.
9. Sniderman AD. Apolipoprotein-B and apolipoprotein-Al as predictors of coronary artery disease. Can J Cardiol 1988; 4 Suppl A:24A-30A.
10. Sandkamp M. Apolipoprotein-Diagnostik: klinische Relevanz. Diagnose & Labor 1990;40:37-43.
11. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5- A2 (ISBN 1-56238-542-9). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
12. Clinical and Laboratory Standards Institute (CLSI). Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition. CLSI document EP9-A2 (ISBN 1-56238-472-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
13. Expected values specification in the Dimension Vista® Software.

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