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1. Clinical and Laboratory Standards Institute (CLSI). Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard-Fifth Edition. CLSI document H3-A5 (ISBN 1-56238-515-1). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
2. Clinical and Laboratory Standards Institute (CLSI). Tubes and additives for venous blood specimen collection; Approved Standard-Fifth Edition. CLSI document H1-A5 (ISBN 1-56238-519-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
3. Clinical and Laboratory Standards Institute (CLSI). Procedures for the handling and processing of blood specimens; Approved Guideline–Third Edition. CLSI document H18-A3 (ISBN 1-56238-555-0). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
4. Clinical and Laboratory Standards Institute (CLSI). Interference testing in clinical chemistry; Approved Guideline-Second Edition. CLSI document EP7-A2 (ISBN 1-56238-584-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
5. Gressner AM, Thomas L. Proteinstoffwechsel. In: Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie. Stuttgart; New York: Schattauer, 1995:225.
6. Rask L, Anundi H, Bohme J, et al. The retinol-binding protein. Scand J Clin Lab Invest Suppl 1980;154:45-61.
7. Haider M, Haider SQ. Assessment of protein-calorie malnutrition. Clin Chem 1984;30:1286-99.
8. Bernstein LH, Leukhardt-Fairfield CJ, Pleban W, Rudolph R. Usefulness of data on albumin and prealbumin concentrations in determining effectiveness of nutritional support. Clin Chem 1989;35:271-4.
9. Clinical and Laboratory Standards Institute (CLSI). Evaluation of precision performance of quantitative measurement methods; Approved Guideline-Second Edition. CLSI document EP5-A2 (ISBN 1-56238-542-9). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
10. Clinical and Laboratory Standards Institute (CLSI). Method comparison and bias estimation using patient samples; Approved Guideline-Second Edition. CLSI document EP9-A2 (ISBN 1-56238-472-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
11. Clinical and Laboratory Standards Institute (CLSI). Protocols for determination of limits of detection and limits of quantitation; Approved Guideline. CLSI document EP17-A (ISBN 1-56238-551-8). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

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