DF115

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1. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fourth Edition. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 2003.
2. Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 2003.
3. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Standard—Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
4. Shaw LM, Holt DW, Oellerich M, et.al. Current issues in therapeutic drug monitoring of mycophenolic acid: report of a roundtable discussion. Ther Drug Monit. 2001; 23: 305
5. Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline. CLSI/NCCLS document EP7-A2 (ISBN 1-56238-584-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA 2002.
6. van Gelder T, LeMeur Y, Shaw L, et.al. Therapeutic Drug Monitoring of Mycophenolate Mofetil in Transplantation. Ther. Drug Monit. 2006; 28 (2): 145-154.
7. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. Clin. Chem 1999; 45:7;942-956.
8. Vaidya HC, Beatty BG, Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')2 conjugate and polyclonal mouse IgG, Clin Chem 1992; 38:1737-1742.
9. Srinivas T, Kaplan B and Meier-Kriesche H., Mycophenolate mofetil in solid-organ transplantation. Expert Opinions in Pharmacotherapy. 2003 4(12): 2325-2345.
10. European Mycophenolate Mofetil Cooperative Study Group. Placebo-controlled study of mycophenolate mofetil combined with cyclosporin and corticosteroids for prevention of acute rejection. Lancet. 1995;345/8961:1321-1325.
11. Sollinger HW for the US Renal Transplant Mycophenolate Mofetil Study Group. Mycophenolate mofetil for the prevention of acute rejection in primary cadaveric renal allograft recipients. Transplantation.1995;60/3:225-232.
12. Kobashigawa J, Miller L, Renlund D, et al. A randomized active-controlled trial of mycophenolate mofetil in heart transplant recipients. Transplantation.1998;66/4:507-515.
13. Allison AC, Kowalski WJ, Muller CD, et al. Mechanisms of action of mycophenolic acid. In: Allison AC, Lafferty K, Fliri H, eds. Immunosuppressive and Antiinflammatory Drugs. New York, NY: The New York Academy of Sciences; 1993;696:63-87.
14. Wu JC. Mycophenolate mofetil: Molecular mechanisms of action. Perspect Drug Discovery Design. 1994;2:185-204.
15. Shaw LM, Nicholls A, Hale M, et al. Therapeutic monitoring of mycophenolic acid: A consensus panel report. Clin Biochem.1998;31/5:317-322.
16. Bullingham RES, Nicholls AJ, Kamm BR. Clinical pharmacokinetics of mycophenolic mofetil. Clin Pharmacokinet.1998;34 /6:429-455.
17. Schütz E, Shipkova M, Armstrong VW, et al. Therapeutic drug monitoring of mycophenolic acid: comparison of HPLC and immunoassay reveals new MPA metabolites. Transplant Proc. 1998;30:1185-1187.
18. Shipkova M, Armstrong VW, Wieland E, Niedmann PD, Schütz E, Brenner-Weiss G, Voihsel M, Braun F, Oellerich M. Identification of glucoside and carboxyl-linked glucuronide conjugates of mycophenolic acid in plasma of transplant recipients treated with mycophenolate mofetil. Br J Pharmacol.1999;126/5:1075-1082.
19. Shaw LM, Korecka M, Venkataramanan R, et.al. Mycophenolic acid pharmacodynamics and pharmacokinetics provide a basis for rational monitoring strategies. Am J Transplant 2003; 3:534.
20. Nowak I, Shaw LM. Mycophenolic acid binding to human serum albumin: Characterization and relation to pharmacodynamics. Clin Chem.1995:41/7:1011-1017.
21. van Gelder T and Shaw LM. The Rationale for and limitations of therapeutic drug monitoring for mycophenolate mofetil in transplantation. Transplantation. 2005; 80 (2S); S244-253.
22. Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004.
23. Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 2002.
24. Clinical and Laboratory Standards Institute/NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. CLSI/NCCLS document EP17-A (ISBN 1-56238-551-8). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA 2004.

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