Atellica® CH
Analyzer

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1. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline−Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29‑A4.
2. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard−Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41‑A6.
3. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard−Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39‑A6.
4. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline−Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44‑A4.
5. Wu AHB. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Saunders; 2006:146.
6. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline−Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7‑A2.
7. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline−Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28‑A3c.
8. Tietz NW. *Clinical Guide to Laboratory Tests*. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:68-71.
9. Brown SA, Hutchinson R, Morrisett J, *et al*. Plasma lipid, lipoprotein cholesterol, and apoprotein distributions in selected US communities. *Arterioscler Thromb*. 1993;13（8）: 1139-1158.
10. Alaupovic P. Structure and function of plasma lipoproteins with particular regard to hyperlipoproteinemias and atherosclerosis. *Ann Biol Clin*. 1980;38（2）:83-93.
11. Naito HK. The clinical significance of apolipoprotein measurements. *J Clin Immunoassay*. 1986;9:11-20.
12. Burtis CA, Ashwood ER, Bruns DE, eds. *Tietz Fundamentals of Clinical Chemistry*. 6th ed. St. Louis, MO: Saunders; 2008:413-427.
13. Riesen WF, Mordasini R, Salzmann C, Theler A, Gurtner HP. Apoproteins and lipids as discriminators of severity of coronary heart disease. *Atherosclerosis*. 37（1）:157-162.
14. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline−Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09‑A3.
15. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline−Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05‑A3.
16. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline−Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17‑A2.
17. Data on file at Siemens Healthcare Diagnostics.

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