Atellica® CH
Analyzer

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1. You K, Bittikofer JA. Quantification of salicylate in serum by use of salicylate hydroxylase. *Clin.Chem*. 1984;30（9）:1549-1551.
2. Hammond PM. A simple colorimetric assay to determine salicylate ingestion utilizing salicylate monooxygenase. Paper presented at: VIIIth International Conference on Enzyme Engineering; 1985.
3. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline−Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29‑A4.
4. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard−Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41‑A6.
5. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard−Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39‑A6.
6. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline−Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44‑A4.
7. World Health Organization. *Use of Anticoagulants in Diagnostic Laboratory Investigations*. Geneva; World Health Organization: 2002;1-62. Government Document WHO/DIL/LAB/99.1 Rev 2.
8. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline−Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7‑A2.
9. Data on file at Siemens Healthcare Diagnostics.
10. Henry JB. *Clinical Diagnosis and Management by Laboratory Methods*. 19th ed. Philidelphia, PA: W.B. Saunders Co.; 1996:1454.
11. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline−Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09‑A3.
12. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline−Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05‑A3.
13. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline−Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17‑A2.

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