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1. Wolde-Mariam W, Peter J. Recent diagnostic advances in cellular immunology. Diagnost Med. 1984;7:25-32.
2. Henry C, Mishell B, Shiigi S. Hemolytic plaque assays. In: Mishell B, Shiigi S, ed. Selected Methods in Cellular Immunology. New York: WH Freeman and Co; 1980:69-123.
3. Grabstein K, Chen Y. Cell-mediated cytolytic responses. In: Mishell B, Shiigi S, ed. Selected Methods in Cellular Immunology. New York: WH Freeman and Co; 1980:124-137.
4. Wofsy L, Henry C, Kimura J, North J. Modification and use of antibodies to label cell surface antigens. In: Mishell B, Shiigi S, ed. Selected Methods in Cellular Immunology. New York: WH Freeman and Co; 1980:287-304.
5. Jackson A, Warner N. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose N, Friedman H, Fahey J, ed. Manual of Clinical Laboratory Immunology. Washington, DC: American Society for Microbiology; 1986:226-235.
6. Jackson A. Basic phenotyping of lymphocytes: selection and testing of reagents and interpretation of data. Clin Immunol Newsletter.

1990;10:43-55.
7. Nicholson J, Jones B. Immunophenotyping by flow cytometry: its use in HIV infection. Labmedica. 1989;6:21-26.
8. Landay A, Muirhead K. Procedural guidelines for performing immunophenotyping by flow cytometry. Clin Immunol Immunopathol. 1989;52:48-60.
9. Schmidt R. Monoclonal antibodies for diagnosis of immunodeficiencies. Blut. 1989;59:200-206.
10. Ohno T, Kanoh T, Suzuki T, et al. Comparative analysis of lymphocyte phenotypes between carriers of human immunodeficiency virus (HIV) and adult patients with primary immunodeficiency using two-color immunofluorescence flow cytometry. J Exp Med.1988;154:157.
11. Antel J, Bania M, Noronha A, Neely S. Defective suppressor cell function mediated by T8+ cell lines from patients with progressive multiple sclerosis. J Immunol. 1986;137:3436-3439.
12. Smolen J, Chused T, Leiserson W, Reeves J, Alling D, Steinberg A. Heterogeneity of immunoregulatory T-cell subsets in systemic lupus erythematosus: Correlation with clinical features. Am J Med. 1982;72:783-790.
13. Legendre C, Schiffrin A, Weitzner G, Colle E, Guttmann R. Two color flow cytometry analysis of activated T-lymphocyte subsets in type 1 diabetes mellitus. Diabetes. 1988;37:792-795.
14. Tomkinson B, Wagner D, Nelson D, Sullivan J. Activated lymphocytes during acute Epstein-Barr virus infection. J Immunol. 1987;139:3802-3807.
15. Bishop G, Hall B, Duggin G, Horvath J, Sheil A, Tiller D. Immunopathology of renal allograft rejection analyzed with monoclonal antibodies to mononuclear cell markers. Kidney Internat. 1986;29:708-717.
16. Lanier L, Le A, Phillips J, Warner N, Babcock G. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. J Immunol. 1983;131:1789-1796.
17. Giorgi J, Hultin L. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. Clin Immunol Newsletter. 1990;10:55-61.
18. Gratama J, Naipal A, Oljans P, et al. T lymphocyte repopulation and differentiation after bone marrow transplantation: Early shifts in the ratio between T4+ and T8+ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. Blood. 1984;63:1416-1423.
19. Stites D, Casavant C, McHugh T, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. Clin Immunol Immunopath. 1986;38:161-177.
20. Lanier L, Le A, Civin C, Loken M, Phillips J. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic Tlymphocytes. J Immunol. 1986;136:4480-4486.
21. Fitzgerald-Bocarlsy P, Herberman R, Hercend T, et al. A definition of natural killer cells. In: Ades E, Lopez C, ed. Natural Killer Cells and Host Defense. Basel: Karger; 1989:1.
22. Aiuti F, Sirianni M, Mezzaroma I, et al. HIV-1 infection: epidemiological features and immunological alterations during the natural history of the disease. Clin Immunol Immunopath. 1989;50:S157-S165.
23. Jandl J. Blood and blood forming tissues. In: Jandl J, ed. Blood Textbook of Hematology. Boston, Mass: Little Brown and Company; 1987:1-48.
24. Beverley P. Production and use of monoclonal antibodies in transplantation immunology. In: Touraine J, Trager J, Betuel H, et al., ed. Transplantation and Clinical Immunology XI. Amsterdam: Excerpta Medica; 1980:87-94.
25. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael A, ed. Leucocyte Typing III. New York, NY: Oxford University Press; 1986:788-803.
26. Rowe D, Beverley P. Characterisation of breast cancer infiltrates using monoclonal antibodies to human leucocyte antigens. Br J Cancer. 1984;49:149-159.
27. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: M2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman S, ed. Leucocyte Typing. Berlin: Springer-Verlag; 1984:82-108.
28. Dimitriu-Bona A, Burmester G, Waters S, Winchester R. Human mononuclear phagocyte differentiation antigens. I. Patterns of antigenic expression on the surface of human monocytes and macrophages defined by monoclonal antibodies. J Immunol. 1983;130:145-152.
29. Herrmann F, Komischke B, Odenwald E, Ludwig W. Use of monoclonal antibodies as a diagnostic tool in human leukemia. I. Acute myeloid leukemia and acute phase of chronic myeloid leukemia. Blut. 1983;47:157-163.
30. Schwinzer R. Cluster report: CD45/CD45R. In: Knapp W, ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:628-634.
31. Goyert S, Tesio L, Ashman L, et al. Report on the CD14 Cluster Workshop. In: Knapp W, ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:789-794.
32. Bernstein I, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leukocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, ed. Leukocyte Typing II: Human Myeloid and Hematopoietic Cells. New York: Springer-
Verlag; 1984:1-25.
33. Jayaram Y, Hogg N. Surface expression of CD14 molecules on human neutrophils. In: Knapp W, Dorken B, Gilks W, et al., ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:796-797.
34. Haynes B. Summary of T cell studies performed during the Second International Workshop and Conference on Human Leukocyte

Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, ed. Leukocyte Typing II: Human T Lymphocytes. New York: Springer-Verlag; 1986:3-30.
35. Ledbetter J, Evans R, Lipinski M, Cunningham-Rundles C, Good R, Herzenberg L. Evolutionary conservation of surface molecules that distinguish T lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. J Exp Med. 1981;153:310-323.
36. Kan E, Wang C, Wang L, Evans R. Non-covalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. J Immunol. 1983;131:536-539.
37. Knowles R. Immunochemical analysis of the T cell-specific antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, ed. Leukocyte Typing II: Human T Lymphocytes. New York: Springer-Verlag; 1986:259-288.
38. Nadler L. B Cell/Leukemia Panel Workshop: Summary and comments. In: Reinherz E, Haynes B, Nadler L, Bernstein I, ed. Leukocyte Typing II: Human B Lymphocytes. New York: Springer-Verlag; 1986:3-43.
39. van Dongen J, Krissansen G, Wolvers-Tettero I, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. Blood. 1988;71:603-612.
40. Lanier L, Federspiel N, Ruitenberg J, et al. The T cell antigen receptor complex expressed on normal peripheral blood CD4-, CD8- T lymphocytes. J Exp Med. 1987;165:1076-1094.
41. Dorken B, Moller P, Pezzutto A, Schwartz-Albiez R, Molderhauer G. B-cell antigens: CD19. In: Knapp W, ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1986:34-36.
42. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausett J, Milstein C, Schlossman S, ed. Leucocyte Typing. Berlin: Springer-Verlag; 1984:25-60.
43. Evans R, Wall D, Platsoucas C, et al. Thymus-dependent membrane antigens in man: Inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. Proc Natl Acad Sci USA. 1981;78:544-548.
44. Maddon P, Dalgleish A, McDougal J, Clapham P, Weiss R, Axel R. The T4 gene encodes the AIDS virus receptor and is expressed in the immune system and the brain. Cell. 1986;47:333-348.
45. Dalgleish A, Beverly P, Clapham P, Crawford D, Greaves M, Weiss R. The CD4 (T4) antigen is an essential component of the receptor for the AIDS virus. Nature. 1984;312:763-767.
46. Kurrle R. Cluster report: CD3. In: Knapp W, ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:290-293.
47. Wood G, Warner N, Warnke R. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. J Immunol. 1983;131:212-216.
48. Anderson P, Blue M-L, Morimoto C, Schlossman S. Cross-linking of T3 (CD3) with T4 (CD4) enhances the proliferation of resting T lymphocytes. J Immunol. 1987;139:678-682.
49. Eichmann K, Johnson J-I, Falk I, Emmrich F. Effective activation of resting mouse Tlymphocytes by cross-linking submitogenic concentrations of the T cell antigen receptor with either Lyt-2 or L3T4. Eur J Immunol. 1987;17:643-650.
50. Gallagher P, Fazekas de St. Groth B, Miller J. CD4 and CD8 molecules can physically associate with the same T-cell receptor. Proc Natl Acad Sci USA. 1989;86:10044-10048.
51. Rudd C, Burgess K, Barber E, Schlossman S. Monoclonal antibodies to the CD4 and CD8 antigens precipitate variable amounts of CD4/CD8-associated p56-lck activity. In: Knapp W, Dorken B, Gilks W, et al., ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:326-327.
52. Terry L, Disanto J, Small T, Flomenberg N. Differential expression of the CD8 and Lyt-3 antigens on a subset of human T-cell receptor γ/δ-bearing lymphocytes. In: Knapp W, Dorken B, Gilks W, et al., ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:345-346.
53. Moebius U. Cluster report: CD8. In: Knapp W, Dorken B, Gilks W, et al., ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:342-343.
54. Lampson L, Levy R. Two populations of Ia-like molecules on a human B cell line. J Immunol. 1980;125:293-299.
55. Brodsky F. A matrix approach to human Class II histocompatibility antigens: reactions of four monoclonal antibodies with the products of nine haplotypes. Immunogenetics. 1984;19:179-194.
56. Robbins P, Evans E, Ding A, Warner N, Brodsky F. Monoclonal antibodies that distinguish between Class II antigens (HLA-DP, DQ, and DR) in 14 haplotypes. Human Immunol. 1987;18:301-313.
57. Terstappen L, Hollander Z, Meiners H, Loken M. Quantitative comparison of myeloid antigens on five lineages of mature peripheral blood cells. J Leuk Biol. 1990;48:138-148.
58. Warnke R, Levy R. Detection of T and B antigens with hybridoma monoclonal antibodies: a biotin-avidin-horseradish peroxidase method. J Histochem Cytochem. 1980;28:771-776.
59. Engleman E, Warnke R, Fox R, Dilley J, Benike C, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. Proc Natl Acad Sci USA. 1981;78:1791-1795.
60. Edwards J, Durant B, Jones D, Evans P, Smith J. Differential expression of HLA Class II antigens in fetal human spleen: relationship of HLA-DP, DQ, and DR to immunoglobulin expression. J Immunol. 1985;137:479-490.
61. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. I. Characterization of the lymphocyte subset reactive with B73.1. J Immunol. 1983;130:2133-2141.
62. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. II. Studies of B73.1 antibody-antigen interaction on the lymphocyte membrane. J Immunol. 1983;130:2142-2148.

63. Perussia B, Trinchieri G, Jackson A, et al. The Fc receptor for IgG on human natural killer cells: phenotypic, functional, and comparative studies with monoclonal antibodies. J Immunol. 1984;133:180-189.
64. Schmidt R. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dorken B, Gilks W, et al., ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:517-542.
65. Schubert J, Lanier L, Schmidt R. Cluster report: CD56. In: Knapp W, Dorken B, Gilks W, et al., ed. Leucocyte Typing IV. New York, NY: Oxford University Press; 1989:699-702.
66. Knapp W, Rieber P, Dörken B, Schmidt R, Stein H, van der Borne A. Towards a better definition of human leucocyte surface molecules. Immunol Today. 1989;10:253-258.
67. Lanier L, Testi R, Bindl J, Phillips J. Identity of Leu-19 (CD56) leukocyte differentiation antigen and neural cell adhesion molecule. J Exp Med. 1989;169:2233-2238.
68. Slockbower J, Belgeri K, Bruck E, et al. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture (H3-A2). Villanova: The National Committee for Clinical Laboratory Standards, 1984.
69. Reichert T, DéBruyere M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. Clin Immunol Immunopath. 1991;60:190-208.
70. Prince H, Hirji K, Waldbeser L, Plaeger-Marshall S, et al. Influence of racial background on the distribution of T cell subsets and Leu-11-positive lymphocytes in healthy blood donors. Diag Immunol. 1985;3:33-37.

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