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1. van Dongen JJM, Krissansen GW, Wolvers-Tettero ILM, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. Blood. 1988;71:603-612.
2. Brenner MB, Groh V, Porcelli SA, et al. Structure and distribution of the human γ/δ T-cell receptor. In: Knapp W, Dörken B, Gilks WR, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens.New York, NY: Oxford University Press; 1989:1049-1053.
3. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. Annu Rev Immunol. 1988;6:629-662.
4. Reichert T, DeBruyére M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. Clin Immunol Immunopath. 1991;60:190-208.
5. Lanier LL, Allison JP, Phillips JH. Correlation of cell surface antigen expression on human thymocytes by multi-color flow cytometric analysis: implications for differentiation. J Immunol. 1986;137:2501-2507.
6. Garson JA, Beverley PCL, Coakham HB, Harper EJ. Monoclonal antibodies against human T lymphocytes label Purkinje neurones of many species. Nature. 1982;298:375-377.
7. Kaneoka H, Perez-Rojas G, Sasasaki T, Benike CJ, Engleman EG. Human T lymphocyte proliferation induced by a pan-T monoclonal antibody (Anti-Leu-4): heterogeneity of response is a function of monocytes. J Immunol. 1983;131:158.
8. Ledbetter JA, Evans RL, Lipinski M, Cunningham-Rundles C, Good RA, Herzenberg LA. Evolutionary conservation of surface molecules that distinguish T-lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. J Exp Med. 1981;153:310-323.
9. Kan EAR, Wang CY, Wang LC, Evans RL. Noncovalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. J Immunol. 1983;131:536-539.
10. Knowles RW. Immunochemical analysis of the T-cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. Leukocyte Typing II: Human T Lymphocytes. New York, NY: Springer-Verlag; 1986;1:259-288.
11. Haynes BF. Summary of T-cell studies performed during the Second International Workshop and Conference on Human Leukocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. Leukocyte Typing II: Human T Lymphocytes. New York, NY: Springer-Verlag; 1986;1:3-30.
12. Allison JP, Lanier LL. Structure, function, and serology of the T-cell antigen receptor complex. Ann Rev Immunol. 1987;5:503.
13. Freidrich W, O'Reilly RJ, Kozenger B, Gebhard DF, Good RA, Evans RL. T lymphocyte reconstitution in recipients of bone marrow transplants: analysis of regulatory T cell imbalances in GVHD. Blood. 1982;59:696.
14. Kurrle R. Cluster report: CD3. In: Knapp W, Dörken B, Gilks WR, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:290-293.
15. Maecker HT, Maino VC. Analyzing T-cell responses to cytomegalovirus by cytokine flow cytometry. Hum Immunol. 2004;65:493-499.
16. Maino VC, Picker LJ. Identification of functional subsets by flow cytometry: intracellular detection of cytokine expression. Cytometry. 1998;34:207-215.
17. Campbell MJ, Scott J, Maecker HT, Park JW, Esserman LJ. Immune dysfunction and micrometastases in women with breast cancer. Breast Cancer Res Treat. 2005;91:163-171.
18. Tu W, Chen S, Sharp M, et al. Persistent and selective deficiency of CD4+ T cell immunity to cytomegalovirus in immunocompetent young children. J Immunol. 2004;172:3260-3267.
19. Wolde-Mariam W, Peter JB. Recent diagnostic advances in cellular immunology. Diag Med. 1984;7:25-32.
20. Schmidt RE. Monoclonal antibodies for diagnosis of immunodeficiencies. Blut. 1989;59:200-206.
21. Giorgi JV, Hultin LE. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. Clin Immunol Newslett. 1990;10:55-61.
22. Bernard A, Boumsell L, Hill C. Joint report of the first international workshop on human leucocyte differentiation antigens by the investigators of the participating laboratories. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman SF, eds. Leucocyte Typing. New York, NY: Springer-Verlag; 1984:9-108.
23. Evans RL, Wall DW, Platsoucas CD, et al. Thymus-dependent membrane antigens in man: inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. Proc Natl Acad Sci USA. 1981;78:544-548.
24. Engleman EG, Benike CJ, Glickman E, Evans RL. Antibodies to membrane structures that distinguish suppressor/cytotoxic and helper T lymphocyte subpopulations block the mixed leukocyte reaction in man. J Exp Med. 1981;153:193-198.
25. Kotzin BK, Benike CJ, Engleman EG. Induction of immunoglobulin secreting cells in the allogeneic mixed leukocyte reaction: regulation by helper and suppressor lymphocyte subsets in man. J Immunol. 1981;127:931-935.
26. Wood GS, Warner NL, Warnke RA. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. J Immunol. 1983;131:212-216.
27. Dalgleish AG, Beverley PCL, Clapham PR, Crawford DH, Greaves MF, Weiss RA. The CD4 (T4) antigen is an essential component of the receptor for the AIDS retrovirus. Nature. 1984;312:763-767.
28. Sattentau QJ, Dalgleish AG, Weiss RA, Beverley PCL. Epitopes of the CD4 antigen and HIV infection. Science. 1986;234:1120-1123.
29. Lewis DE, Puck JM, Babcock GF, Rich R. Disproportionate expansion of a minor T cell subset in patients with lymphadenopathy syndrome and acquired immunodeficiency syndrome. J Infect Dis. 1985;151:555-559.
30. Ohno T, Kanoh T, Suzuki T, et al. Comparative analysis of lymphocyte phenotypes between carriers of human immunodeficiency

virus (HIV) and adult patients with primary immunodeficiency using two-color immunofluorescence flow cytometry. J Exp Med. 1988;154:157-172.
31. Stites DP, Casavant CH, McHugh TM, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. Clin Immunol Immunopathol. 1986;38:161-177.
32. Gatenby PA, Kansas GS, Xian CY, Evans RL, Engleman EG. Dissection of immunoregulatory subpopulations of T lymphocytes within the helper and suppressor sublineages in man. J Immunol. 1982;129:1997-2000.
33. Picker LJ, Weiss LM, Mederios LJ, Wood GS, Warnke RA. Immunophenotypic criteria for the diagnosis of non-Hodgkin's lymphoma. Am J Pathol. 1987;128:181-201.
34. Maecker HT, Maino VC, Picker LJ. Immunofluorescence analysis of T-cell responses in health and disease. J Clin Immunol. 2000;20:391399.
35. Knowles RW. Immunochemical analysis of the T-cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. Leukocyte Typing II: Human T Lymphocytes. New York, NY: Springer-Verlag; 1986;1:259-288.
36. Engleman EG, Warnke R, Fox RI, Dilley J, Benike CJ, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. Proc Natl Acad Sci USA. 1981;78:1791-1795.
37. Ledbetter JA, Frankel AE, Herzenberg LA. Human Leu T-cell differentiation antigens: quantitative expression on normal lymphoid cells and cell lines. In: Hämmerling G, Hämmerling U, Kearney J, eds. Monoclonal Antibodies and T Cell Hybridomas: Perspectives and Technical Notes. New York, NY: Elsevier/North Holland; 1981:16-22.
38. Warnke RA, Levy R. Detection of T and B antigens with hybridoma monoclonal antibodies: a biotinavidin-horseradish peroxidase method. J Histochem Cytochem. 1980;28:771-776.
39. Warnke R, Miller R, Grogan T, Pederson M, Dilley J, Levy R. Immunologic phenotype in 30 patients with diffuse large-cell lymphoma. N Eng J Med. 1980;303:293-300.
40. Zipf RF, Fox R, Dilley J, Levy R. Definition of the high risk ALL patient by immunologic phenotyping with monoclonal antibodies. Cancer Res. 1981;41:4786.
41. Gadol N, Ault KA. Phenotypic and functional characterization of human Leu-1 (CD5) B cells. Immunol Rev. 1986;93:23.
42. Royston I, Majda JA, Baird SM, Meserve BL, Griffiths JC. Human T-cell antigens defined by monoclonal antibodies: the 65,000-dalton antigen of T cells (T65) is also found on chronic lymphocytic leukemia cells bearing surface immunoglobulin. J Immunol. 1980;125:725.
43. Antin JH, Emerson SG, Martin P, Gadol N, Ault KA. Leu-1+ (CD5) B cells. A major lymphoid subpopulation in human fetal spleen: Phenotypic and functional studies. J Immunol. 1986;136:505.
44. Hardy RR, Hayakawa K, Shimizu M, Yamasaki K, Kishimoto T. Rheumatoid factor secretion from human Leu-1+ B cells. Science. 1987;236:81-83.
45. Lanier LL, Le AM, Civin CI, Loken MR, Phillips JH. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. J Immunol. 1986;136:4480-4486.
46. Lanier LL, Kipps T, Phillips JH. Functional properties of a unique subset of cytotoxic CD3+ T lymphocytes that express Fc receptors for IgG (CD16/Leu-11 antigen). J Exp Med. 1985;162:2089.
47. Lanier LL, Le AM, Phillips JH, Warner NL, Babcock GF. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. J Immunol. 1983;131:1789-1796.
48. Perussia B, Trinchieri G, Jackson A, et al. The Fc receptor for IgG on human natural killer cells: phenotypic, functional, and comparative studies with monoclonal antibodies. J Immunol. 1984;133:180-189.
49. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. I. Characterization of the lymphocyte subset reactive with B73.1. J Immunol. 1983;130:2133-2141.
50. Gerosa F, Baldani-Guerra B, Nisii C, Marchesini V, Carra G, Trinchieri G. Reciprocal activating interaction between natural killer cells and dendritic cells. J Exp Med. 2002;195:327-333.
51. Phillips JH, Babcock GF. A monoclonal antibody reactive against purified human natural killer cells and granuocytes. Immunol Letters. 1983;6:143.
52. Thompson J, Goeken N, Brown S, Rhodes J. Phenotypic definition of human monocyte/macrophage subpopulations by monoclonal antibodies. II. Distinction of non-T, non-B natural killer cells (NK) from antigen-presenting stimulating cells in the mixed lymphocyte response (MLR). Presented at the American Association for Clinical Histocompatibility Testing;
53. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. II. Studies of B73.1 antibody-antigen interaction on the lymphocyt membrane. J Immunol. 1983;130:2142-2148.
54. Schmidt RE. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dörken B, Gilks WR, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:517-542.
55. Cooper MD, Fehniger TA, Caligiuri MA. The biology of human natural killer-cell subsets. Trends in Immunol. 2001;22:633-640.
56. Campbell JJ, Qin S, Unutmaz D, et al. Unique subpopulations of CD56+ NK and NK-T peripheral blood lymphocytes identified by chemokine receptor expression repertoire. J Immunol. 2001;166:6477-6482.
57. Galandrini R, Tassi I, Mattia G. SH2-containing inositol phosphatase (SHIP-1) transiently translocates to raft domains and modulates CD16-mediated cytotoxicity in human NK cells. Blood. 2002;100:4581-4589.
58. Fröland S, Natvig J, Berdal P. Surface-bound immunoglobulin as a marker of B lymphocytes in man. Nature. 1971;234:251-252.
59. Warner N. Membrane immunoglobulins and antigen receptors on B and T lymphocytes. Advances Immunol. 1974;19:67.
60. Kohler G, Milstein C. Continuous cultures of fused cells secreting antibody of predefined specificity. Nature. 1975;256:495.

61. Dörken B, Möller P, Pezzutto A, Schwartz-Albiez R, Moldenhauer G. B-cell antigens: CD19. In: Knapp W, Dörken B, Gilks WR, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford: Oxford University Press; 1989:34-36.
62. Loken M, Shah V, Dattilio K, Civin C. Flow cytometric analysis of human bone marrow. II. Normal B-lymphocyte development. Blood. 1987;70:1316-1324.
63. Ryan D, Kossover S, Mitchell S, Frantz C, Hennessy L, Cohen H. Subpopulations of common acute lymphoblastic leukemia antigen-positive lymphoid cells in normal bone marrow identified by hematopoietic differentiation antigens. Blood. 1986;68:417-425.
64. Foucar K, Goeken J. Clinical application of immunologic techniques to the diagnosis of lymphoproliferative and immunodeficiency disorders. Lab Medicine. 1982;13:403-413.
65. Meeker T, Miller R, Link M, Bindl J, Warnke R, Levy R. A unique human B lymphocyte antigen defined by a monoclonal antibody. Hybridoma. 1984;3:305-320.
66. Clark P, Normansell D, Innes D, Hess C. Lymphocyte subsets in normal bone marrow. Blood. 1986;67:1600.
67. Nadler L. B cell/leukemia panel workshop: Summary and comments. In: Reinherz E, Haynes B, Nadler L, ID B, eds. Leukocyte Typing II: Human B Lymphocytes. New York: Springer-Verlag; 1986:3-43.
68. Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens on health-care settings. MMWR. 1988; 37(24).
69. National Committee for Clinical Laboratory Standards. Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue: Tentative Guideline (M29-T2). Villanova, PA: National Committee for Clinical Laboratory Standards; 1991.
70. Centers for Disease Control. Guidelines for the performance of CD4+ T-cell determinations in persons with human immunodeficiency virus infection. MMWR. 1992; 41(No. RR-8).
71. National Committee for Clinical Laboratory Standards. Clinical Applications of Flow Cytometry: Quality Assurance and Immunophenotyping of Peripheral Blood Lymphocytes (H42-T). Villanova, PA: National Committee for Clinical Laboratory Standards; 1992.
72. Giorgi J. Lymphocyte subset measurements: Significance in clinical medicine. In: Rose N, Friedman H, Fahey J, eds. Manual of Clinical Laboratory Standards. Washington DC: American Society of Microbiology; 1986:236-246.
73. Mishell B, Shiigi S, Henry C, et al. Preparation of mouse cell suspensions. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York: WH Freeman and Co; 1980:16-17.
74. Jackson A, Warner N. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose N, Friedman H, Fahey J, eds. Manual of Clinical Laboratory Immunology. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
75. Tursz T, Brouet J-C, Flandrin G, Danon F, Clauvel J-P, Seligmann M. Clinical and pathologic features of Waldenström's macroglobulinemia in seven patients with serum monoclonal IgG or IgA. Am J Med. 1977;63:499-502.
76. Brown G, Greaves M. Enumeration of absolute numbers of T and B cells in human blood. Scand J Immunol. 1974;3:161.
77. Dwyer J, Finklestein F, Mangi R, Fisher K, Hendler E. Assessment of the adequacy of immunosuppressive therapy using microscopy techniques to study immunologic competence. Transplant Proc. 1975;7:785.
78. Dörken B, Möller P, Pezzutto A, Schwartz-Albiez R, Moldenhauer G. B-cell antigens: CD20. In: Knapp W, Dörken B, Gilks WR, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. 1 ed. New York, NY: Oxford University Press; 1989:46-48.
79. Hultin LE, Hausner MA, Hultin PM, Giorgi JV. CD20 (Pan-B cell) antigen is expressed at a low level on a subpopulation of human T lymphocytes. Cytometry. 1993;14:196-204.
80. Ling NR, Maclennan ICM, Mason DY. B-cell and plasma cell antigens: new and previously defined clusters. In: McMichael AJ, ed. Leucocyte Typing III: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1987:303-334.
81. Clark EA, Shu G, Ledbetter JA. Role of the B p35 cell surface polypeptide in human B cell activation. Proc Natl Acad Sci USA. 1985;82:1766.
82. Abts H, Emmerich M, Miltenyi S, Radbruch A, Tesch H. CD20 positive human B lymphocytes separated with the magnetic cell sorter (MACS) can be induced to proliferation and antibody secretion in vitro. J Immunol Methods. 1989;125:19-28.
83. Marti GE, Faguet G, Bertin P, et al. CD20 and CD5 expression in B-chronic lymphocytic leukemia. Ann NY Acad Sci. 1992;651:480-483.
84. Deans JP, Li H, Polyak MJ. CD20-mediated apoptosis: signalling through lipid rafts. Immunology. 2002;107:176-182.
85. Lanier LL, Testi R, Bindl J, Phillips JH. Identity of the Leu-19 (CD56) leukocyte differentiation antigen and neural cell adhesion molecule. J Exp Med. 1989;169:2233-2238.
86. Lanier LL, Chang C, Azuma M, Ruitenberg JJ, Hemperly JJ, Phillips JH. Molecular and functional analysis of human natural killer cell-associated neural cell adhesion molecule (N-CAM/CD56). J Immunol. 1991;146:4421-4426.
87. Schubert J, Lanier LL, Schmidt RE. Cluster report: CD56. In: Knapp W, Dörken B, Gilks WR, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York: Oxford University Press; 1989:699-702.
88. Hercend T, Griffin J, Benussan A, et al. Generation of monoclonal antibodies to a human natural killer clone. Characterization of two natural killer-associated antigens NKH1 and NKH2 expressed on subsets of large granular lymphocytes. J Clin Invest. 1985;75:932.
89. Edelman GM. Cell adhesion molecules. Science. 1983;219:450-457.
90. Phillips JH, Lanier LL. Dissection of the lymphokine-activated killer phenomenon: Relative contribution of peripheral blood natural killer cells and T lymphocytes to cytolysis. J Exp Med. 1986;164:814-825.
91. Abo T, Balch CMA. A differentiation antigen of human NK and K cells identified by a monoclonal antibody (HNK-1). J Immunol. 1981;127:1024-1029.

92. Manara DC, de Panfilis B, Ferrari C. Ultrastructural characterization of human large granular lymphocyte subsets defined by the expression of HNK-1 (Leu-7), Leu-11, or both HNK-1 and Leu-11 antigens. J Histochem Cytochem. 1985;33:1129-1133.
93. Mechtersheimer G, Staudter M, Möller P. Expression of the natural killer cell-associated antigens CD56 and CD57 in human neural and striated muscle cells and in their tumors. Cancer Res. 1991;51:1300-1307.
94. Kruse J, Mailhammer R, Wennecke H, et al. Neural cell-adhesion molecules and MAG share a common carbohydrate moiety recognized by monoclonal antibodies L2 and HNK-1. Nature. 1984;311:153-155.
95. McGarry RC, Helfand SL, Quarles RH, Roder JC. Recognition of myelinassociated glycoprotein by the monoclonal antibody HNK-1. Nature. 1983;306:376.
96. Bunn PA, Gazdar AF, Carney DN, Minna J. Small cell lung carcinoma and natural killer cells share an antigenic determinant, Leu-7. Clin Res. 1984;32:416A.
97. Schubert J, Lanier LL, Schmidt RE. Cluster report: CD57. In: Knapp W, Dörken B, Gilks WR, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford: Oxford University Press; 1989:711-714.
98. Lewis DE, Puck JM, Babcock GM, Rich RR. Disproportionate expansion of a minor T cell subset in patients with lymphadenopathy syndrome and acquired immunodeficiency syndrome. J Infect Disease. 1985;151:555.
99. Prince HE, Kreiss JK, Kasper CK, et al. Distinctive lymphocyte subpopulation abnormalities in patients with congenital coagulation disorders who exhibit lymph node enlargement. Blood. 1985;66:64.
100. Stites D, Casavant C, McHugh T, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. Clin Immunol Immunopathol. 1986;38:161-177.
101. Wood GS, Warner N, Warnke R. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. J Immunol. 1983;131(1):212-216.
102. Giorgi J, Hultin L. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. Clin Immunol Newslett. 1990;10(4):55-62.
103. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman S, eds. Leucocyte Typing. Berlin: Springer-Verlag; 1984:25-60.
104. Picker LF, Weiss LM, Medeiros LJ, Wood GS, Warnke RA. Immunophenotypic criteria for the diagnosis of non-Hodgkin's lymphoma. Am J Pathol. 1987;128:181-201
105. Schwinzer R. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks WR, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:628-634.
106. Jackson A. Basic phenotyping of lymphocytes: selection and testing of reagents and interpretation of data. Clin Immunol Newslett. 1990;10:43-55.
107. Terstappen LWMM, Levin J. Bone marrow cell differential counts obtained by multidimensional flow cytometry. Blood Cells. 1992;18:311-330.
108. Loken MR, Brosnan J, Bach B, Ault K. Establishing optimal lymphocyte gates for immunophenotyping by flow cytometry. Cytometry. 1990;11:453-459.
109. Hermiston ML, Xu Z, Weiss A. CD45: a critical regulator of signaling thresholds in immune cells. Annu Rev Immunol. 2003;21:107-137.
110. Pellat-Deceunynck C, Bataille R. Normal and malignant human plasma cells: proliferation, differentiation, and expansions in relation to CD45 expression. Blood Cells Mol Dis. 2004;32:293-301.
111. Link M, Warnke R, Finlay J, et al. A single monoclonal antibody identifies T-cell lineage of childhood lymphoid malignancies. Blood. 1983;62:722-728.
112. Palker TJ, Scearce RM, Hensley LL, Ho W, Haynes BF. Comparison of the CD7 (3A1) group of T cell workshop antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. Leukocyte Typing II: Human T Lymphocytes. New York, NY: Springer-Verlag; 1986:303-313.
113. Weiss LM, Crabtree GS, Rouse RV, Warnke RA. Morphologic and immunologic characterization of 50 peripheral T-cell lymphomas. Am J Pathol. 1985;118:316-324.
114. Weiss LM, Wood GS, Warnke RA. Immunophenotypic differences between dermatopathic lymphadenopathy and lymph node involvement in mycosis fungoides. Am J Pathol. 1985;120:179-185.
115. Foon KA, Todd RF. Immunologic classification of leukemia and lymphoma. Blood. 1986;68:1-31.
116. Wood GS, Abel EA, Hoppe RT, Warnke RA. Leu-8 and Leu-9 antigen phenotypes: immunologic criteria for the distinction of mycosis fungoides from cutaneous inflammation. J Am Acad Dermatol. 1986;14:1006-1013.
117. Lazarovits AI, White MJ, Karsh J. CD7- T cells in rheumatoid arthritis. Arthritis Rheum. 1992;35:615-624.
118. Lazarovits AI, Karsh J. Decreased expression of CD7 occurs in rheumatoid arthritis. Clin Exp Immunol. 1988;72:470-475.
119. Greaves MF. Monoclonal antibodies as probes for leukemic heterogeneity and hematopoietic differentiation. In: Knapp W, ed. Leukemia Markers. New York, NY: Academic Press; 1981:19.
120. Zola H. CD10 Workshop Panel report. In: Schlossman SF, Boumsell L, Gilks W, et al, eds. Leucocyte Typing V: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1995:505-507.
121. Letarte M, Vera S, Tran R, et al. Common acute lymphocytic leukemia antigen is identical to neutral endopeptidase. J Exp Med. 1988;168:1247-1253.
122. Le Bien TW, McCormack RT. The common acute lymphoblastic leukemia antigen (CD10): emancipation from a functional enigma.

Blood. 1989;73:625-635.
123. Consolini R, Legitimo A, Rondelli R, et al. Clinical relevance of CD10 expression in childhood ALL. Haematologica. 1998;83:967-973.
124. Nadler LM, Korsmeyer SJ, Anderson KC, et al. B cell origin of non-T cell acute lymphoblastic leukemia: a model for discrete stages of neoplastic and normal pre-B cell differentiation. J Clin Invest. 1984;74:332-340.
125. Hann IM, Richards SM, Eden OB, Hill FGH. Analysis of the immunophenotype of children treated on the Medical Research Council United Kingdom Acute Lymphoblastic Leukaemia Trial XI (MRC UKALLXI). Leukemia. 1998;12:1249-1255.

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