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1. Fitzgerald-Bocarlsly P, Herberman R, Hercend T, et al. A definition of natural killer cells. In: Ades E, Lopez C, eds. Natural Killer Cells and Host Defense. Basel: Karger; 1989:1.
2. Brain P, Gordon J, Willets W. Rosette formation by peripheral lymphocytes. Clin Exp Immunol. 1970;6:681-688.

3. Jondal M, Holm G, Wigzell H. Surface markers on human T and B lymphocytes. I. A large population of lymphocytes forming nonimmune rosettes with sheep red blood cells. J Exp Med. 1972;136:207-215.
4. Wolde-Mariam W, Peter J. Recent diagnostic advances in cellular immunology. Diagnost Med. 1984;7:25-32.
5. Froland S, Natvig J, Berdal P. Surface-bound immunoglobulin as a marker of B lymphocytes in man. Nature. 1971;234:251-252.
6. Bellanti J, Peters S. Diagnostic applications of immunology. In: Bellanti J, ed. Immunology III. Philadelphia: WB Saunders Co; 1985:558-569.
7. Winchester R, Fu S, Hoffman H, Kunkel H. IgG on lymphocyte surfaces: technical problems and the significance of a third cell population. J Immunol. 1975;114:1210-1212.
8. Wofsy L, Henry C, Kimura J, North J. Modification and use of antibodies to label cell surface antigens. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York, NY: WH Freeman and Co; 1980:287-304.
9. Henry C, Mishell B, Shiigi S. Hemolytic plaque assays. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York, NY: WH Freeman and Co; 1980:69-123.
10. Grabstein K, Chen Y. Cell-mediated cytolytic responses. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York, NY: WH Freeman and Co; 1980:124-137.
11. Oi V, Herzenberg L. Immunoglobulin-producing hybrid cell lines. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York, NY: WH Freeman and Co; 1980:351-372.
12. Jackson A. Basic phenotyping of lymphocytes: selection and testing of reagents and interpretation of data. Clin Immunol Newslett. 1990;10:43-55.
13. Jackson A, Warner N. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose N, Friedman H, Fahey J, eds. Manual of Clinical Laboratory Immunology. 3rd ed. Washington DC: American Society for Microbiology; 1986:226-235.
14. Nicholson J, Jones B. Immunophenotyping by flow cytometry: its use in HIV infection. Labmedica. 1989;6:21-26.
15. Kidd P, Vogt R. Report of the workshop on the evaluation of T-cell subsets during HIV infection and AIDS. Clin Immunol Immunopathol. 1989;52:3-9.
16. Landay A, Muirhead K. Procedural guidelines for performing immunophenotyping by flow cytometry. Clin Immunol Immunopathol. 1989;52:48-60.
17. Schmidt RE. Monoclonal antibodies for diagnosis of immunodefiencies. Blut. 1989;59:200-206.
18. Nicholson J. Use of flow cytometry in the evaluation and diagnosis of primary and secondary immunodeficiency diseases. Arch Pathol Lab Med. 1989;113:598-605.
19. Foucar K, Gueken JA. Clinical application of immunologic techniques to the diagnosis of lymphoproliferative and immunodeficiency disorders. Lab Medicine. 1982;13:403-413.
20. Cohen S, Weetman A. Activated interstitial and intraepithelial thyroid lymphocytes in autoimmune thyroid disease. Acta Endocr. 1988;119:161-166.
21. Smolen J, Chused T, Leiserson W, Reeves J, Alling D, Steinberg A. Heterogeneity of immunoregulatory T-cell subsets in systemic lupus erythematosus. Am J Med. 1982;72:783-790.
22. Giorgi J, Hultin L. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. Clin Immunol Newslett. 1990;10:55-62.
23. Antel J, Bania M, Noronha A, Neely S. Defective suppressor cell function mediated by T8+ cell lines from patients with progressive multiple sclerosis. J Immunol. 1986;137:3436-3439.
24. Gratama J, Naipal A, Oljans P, et al. T lymphocyte repopulation and differentiation after bone marrow transplantation: early shifts in the ratio between T4+ and T8+ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. Blood. 1984;63:1416-1423.
25. Bishop G, Hall B, Duggin G, Horvath J, Sheil A, Tiller D. Immunopathology of renal allograft rejection analyzed with monoclonal antibodies to mononuclear cell markers. Kidney Internat. 1986;29:708-717.
26. Lanier L, Le A, Phillips J, Warner N, Babcock G. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. J Immunol. 1983;131:1789-1796.
27. Landay A, Ohlsson-Wilhelm B, Giorgi J. Application of flow cytometry to the study of HIV infection. AIDS. 1990;4:479-497.
28. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael A, ed. Leucocyte Typing III: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1986:788-803.
29. Beverley P. Production and use of monoclonal antibodies in transplantation immunology. In: Touraine J, Trager J, Betuel H, eds. Transplantation and Clinical Immunology XI. Amsterdam: Excerpta Medica; 1980:87-94.
30. Rowe D, Beverley P. Characterisation of breast cancer infiltrates using monoclonal antibodies to human leucocyte antigens. Br J Cancer. 1984;49:149-159.
31. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman S, eds. Leucocyte Typing. Berlin: Springer-Verlag; 1984:25-60.
32. Dimitriu-Bona A, Burmester G, Waters S, Winchester R. Human mononuclear phagocyte differentiation antigens. I. Patterns of antigenic expression on the surface of human monocytes and macrophages defined by monoclonal antibodies. J Immunol. 1983;130(1):145-152.

33. Herrmann F, Komischke B, Odenwald E, Ludwig W. Use of monoclonal antibodies as a diagnostic tool in human leukemia. I. Acute myeloid leukemia and acute phase of chronic myeloid leukemia. Blut. 1983;47:157-163.
34. Schwinzer R. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:628-634.
35. Goyert S, Tesio L, Ashman L, et al. Report on the CD14 Cluster Workshop. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:789-794.
36. Bernstein I, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leukocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. Leukocyte Typing II: Human Myeloid and Hematopoietic Cells. New York, NY: Springer-Verlag; 1984:3:1-25.
37. Jayaram Y, Hogg N. Surface expression of CD14 molecules on human neutrophils. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:796-797.
38. Haynes B. Summary of T cell studies performed during the Second International Workshop and Conference on Human Leukocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. Leukocyte Typing II: Human T Lymphocytes. New York, NY: Springer-Verlag; 1986:1:3-30.
39. Ledbetter J, Evans R, Lipinski M, Cunningham-Rundles C, Good R, Herzenberg L. Evolutionary conservation of surface molecules that distinguish T lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. J Exp Med. 1981;153:310-323.
40. Kan E, Wang C, Wang L, Evans R. Non-covalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. J Immunol. 1983;131:536-539.
41. Knowles R. Immunochemical analysis of the T cell-specific antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. Leukocyte Typing II: Human T Lymphocytes. New York, NY: Springer-Verlag; 1986:259-288.
42. Nadler L. B cell/leukemia panel workshop: summary and comments. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. Leukocyte Typing II: Human B Lymphocytes. New York, NY: Springer-Verlag; 1986:2:3-43.
43. van Dongen J, Krissansen G, Wolvers-Tettero I, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. Blood. 1988;71:603-612.
44. Brenner M, Groh V, Porcelli F, Hochstenbach H, et al. Structure and distribution of the human γ δ T-cell receptor. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:1049-1053.
45. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. Annual Rev Immunol. 1988;6:629-662.
46. Dörken B, Moller P, Pezzutto A, Schwartz-Albiez R, Molderhauer G. B-cell antigens: CD19. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1986:34-36.
47. Evans R, Wall D, Platsoucas C, et al. Thymus-dependent membrane antigens in man: inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. Proc Natl Acad Sci USA. 1981;78(1):544-548.
48. Maddon P, Dalgleish A, McDougal J, Clapham P, Weiss R, Axel R. The T4 gene encodes the AIDS virus receptor and is expressed in the immune system and the brain. Cell. 1986;47:333-348.
49. Dalgleish A, Beverly P, Clapham P, Crawford D, Greaves M, Weiss R. The CD4 (T4) antigen is an essential component of the receptor for the AIDS virus. Nature. 1984;312:763-767.
50. Kurrle R. Cluster report: CD3. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:290-293.
51. Wood F, Warner N, Warnke R. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. J Immunol.1983;131:212-216.
52. Anderson P, Blue M-L, Morimoto C, Schlossman S. Cross-linking of T3 (CD3) with T4 (CD4) enhances the proliferation of resting T lymphocytes. J Immunol. 1987;139:678-682.
53. Eichmann K, Johnson J-I, Falk I, Emmrich F. Effective activation of resting mouse T lymphocytes by cross-linking submitogenic concentrations of the T cell antigen receptor with either Lyt-2 or L3T4. Eur J Immunol. 1987;17:643-650.
54. Gallagher P, Fazekas de St. Groth B, Miller J. CD4 and CD8 molecules can physically associate with the same T-cell receptor. Proc Natl Acad Sci USA. 1989;86:10044-10048.
55. Rudd C, Burgess K, Barber E, Schlossman S. Monoclonal antibodies to the CD4 and CD8 antigens precipitate variable amounts of CD4/CD8-associated p56-lck activity. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989: 326-327.
56. Terry L, Disanto J, Small T, Flomenberg N. Differential expression of the CD8 and Lyt-3 antigens on a subset of human T-cell receptor γ/δ-bearing lymphocytes. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:345-346.
57. Moebius U. Cluster report: CD8. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigen. New York, NY: Oxford University Press; 1989:342-343.
58. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. I. Characterization of the lymphocyte subset reactive with B73.1. J Immunol. 1983;130(5):2133-2141.
59. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. II. Studies of B73.1 antibody-antigen interaction on the lymphocyte membrane. J Immunol. 1983;130:2142-2148.

60. Perussia B, Trinchieri G, Jackson A, et al. The Fc receptor for IgG on human natural killer cells: phenotypic, functional, and comparative studies with monoclonal antibodies. J Immunol. 1984;133:180-189.
61. Schmidt R. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:517-542.
62. Schubert J, Lanier L, Schmidt R. Cluster report: CD56. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:699-702.
63. Lanier L, Testi R, Bindl J, Phillips J. Identity of Leu-19 (CD56) leukocyte differentiation antigen and neural cell adhesion molecule. J Exp Med. 1989;169:2233-2238.
64. Lanier L, Chang C, Azuma M, Ruitenberg J, Hemperly J, Phillips J. Molecular and functional analysis of human natural killer cell-associated neural cell adhesion molecule (N-CAM/CD56). J Immunol. 1991;146:4421-4426.
65. Lanier L, Le A, Civin C, Loken M, Phillips J. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. J Immunol. 1986;136:4480-4486.
66. Knapp W, Rieber P, Dörken B, Schmidt R, Stein H, van dem Borne A. Towards a better definition of human leucocyte surface molecules. Immunol Today. 1989;10:253-258.
67. Phillips J, Lanier L. Dissection of the lymphokine-activated killer phenomenon: relative contribution of peripheral blood natural killer cells and T lymphocytes to cytolysis. J Exp Med. 1986;164:814-825.
68. National Committee for Clinical Laboratory Standards. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture (H3-A3).Villanova, PA: The National Committee for Clinical Laboratory Standards; 1991.
69. National Committee for Clinical Laboratory Standards. Clinical Applications of Flow Cytometry: Quality Assurance and Immunophenotyping of Peripheral Blood Lymphocytes; Tentative Guideline. (H42-T). Villanova, PA: National Committee for Clinical Laboratory Standards; 1992.
70. Giorgi J. Lymphocyte subset measurements: significance in clinical medicine. In: Rose N, Friedman H, Fahey J, eds. Manual of Clinical Laboratory Immunology. 3rd ed. Washington DC: American Society of Microbiology. 1986:236-246.
71. Angadi C. Lack of Leu-3a epitope on T-helper (CD4) lymphocytes. J Clin Lab Anal. 1990;4:193-195.
72. Reichert T, DeBruyére M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. Clin Immunol Immunopathol. 1991;60:190-208.
73. Prince H, Hirji K, Waldbeser L, Plaeger-Marshall S, Kleinman S, Lanier L. Influence of racial background on the distribution of T cell subsets and Leu-11-positive lymphocytes in healthy blood donors. Diag Immunol. 1985;3:33-37.
74. Hicks C. Fundamental Concepts in the Design of Experiments. 3rd Edition. New York, NY: Holt, Rinehart, and Winston; 1982:222-223.

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