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Warranty Statement

These products are warranted to perform as described in their labeling and MicroScan literature when used in accordance with their instructions. THERE ARE NO WARRANTIES WHICH EXTEND BEYOND THIS EXPRESS WARRANTY AND BECKMAN COULTER DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. The sole obligation of Beckman Coulter and the purchaser's exclusive remedy for breach of this warranty shall be, at the option of Beckman Coulter to repair or replace the products. In no event shall Beckman Coulter be liable for any proximate, incidental or consequential damages in connection with the products.

References

1) Bobey, D.G., and G.M. Ederer. 1981. Rapid detection of year enzymes by using 4-methylumberyl substrates. J. Clin. Microbiol. 13:393.
2) Callaway, C.J., and L.D. Haley. 1978. Laboratory Methods in Medical Mycology, 4th ed. Centers for Disease control, Atlanta, GA.
3) Emmons, C.W., C.H. Binford, J.P. Utz, and K.J. Kwon- Chung. 1977. Medical Mycology, 3rd ed. Lea and Febiger, Philadelphia, PA.
4) Koneman, E.W., G.D. Roberts, and S.F Wright. 1985. Practical Laboratory Mycology, 3rd ed. The Williams & Wilkins Co., Baltimore, MD.
5) Kreger-van Rij, N.J.W (ed). 1984. The Yeast-A Taxonomic Study, 3rd ed. Elsevier Science Publishers B.V., Amsterdam.
6) LaRone, D.H.1987. Medically Important Fungi. A Guide to Identification. 2nd ed. Elsevier Science Publishing Co., New York, NY.
7) Lee, K.L., M.E. Reca, R.R. Watson, and C.C. Campbell. 1975, Identification of yeast phase of pathogenic fungi by the specifity of their aminopeptidase(S). Sabouraudia 13:132.
8) McGinnis, M.R. 1980. Laboratory Handbook of Medical Mycology. Academic Press New York.
9) Merz, W.G. and G.D. Roberts. 1996. Detection and Recovery of Fungi from Clinical Specimens. P.709-722. In P. R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H.Yolken (eds.), Manual of Clinical Microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
10) Pore, R.S. 1985, Prototheca Taxonomy. Mycopathol. 90:129-139.
11) Roberts, G.D. 1977. Laboratory identification of fungi, p. 147-160. In E.W. Koneman, and M.S. Britt (ed.), Clinical Microbiology. Health Education Resources, Inc., Bethesda, MD.
12) Warren, N.G., and K. C Hunter. 1996. Candida, Cryptococcus, and Other Yeasts of Medical Importance. P. 723-738. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Yolken (eds.), Manual of Clinical Microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
13) Watson, R.R. 1976. Substrates Specificities of Aminopeptidase: A Specific Method of Microbial Differentiation, p. 1-14. In J.R. Norris, and D.W. Ribbons (ed.), Methods in Microbiology, Vol. 9, Academic Press, Inc., New York, NY.

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Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

Made in USA
Beckman Coulter, Inc. EC REP
250 S. Kraemer Blvd.
Brea, CA 92821 United States
www.beckmancoulter.com